



The future of mobility is here.

**Jazzy**  
**air**

**E M P O W E R S**

The Jazzy Air® is a new era in power mobility. Whether at home, or out on the town, the unique appeal of Jazzy Air will elevate your independence and social confidence in style.

- Elevate up to 25.4 cm (10 inches) in 16 seconds
- Enjoy face-to-face social engagement
- Safely drive elevated at walking speed (5.6 km/h (3.5 mph))



[pridemobility.com.au](http://pridemobility.com.au)

# Specifications

<b>Weight Capacity</b>	136 kg (300 lbs.)
<b>Maximum Speed<sup>1</sup></b>	Up to 6.44 km/h (4.0 mph) (non-elevated position) Up to 5.6 km/h (3.5 mph) (elevated position)
<b>Ground Clearance<sup>2</sup></b>	6.15 cm (2.42") at footplate 2.54 cm (1.00") at motor bracket
<b>High-Back Seat to Floor Range</b>	58.4 cm - 83.6 cm (23" - 32.9")
<b>Turning Radius<sup>2</sup></b>	57.15 cm (22.5")
<b>Overall Size<sup>2,6</sup></b>	Length: 106.5 cm (42") with foot rigging; 96.5 cm (38") w/out foot rigging Width: 59.5 cm (23.375") at the wheels 64.8 cm (25.5") with armrests
<b>Drive Controller</b>	PG VR2
<b>Drive Wheels</b>	25.4 cm (10") flat free
<b>Caster Wheels</b>	15.2 cm (6") solid caster
<b>Suspension</b>	Active-Trac <sup>®</sup> with Mid-Wheel 6 <sup>®</sup> Technology
<b>Range<sup>1,4</sup></b>	Up to 24.6 km (15.3 miles) w/ U1 batteries Up to 29.64 km (18.4 miles) w/ 40 Ah batteries
<b>Component Weights<sup>5</sup></b>	Base: 181 lbs. (82 kg) complete assembly (without batteries) Batteries: U1-11.1 kg (24.5 lbs.) each 40 Ah-14.6 kg (32.2 lbs.) each
<b>Drivetrain</b>	Two motor, Mid-Wheel 6 <sup>®</sup> drive
<b>Brakes</b>	"Intelligent Braking" (electronic regenerative, disc park brakes)
<b>Batteries<sup>3</sup></b>	(2) 12 V, U1 (2) 40 Ah
<b>Battery Charger</b>	3.5 amp, off-board
<b>Max. Allowable Incline</b>	6°
<b>Max. Angle Capacity</b>	6°
<b>Maximum Obstacle Climbing Ability</b>	Lowered Position: 5.6 cm (2.2") Elevated Position: 4.6 cm (1.8")

- 1) Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tire type, and tire condition. This specification is subject to a variance of +10%, -5%.
- 2) Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 3) AGM or Gel-Cell type required. See VI. "Batteries and Charging."
- 4) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5) Battery weight may vary based on manufacturer.
- 6) Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA, ISO 7176 Series, and EN12184 standards. All specifications subject to change without notice.

**Pride<sup>®</sup>**  
Mobility Products  
Australia Pty. Ltd.  
— LIVE YOUR BEST —

(Tel) 03 8770 9600 • (Fax) 03 9703 2720  
20-24 Apollo Drive Hallam, Victoria 3803  
ACN 088 609 661 ABN 66 088 609 661  
[www.pridemobility.com.au](http://www.pridemobility.com.au)

# Features

- Achieve 25.4 cm (10") of power adjustable seat height in 16 seconds
- Elevate or lower seat height while on the move with the "Air" button or PG VR2 joystick drive controller
- Seat and integrated footplate design provides maximum comfort in both elevated and non-elevated positions
- Angle adjustable footplate with three height adjustments
- Choice of 41 x 41-46 cm (16" x 16"-18") depth adjustable, high-back, reclining sport seat
- Choice of 46 x 46 - 51 cm (18" x 18"-20") or 51 x 51-56 cm (20" x 20"-22") depth adjustable, fold-flat, high-back reclining sport seat
- Choice of 41 x 41-46 cm (16" x 16"-18"), 46 x 46 - 51 cm (18" x 18"-20") and 51 x 51-56 cm (20" x 20"-22") fold-flat solid seat pan
- Angle, height and depth adjustable armrest
- Bright LED marker lights
- 152 cm (60") lap belt

# Colours



Citrine  
Yellow



Sapphire  
Blue



Pink  
Topaz



Onyx  
Black



Ruby  
Red



Amber  
Orange



Pearl  
White



Gold



Silver



Rose  
Quartz

# Accessories

- Rear basket
- Cup holder
- Flag kit
- Cell phone holder
- Cane holder
- Walker holder
- Oxygen holder
- Saddle bag
- Replacement shroud (ten available colors)
- Alternative size lap belt: 127 cm, 178 cm, 203 cm, 229 cm (50", 70", 80", 90")
- XLR USB charger

# Options

- 40 amp AGM batteries
- Swing-away joystick
- Synergy<sup>®</sup> Cushions
  - Solution<sup>®</sup>
  - Solution<sup>®</sup> 1
  - Spectrum<sup>®</sup> Gel
  - Spectrum<sup>®</sup> Foam



Authorized Pride<sup>®</sup> Provider