

Jazzy[®]

SELECT 6

The Jazzy Select[®] 6 incorporates sporty design and performance features for a truly unique ride experience. Pride's patented Active-Trac[®] Suspension combined with six wheels on the ground and an innovative articulating beam allow for superior performance indoors and out, with added stability over uneven terrain and ramp transitions. Mid-wheel drive provides the tightest turning radiuses in their class and caster wheels absorb the shock of obstacles, offering outstanding handling and control.



*Jazzy Select[®] 6
shown in Jazzy Red with
high-back seat and Jazzy[®]
Armor[®] (standard)*

SPECIFICATIONS

Weight Capacity	136 kg (300 lbs.)
Maximum Speed	Up to 6.8 km/h (4.2 mph)
Ground Clearance	4.45 cm (1.75") at motors and at center frame
Turning Radius	54.6 cm (21.5")
Overall Length	90.2 cm (35.5")
Overall Width	59.7 cm (23.5")
Drive Wheels	25.4 cm (10")
Front Wheels	15.24 cm (6") solid casters
Rear Wheels	15.24 cm (6") solid casters
Suspension	Active-Trac® (ATX)
Drivetrain	Two-motor, in-line, mid-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 20.8 km (12.92 miles)
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	11 kg (24.5 lbs.) each
Base Weight	45.4 kg (100 lbs.)
Standard Seat Weight	19.5 kg (43 lbs.) (comfort, high-back)
Standard Seat	46 x 46-20 cm (18" x 18-20")
Max. Seat Dimensions	51 x 51 cm (20" x 20")

FEATURES

- Jazzy® Armor® includes shroud and controller guards to protect against daily wear and tear
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 152 cm (60") lap belt

ACCESSORIES*

- Pride® rear bag (includes crutch/ walking stick holder)
- Walker holder (must have Pride® rear bag to attach)
- Oxygen tank holder (must have Pride® rear bag to attach)
- Saddlebag
- Cup holder
- Mobile phone holder

* Please note, accessory bracket may be needed for attaching rear bag.

*Jazzy Select® 6
shown with power
elevating seat option*



COLOUR CHOICES



Jazzy
Red



Jazzy
Blue

OPTIONAL SEATING

- Synergy® Seat
- Solid seat pan

Pride
Mobility Products
Australia Pty. Ltd.

20-24 Apollo Drive Hallam, Victoria 3803
ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au



The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.



AU-JZSelect6-111017