

JET™ 3

U L T R A

The Jet™ 3 Ultra delivers a broad range of standard performance features like Active-Trac® Suspension and a new, compact 50 amp, PG VSI controller on a very manoeuvrable frame. Plus, the style is fresh, exciting and contemporary. These things combine with a truly economical price to make the Jet 3 Ultra an excellent overall value in personal mobility.



Jet™ 3 Ultra
shown in Red

21 Healey Rd., Dandenong, Victoria 3175

ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

Pride[®]
Mobility Products
Australia Pty. Ltd.

— LIVE YOUR BEST® —

SPECIFICATIONS

Weight capacity	136 kg (300 lbs.)
Maximum speed	Up to 6.44 km/h (4 mph)
Ground clearance	7 cm (2.75")
Turning radius	50.80 cm (20")
Overall length	102.87 cm (40.5")
Overall width	57.15 cm (22.5")
Wheels	Drive: 25.4 cm (10") pneumatic or solid
Front anti-tips	15.24 cm (6") solid
Rear castors	20.32 cm (8") solid
Suspension	Active-Trac®
Drive train	Two-motor, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc park brake)
Standard electronics	50 amp, PG VSI controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp (standard)
Per-charge range	Up to 17.7 km (11 miles)
Battery requirements	(2) 12 volt, deep cycle
Battery size	32 AH
Battery weight	11 kg (24.5 lbs.) each
Base weight	41.74 kg (92 lbs.)
Standard seat weight	19 kg (42.5 lbs.) (high-back)
Std. seat dimensions	46 cm x 46 cm (18" x 18")
Warranty	2-year limited

COLOUR CHOICES



Jet Red



Jet Blue

FEATURES

- Compact size for superior handling in tight quarters
- Standard Active-Trac® Suspension
- 136 kg (300 lbs.) weight capacity
- Utilizes the trapeze bar seat mount for quick adjustability and easy interchange of a broad range of seats
- Convenient front battery access
- Curved rear castor beam provides excellent stability
- Maximum speed up to 6.44 km/h (4 mph)
- Up to 17.7 km (11 miles) of range per charge
- Standard 50 amp PG VSI controller
- 50.80 cm (20") turning radius

OPTIONS

- Power seat
- Cane/crutch holder
- Cup holder
- Walker holder
- Heavy-duty dust cover
- Swing-away footrests
- Rear basket
- Elevating leg rests
- Oxygen tank holder
- Heel loops
- Contoured seat with solid seat pan
- Lap belt

Jet™ 3 Ultra shown in Red



21 Healey Rd., Dandenong, Victoria 3175

ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

The information contained herein is correct at the time of publication; due to the Pride's commitment to constant improvement and development, We reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.



PRIDE
PERFORMANCE MOBILITY

— LIVE YOUR BEST® —