

SPECIFICATIONS

	C1	C5²	C6²	C-101	LC-101 Euro Leather	D30
Back Width:	51 cm (20")	71 cm (28")	71 cm (28")	72 cm (28.5")	58 cm (23")	74 cm (29")
Seat-To-Floor:	41 cm (16")	46 cm (18")	48 cm (19")	47 cm (18.5")	47 cm (18.5")	48 cm (19")
Seat Depth:	51 cm (20")	46 - 51 cm (18 - 20")	53 cm (21")	51 cm (20")	51 cm (20")	56 cm (22")
Seat Width:	46 cm (18")	56 cm (22")	56 cm (22")	51 cm (20")	51 cm (20") ¹	56 cm (22")
Top of Back to Seat:	66 cm (26")	66 - 71 cm (26 - 28")	66 - 71 cm (26-28")	70 cm (27.5")	70 cm (27.5")	66 cm (26")
Weight Capacity:	113 kg (250 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	147 kg (325 lbs.)	147 kg (325 lbs.)	170 kg (375 lbs.)
Colour Options:	Arctic Blue Chocolate Oatmeal	Arctic Blue Chocolate Oatmeal	Arctic Blue Chocolate Oatmeal Levante Vinyl: Royal Blue	Cocoa Fern Sandal	Leather: Beige Brown Black Red	Ocean Burgundy Spruce Latte Levante Vinyl: Burgundy Royal Blue Pine Mountain

	560	660	LL-805	T-3	LC-440 Power Stand Assist
Back Width:	71 cm (28")	71 cm (28")	77.5 cm (30.5")	69 cm (27")	53.5 cm (21")
Seat-To-Floor:	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	49 cm (19")
Seat Depth:	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	46 cm (18")
Seat Width:	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	54.5 cm (21.5")
Top of Back to Seat:	63.5 cm (25")	63.5 cm (25")	68.5 cm (27")	68.5 cm (27")	69 cm (27")
Weight Capacity:	170 kg (375 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	170 kg (375 lbs.)	125 kg (275 lbs.)
Colour Options:	Ocean Burgundy Spruce	Latte	Cocoa Fern Sandal	Black (Leather)	Cocoa Fern Sandal

BENEFITS

- Engineered furniture laminate/hardwood frame
- Patented, quiet & smooth lift system
- Large side pocket ¹
- Hand control can be easily moved to left-hand side
- Integrated emergency battery backup (batteries not included)
- External transformer for easy access to battery backup

Note: Measurements may vary +/- 2.5 cm (1"). Specifications subject to change without prior notice.

¹ Large side pocket not available on LC-440.

² C5 and C6 also available in Pride's custom fabric selections range. Please see "Custom Fabric Selections" brochure for further details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



Pride® Lift Chairs
Experience the Luxury

Pride Lift Chairs

Experience the Luxury

Pride® Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance and are available in several beautiful fabrics and styles to fit any décor. The state-of-the-art lift system utilised in Pride Lift Chairs will recline you further back than any other lift chair. Plus, Pride Lift Chairs are equipped with an emergency battery backup system should the power fail in your home.

LIFT CHAIR DESCRIPTIONS

C1
3-position, petite to small chaise lounger

C5
3-position, medium to large chaise lounger with fully removable upholstery components

C6
Infinite-position, medium to large chaise lounger with fully removable upholstery components and dual motors which enable independent recline/footrest function

C-101
3-position, medium to large lift chair

LC-101 Euro Leather
3-position, medium to large lift chair

D30
3-position, medium to large chaise lounger

560
3-position, small to medium chaise lounger



660
Infinite-position, small to medium chaise lounger with dual motors that enable independent recline/footrest function

LL-805
2-position, medium to large wall hugger chaise lounger which requires only 10 cm (4") of space between wall and chair

T-3
Infinite-position, medium to large lift chair featuring Trendelenburg position (elevates legs above head for greater relaxation) and zero gravity position

LC-440 Power Stand Assist³
Power seat lift

AVAILABLE FABRICS ¹



LIFT CHAIR MODELS



C1

Accented Pillow-back Design
Shown in Arctic Blue



C5

Sewn Pillow-back Design
Shown in Arctic Blue



C6

Sewn Pillow-back Design
Shown in Chocolate



C-101

Button-back Design
Shown in Sandal



LC-101 Euro Leather

Split-back Design
Shown in Beige Leather



D30

Split-back Design
Shown in Burgundy



560

Pillow-back Design
Shown in Spruce



660

Pillow-back Design
Shown in Latte



805

Sewn Split-back Design
Shown in Cocoa



T-3

Mongolian-back Design
Shown in Black (Leather)



LC-440 Power Stand Assist²

Shown in Sandal



PLEASE NOTE: 2-position is a TV recline.
3-position is a full recline.
Infinite-position is a full recline with independent leg rest function.

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

1. Not all fabrics are available for every model. Please see Specifications for details.
2. This chair does not recline.