

# Pride<sup>®</sup> POWER RECLINERS SPECIFICATIONS



	Petite to Small	Medium to Large	Medium to Large	Medium to Large	Medium to Large	Medium	Medium to Large	Small to Medium	Small to Medium	Medium to Large	Medium to Large	Medium	Medium	Heavy Duty
	C1	C5	C6	C-101	LC-101 Euro Leather	LC-107	D30	560	660	805	T-3	SR-358M	SR-525M	LC-358XL
	3-Position Lift Chair with fully removable upholstery components <b>163 cm &amp; below (5'4" &amp; below)</b>	3-Position Chaise Lounger with fully removable upholstery components <b>163 cm - 183 cm (5'4" - 6'0")</b>	Infinite-Position Chaise Lounger with fully removable upholstery components <b>163 cm - 183 cm (5'4" - 6'0")</b>	3-Position Chaise Lounger <b>163 cm - 183 cm (5'4" - 6'0")</b>	3-Position Chaise Lounger with European style rounded arm design <b>163 cm - 183 cm (5'4" - 6'0")</b>	Infinite-Position Chaise Lounger <b>5'4" - 6'0"</b>	3-Position Chaise Lounger <b>163 cm - 183 cm (5'4" - 6'0")</b>	3-Position Chaise Lounger <b>163 cm (5'4") &amp; below</b>	Infinite-Position Chaise Lounger <b>163 cm (5'4") &amp; below</b>	2-Position Chaise Lounger Wall Hugger that requires only 10 cm (4") of space between the wall and chair <b>163 cm - 183 cm (5'4" - 6'0")</b>	Infinite-Position Chaise Lounger featuring Trendelenburg position (elevates legs above head for greater relaxation) and zero-gravity position <b>163 cm - 183 cm (5'4" - 6'0")</b>	3-Position Chaise Lounger <b>5'4" - 6'0"</b>	Infinite-Position Zero-gravity Trendelenburg Chaise Lounger <b>5'4" - 6'0"</b>	3-Position Chaise Lounger <b>5'6" - 6'1"</b>
Weight Capacity	113 kg (250 lbs.)	170 kg (375 lbs.)	136 kg (300 lbs.)	147 kg (325 lbs.)	147 kg (325 lbs.)	147 kg (325 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	136 kg (300lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	170 kg (375 lbs.)	227 kg (500 lbs.)
Back Width	51 cm (20")	71 cm (28")	71 cm (28")	72 cm (28.5")	58 cm (23")	74 cm (29")	74 cm (29")	71 cm (28")	71 cm (28")	76 cm (30")	69 cm (27")	71 cm (28")	69 cm (27")	91 cm (36")
Seat-to-Floor	41 cm (16")	46 cm (18")	48 cm (19")	47 cm (18.5")	47 cm (18.5")	51 cm (20")	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")
Seat Depth	51 cm (20")	46 cm - 51 cm (18-20")	53 cm (21")	51 cm (20")	51 cm (20")	53 cm (21")	56 cm (22")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	51 cm (20")	51 cm (20")	56 cm (22")
Seat Width	46 cm (18")	56 cm (22")	56 cm (22")	51 cm (20")	51 cm (20")	51 cm (20")	56 cm (22")	48 cm (19")	48 cm (19")	51 cm (20")	51 cm (20")	51 cm (20")	51 cm (20")	66 cm (26")
Top of Back to Seat	66 cm (26")	66 cm - 71 cm (26-28")	66 cm - 71 cm (26-28")	70 cm (27.5")	70 cm (27.5")	70 cm (27.5")	66 cm (22")	63.5 cm (25")	63.5 cm (25")	69 cm (27")	68.5 cm (27")	62 cm (24.5")	67 cm (26.5")	74 cm (29")
Back Style	Blown Fibre Accented Pillow-back	Blown Fibre Sewn Pillow-back	Blown Fibre Sewn Pillow-back	Blown Fibre Button-back	Blown Fibre Split-back	Blown Fibre Sewn Split-Back	Blown Fibre Split-back	Blown Fibre Pillow Back	Blown Fibre Pillow Back	Blown Fibre Sewn Split-Back	Blown Fibre Mongolian-back	Blown Fibre Button-back	Blown Fibre Split-back	Blown Fibre Button-back
Fully Padded Chaise	Yes <sup>1,2</sup>	Yes <sup>1,2</sup>	Yes <sup>1,2</sup>	Yes <sup>1,2</sup>	Yes <sup>1,2</sup>	Yes	Yes <sup>1</sup>	Yes <sup>1</sup>	Yes <sup>1</sup>	Yes	Yes	Yes	Yes	Yes
Integrated Battery Backup <sup>3</sup>	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Recline Positions	3	3	Infinite	3	3	Infinite	3	3	Infinite	2	Infinite	3	Infinite	3
Standard Arm Covers	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Standard Head Cover	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Colour Options	Arctic Blue, Chocolate, Oatmeal	Arctic Blue, Chocolate, Oatmeal	Arctic Blue, Chocolate, Oatmeal <b>Levante (vinyl):</b> Royal Blue	Cocoa, Fern, Sandal	<b>Leather:</b> Beige, Brown, Black, Red	Pecan, Latte, Spruce, Burgundy, Ocean	Ocean, Burgundy, Spruce, Latte <b>Levante (vinyl):</b> Burgundy, Royal Blue, Pine Mountain	Ocean, Burgundy, Spruce	Latte	Cocoa, Fern, Sandal	Black (Leather)	Ember, Timber	Ember, Timber	Burgundy (Vinyl), Latte

## ACCESSORIES



Non-Skid Leg Levelers



Head Pillow



Lumbar Pillow



Footrest Extension



Heat and Massage



Protective Pad



Cleaning Kit

1. Removable fully padded chaise.
2. All upholstery components are fully removable.
3. Batteries are not included.
4. This chair does not recline.

Wood Frame: Engineered Furniture Laminate/Hardwood Frame\*

Not all accessories are available for all products. Please speak with your Provider for details and availability.

- \* C5 & C6 frames are made of steel, not Laminate/Hardwood.
- \*\* Not all fabrics are available for every model. Please see Specifications for details.

**Note:** Measurements may vary +/- 2.5 cm (1"). Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

## LIFT CHAIR POSITION GUIDE

The lifted, seated and reading positions are achievable on all Pride<sup>®</sup> Lift Chairs

	Lifted	Seated	Reading	Napping	Relaxing	Sleeping	Trendelenburg	Zero Gravity
2-Position	The 2-position TV recline is achievable on all Pride <sup>®</sup> Lift Chairs. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.							
3-Position	The 3-position full recline is available on most Pride <sup>®</sup> Lift Chair models and puts you in a position ideal for napping or relaxing.							
Infinite Position	The infinite-position sleep recline is perfect for taking a nap - whether quick or extended - or even for a full night's sleep. Infinite-position Pride <sup>®</sup> Lift Chairs enable the backrest and footrest to move independently from one another so that a variety of positions can be reached.							
Trendelenburg	The Trendelenburg position, elevates one's feet higher than his or her head for increased comfort and relaxation.							
Zero gravity	The zero gravity position, which can be reached on infinite-position Pride <sup>®</sup> Lift Chair models, positions one's thighs and lower legs higher than his or her torso.							

## FABRIC SELECTIONS\*\*



### SUEDE




### LEVANTE VINYL



### LEATHER



### DURASOFT FABRICS



20-24 Apollo Drive Hallam, Victoria 3803  
ACN 088 609 661 ABN 66 088 609 661

Tel: 03 8770 9600 • Fax: 03 9703 2720  
[pridemobility.com.au](http://pridemobility.com.au)

