

## SPECIFICATIONS

**Product Weight:** 12 kg (26.5 lbs.)<sup>1</sup>

**Weight Capacity:** 113 kg (250 lbs.) (159 kg (350 lbs.) with HD configuration)

**Seat-to-Floor Heights:** 38 cm - 49.5 cm (15"-19.5") front  
36.8 cm - 48.3 cm (14.5"-19") rear  
31.8 cm - 42 cm (12.5"-16.5") front (Ultra Hemi)  
31.8 cm - 48.3 cm (12.5"-19") rear (Ultra Hemi)

**Seat Sizes:** Width: 35.5 cm - 53.3 cm (14"-21")  
(HD: 45.7 cm - 56 cm (18" - 22"))  
Depth: 35.5 cm - 50.8 cm (14"-20")  
(HD: 45.7 cm - 50.8 cm (18" - 20"))  
(Available in 2.5 cm (1") increments)

## COLOURS

### METALLIC



Cosmic Silver Pacific Blue Plum Crazy



Purple Haze Rainforest Green

### NON-METALLIC



Yellow Sunshine Blue Sky Green Mantis Pink Flamingo



Red Hot Black Abyss Orange Crush Arctic White

## OPTIONS

**REAR WHEELS:** 51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" & 26") Spinergy wire spoke  
56 cm, 61 cm, 63.5 cm & 66 cm (22", 24" 25" & 26") Spinergy Light Extreme (available in red, yellow, blue, white, pink, green, orange and black)  
61 cm, 63.5 cm & 66 cm (24", 25" & 26") Spinergy Sport Light Extreme (available in red, yellow, blue, white, pink, green, orange and black)  
51 cm, 56 cm & 61 cm (20", 22" & 24") Composite mag  
61 cm (24") heavy duty Skyway Tuff composite mag wheels  
Magic Wheels  
One-arm drive system

**TYRES:** 56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") 110 PSI Primo tyre  
61 cm & 66 cm (24" & 26") Aggressive knobby tyre  
51 cm, 56 cm, 61 cm & 66 cm (20", 22", 24" & 26") Everyday grey tyre  
51 cm, 56 cm, 61 cm & 63.5 cm (20", 22", 24" & 25") Polyurethane tyre (61 cm (24") also available in light grey)  
66 cm (26") Kenda mountain bike knobby tyres (available in black or grey)  
61 cm & 66 cm (24" & 26") 110 PSI Primo coloured tyres (available in red, blue or yellow)  
56 cm, 61 cm & 63.5 cm (22", 24", 25" & 26") Schwalbe Marathon Plus and Marathon Plus Evolution

**HANDRIMS:** 51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" & 26") Aluminum  
56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Vinyl-coated  
51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" and 26") Projection (vertical or 45° oblique)  
56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Natural-Fit (available in standard or super grip)  
61 cm (24") Natural-Fit LT (available in standard or super grip)  
61 cm (24") Spinergy FlexRims  
51 cm, 56 cm, 61 cm (20", 22" & 24") Chrome

**AXLE PLATES:** Standard (0°, 1.5°, 3° camber options)  
Amputee (0°, 1.5°, 3° camber options)

**CASTERS:** 10 cm, 13 cm, 15 cm & 20 cm (4", 5", 6" & 8") composite  
10 cm, 13 cm & 15 cm (4", 5" & 6") soft roll  
13 cm (5") aluminum  
10 cm (4") light-up casters  
10 cm (4") coloured casters (available in red, blue, green, orange, yellow or clear)  
10 cm, 13 cm, 15 cm & 20 cm (4", 5", 6" & 8") black, smooth 3-spoke tyres  
15 cm (6") pneumatic and 20 cm (8") pneumatic or flat-free polyurethane casters

**FORKS:** 13 cm (5") or 20 cm (8") quick release forks  
13 cm (5") or 20 cm (8") bolt-on forks  
Frog Legs suspension forks  
60 mm forks (for use with 20 cm (8") pneumatic casters only)  
Caster pin locks

**WHEEL LOCKS:** Pull-to-lock Push-to-lock Scissor locks Pull-to-lock with extensions Push-to-lock with extensions Hill holder Cable-linked, unilateral wheel locks

**ARMRESTS:** 30.5 cm (12") tubular swing-away (desk length)  
Height-adjustable drop-in (desk or full length)  
Locking or non-locking, flip-up (depth adjustable only)  
Height-adjustable, flip-back, dual post armrests (desk or full length)

**BACK CANES:** Straight, with or without push handles  
Angle-adjustable, with or without push handles  
8 degree with push handles  
Straight reinforced  
Depth-adjustable frame option with fold forward back canes

**BACKS:** Standard nylon sling with convenience pocket  
Tension-adjustable upholstery  
Synergy® Shape Back

**FRONT RIGGINGS:** Aluminum, bi-directional swing-away leg rests, 60°, 70°, 70° tapered, 80° and 90°  
Elevating leg rests with calf pads  
Comfort adjust footplates  
Angle-adjustable foot plates  
Nylon heel loops

**ACCESSORIES:** Height-adjustable, removable anti-tip wheels  
Composite side guards  
Aluminum side guards  
Lap belt  
Quantum® backpack  
Spoke protectors  
Transit securement package  
Gel tubular arm cover  
Case Logic underseat organizer bag  
Combat, cane and crutch holder  
Case Logic cup holder

<sup>1</sup> Standard configuration: 40.6 cm x 40.6 cm (16" x 16"), 48.3 cm (19") backrest, sling upholstery, push-to-lock brakes, 61 cm (24") spoke wheels, 3.5 cm (1-3/8") pneumatic tyres, 15.2 cm (6") casters (not including leg rests)

**NOTE:** Not all options available for all configurations. Please check order form for compatibility.

# Litestream<sup>®</sup> XF



“Experience the performance and durability that I demand.”

– Bryan Anderson  
Quantum Rehab<sup>®</sup> Spokesperson



Quality System  
Quality Endorsed Company  
ISO 9001 QIC20029  
SAI Global

20-24 Apollo Drive, Hallam, Victoria 3803

ACN 088 609 661 ABN 66 088 609 661

Tel (03) 8770 9600 • Fax (03) 9703 2720 • [www.pridemobility.com.au](http://www.pridemobility.com.au)

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

**QUANTUM**  
LIFE BEYOND LIMITS™

[www.quantumrehab.com.au](http://www.quantumrehab.com.au)



HD configuration

Ultra-hemi configuration

## The Litestream® XF

Ultralight manual wheelchair uses a unique, innovative, rigidising double X folding frame that delivers the versatility and performance of a rigid frame while folding to a width of only 23 cm (9") for easy transport. The Litestream XF accepts a broad selection of accessories to meet a wide variety of needs and personal preferences while remaining an excellent overall value. The Litestream XF is also available in an HD configuration that weighs just 18 kg (39 lbs.) and features a cane-to-cane rigidising bar for added stability.

### REAR WHEELS



51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" & 26") high-performance Spinerger wheels feature double-wall rims with straight-pull spokes for reduced weight and optimal durability

56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Spinerger Light Extreme 12 spoke PBO Fiber wheels

AVAILABLE COLOURS:



61 cm, 63.5 cm & 66 cm (24", 25" & 26") Spinerger Sport Light Extreme 24 spoke PBO Fiber wheels

AVAILABLE COLOURS:



51 cm, 56 cm & 61 cm (20", 22" & 24") Skyway Tuff composite mag wheels are exceptionally durable and low maintenance

Also available:

61 cm (24") heavy duty Skyway Tuff composite mag wheels with aluminum handrim



Magic Wheels feature two gears to address real-life situations – 1st gear, for level ground and 2nd gear, providing a 100-percent increase in hill climbing torque



61 cm (24") lightweight, quick-release, one-arm drive system with aluminum handrim

### CASTERS



10 cm, 13 cm, 15 cm & 20 cm (4", 5", 6" & 8") composite casters feature precision bearings and long-life poly tyres



10 cm, 13 cm & 15 cm (4", 5" & 6") soft-roll casters smooth lumps and bumps in everyday use



13 cm (5") aluminum casters with precision bearings and polyurethane tyres



10 cm (4") light-up casters



10 cm (4") coloured casters available in red, blue, green, yellow, orange or clear



10 cm, 13 cm, 15 cm & 20 cm (4", 5", 6" & 8") black, smooth 3-spoke tyres



15 cm (6") pneumatic casters



20 cm (8") pneumatic or flat-free polyurethane casters

### HANDRIMS



51 cm, 56 cm, 61 cm, 63.5 cm & 66 cm (20", 22", 24", 25" & 26") Projection Handrim (vertical or 45° oblique)



61 cm (24") Spinerger FlexRims



56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" and 26") Natural-Fit and 61 cm (24") Natural-Fit LT Handrim (available in standard or super grip)

(Aluminum, vinyl-coated & chrome handrims also available but not shown)

### TYRES



61 cm & 66 cm (24" & 26") Aggressive Knobby Tyre



51 cm, 56 cm, 61 cm & 66 cm (20", 22", 24" & 26") Everyday Grey Tyre



51 cm, 56 cm, 61 cm & 63.5 cm (20", 22", 24" & 25") Polyurethane Tyre (61 cm (24") also available in light grey)



56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") 110 PSI Primo Tyre



61 cm & 66 cm (24" & 26") 110 PSI Primo Coloured Tyre



66 cm (26") Kenda Mountain Bike Knobby Tyres



56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Schwalbe Marathon Plus Evolution



56 cm, 61 cm, 63.5 cm & 66 cm (22", 24", 25" & 26") Schwalbe Marathon Plus

### FORKS



Frog Legs suspension fork



20 cm (8") fork



13 cm (5") fork



Quick-release 60 mm fork

(for use with 20 cm (8") pneumatic or flat-free polyurethane casters only)



Caster pin lock (available as bolt-on only)

### WHEEL LOCKS



Pull-to-lock



Scissor lock



Pull-to-lock with extension



Push-to-lock with extension



Hill holder



Cable-linked, unilateral wheel locks (for use with one-arm drive system)

### LEG RESTS



Elevating bi-directional leg rests with swing-away calf pads



60°, 70°, 70° tapered, 80° and 90° aluminum bi-directional swing-away leg rests



10-degree comfort adjust footplates in 11 cm, 13 cm & 15 cm (4.25", 5" & 6") widths



Angle-adjustable footplates in 14 cm, 16.5 cm & 19 cm (5.5", 6.5" & 7.5") widths



Nylon heel loops

### BACK OPTIONS



Depth-adjustable frame option with fold forward back canes



Tension-adjustable upholstery



Standard nylon sling



Synergy® SHAPE back

### TRANSIT SECUREMENT



The Transit Securement Package allows convenient securement of the Litestream® XF for safe transport.



### ACCESSORIES AND PERSONALISATION



Height-adjustable anti-tips



Composite side guards



Aluminum side guards



Spoke protector



Combat, cane & crutch holder



Case Logic cup holder



Quantum® backpack



Case Logic underseat organizer bag

### ARMRESTS



30.5 cm (12") tubular, swing-away, desk-length armrests (shown with optional gel arm covers)



Height-adjustable, flip-back, dual post armrests (desk or full length)



Locking or non-locking, flip-up armrests (for use with depth-adjustable option only)



Height-adjustable drop-in armrests (desk or full length)