



NEW!

Pride Lift Chairs

MEDIA Chair

Experience the Luxury

Pride
Mobility Products
Australia Pty. Ltd.

PRIDE
PERFORMANCE MOBILITY
— LIVE YOUR BEST™ —

NEW!

MEDIA Chair

Pride® Lift Chairs | *Experience the Luxury*

The LC-900 Media Chair brings the excitement of the theatre to the comfort of your own home. The independent accessory wedge provides two cup holders and ample storage space. The LC-900 incorporates a patented power lift and recline system for maximum comfort. Available in stylish Lexis Black Vinyl, the LC-900 is a welcome addition to any family room, home theatre or living room.


SPECIFICATIONS

	LC-900^{1,3,4}
	3-Position Chaise Lounger 5'5" - 6'0"
Weight Capacity	170 kg (375 lbs.)
Seat to Floor	52 cm (20.5")
Seat Depth	56 cm (22")
Seat Width	56 cm (22")
Top of Back to Seat	58.5 cm (23")
Back Style	Blown Fibre Sewn Pillow Back
Fully Padded Chaise	Yes
Distance From Wall	43 cm (17")
Overall Chair Width	86 cm (34")
Chair Weight	54 kg (120 lbs.) ²
Warranty	2-Year Limited



FABRIC

The LC-900 Media Chair is only available in **LEXIS BLACK VINYL**.



BLACK

1. Media wedge dimensions: 30.5 cm (12") W x 67 cm (26.5") H x 81 cm (32") L.
2. Media wedge weight: 17 kg (38 lbs.)
3. No side pockets.
4. Media wedge sold separately and not included with purchase of the LC-900.



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

20 Apollo Drive
Hallam, Victoria 3803
www.pridemobility.com

Tel: 03 9706 4611
Fax: 03 9706 4622

Pride®
Mobility Products
Australia Pty. Ltd.

PRIDE®
PERFORMANCE MOBILITY
— LIVE YOUR BEST™ —