


Q-Logic

Drive Control System

The fully programmable Q-Logic Drive Control System is packed with unique features including an innovative display that provides users with more information than ever before. Speed, battery condition, odometer information, and time are presented in an easy-to-read colour display. In addition, Q-Logic offers more programming options allowing for a new level of customisation.



www.quantumrehab.com

Pride[®]
Mobility Products
Australia Pty. Ltd.

AVAILABLE MODULES

- Hand Control
- Stand Alone Joystick
- Attendant Control
- Enhanced Display
- Power Base Module
- Sip-N-Puff Module
- ECU Module
- Advanced Actuator Module
- Hand-held Programmer
- PC Based Programmer



Stand Alone Joystick/
Attendant Control



Enhanced
Display



Power Base
Module



Hand-held
Programmer

FULL-COLOUR DISPLAY SCREENS



Drive Screen

- Speedometer
- Battery Charge Indicators
- Profile Indicator
- Real Time Clock
- Speed Knob Position Indicator
- Odometer
- Trip Odometer
- Drive Status Indicator



Seating Screen

- Angle Indication
- Illuminated Panels Indicate Active Seating Functions
- Drive Status Indicator
- Battery Indicator
- Profile Indicator
- Real Time Clock



Speed Display



Actuator Display



Latched Screen



Auxiliary Screen

"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".

QUANTUM
REHAB
INNOVATIVE REHAB SOLUTIONS

21 Healey Rd.
Dandenong, Victoria 3175
www.pridemobility.com

ACN 088 609 661 ABN 66 088 609 661

Tel: (03) 9706 4611
Fax: (03) 9706 4622

