

EDGE 3

EDGE 3 shown in Blue Jean
with iLevel® power seat elevation
and Q-Logic 3 with Bluetooth®

AVAILABLE WITH
30.5 CM (12") iLEVEL®
POWER SEAT ELEVATION AT
7.24 KM/H
(4.5 MPH)



THE EDGE 3 IS PROUDLY
MADE IN THE
USA
OF FOREIGN AND DOMESTIC PARTS


QUANTUM®

quantumrehab.com



Specifications

DRIVE WHEELS	Standard: 35.6 cm (14") Solid Optional: 35.6 cm (14") Pneumatic
CASTOR WHEELS	Front: 15.2 cm (6") Solid Rear: 15.2 cm (6") Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 9.6 km/h (6 mph) Up to 7.24 km/h (4.5 mph) w/iLevel activated
GROUND CLEARANCE	6.9 cm (2.7") (battery tray)
TURNING RADIUS ²	52.1 cm (20.5")
OVERALL LENGTH ²	82 cm (32.23") without front riggings
BASE WIDTH	62.2 cm (24.4")
SEATING SIZES	
Synergy®/TRU-Balance® 3 Static Seat	Width: 25.4 - 55.9 cm (10-22") Depth: 25.4 - 55.9 cm (10-22")
TRU-Balance® 3 Power Tilt	Width: 30.5 - 55.9 cm (12-22") Depth: 30.5 - 55.9 cm (12-22")
TRU-Balance® 3 Tilt and iLevel®	Width: 30.5 - 55.9 cm (12-22") Depth: 30.5 - 55.9 cm (12-22")
TRU-Balance® Flex Seating	Width: 38.1 - 55.9 cm (15-22") Depth: 30.5 - 50.8 cm (12-20")
SEAT-TO-FLOOR HEIGHTS¹⁰	
Synergy®/TRU-Balance® 3 Static Seat ⁵	41.9 - 47 cm (16.5-18.5")
TRU-Balance® 3 Power Tilt ⁵	44.5 - 49.5 cm (17.5-19.5")
TRU-Balance® 3 Tilt & iLevel® ⁵	44.5 - 49.5 cm (17.5-19.5")
TRU-Balance® Flex Seating	41.9 - 49.5 cm (16.5" - 19.5")
MANUAL RECLINE ⁶	Yes
BATTERY SIZE ³	NF-22 (2)
BATTERY WEIGHT ⁷	17.24 kg (38 lbs.) each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel 6®
WEIGHT CAPACITY ⁴	136 kg (300 lbs.)
BASE WEIGHT	63 kg (138.6 lbs.)


1) Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 3 mm (1/8") with 46 x 46 cm (18" x 18") seat, fixed towers, 136 kg (300 lbs.) weight capacity and solid tires (add 10 mm (3/8") with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
 7) Battery weight may vary +/- 1.4 kg (3 lbs.) based on manufacturer.
 8) All Quantum Rehab® units can be ordered with specialty controls.
 9) Call for availability
 10) 41.9 cm (16.5") seat-to-floor height is not available with power tilt or iLevel®. 49.5 cm (19.5") seat-to-floor height is not available with static seat or power recline.

Features


- Upgraded SRS (Smooth Ride Suspension)
 - Gas-charged monotube shocks with springs mounted in a coil-over assembly for a smooth, comfortable ride and advanced stability
 - The comfort of soft springs with the control of real dampening shocks
 - All shocks function independently of each other
 - Same springs and shocks as the 4Front™
 - Close simulation to automotive suspension
- 7.24 km/h (4.5 mph) with iLevel®
- Fender lights and USB charger included
- Bold new colours
- Made in the USA of foreign and domestic parts



7 New Colours*



HYDROGRAPHICS



CARBON FIBER

* Visit quantumrehab.com to see all 18 available colours.