

Compact Design with Maximum Stability

FLASH

Pride Mobility Products

Foldable Power Wheelchair



It folds in a flash! The innovative Flash Power Wheelchair offers easy portability while meeting the demands of daily life. Rear anti-tip wheels combine with a compact frame to ensure stable outdoor performance and smooth indoor manoeuvring. The Flash is the ultimate economical folding power chair.



Folds with Ease

Pride[®]
Mobility Products
Australia Pty. Ltd.

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au



SPECIFICATIONS

Weight capacity	113 kg (250 lbs.)
Maximum speed¹	Up to 6.4 km/h (4 mph)
Ground clearance²	10 cm (4") (at motor)
Turning radius²	81.3 cm (32")
Overall length²	82.5 cm (32.5")
Overall width²	61.5 cm (24.25")
Drive wheels	31.75 cm x 5.7 cm (12.5" x 2.25") (pneumatic)
Anti-tips	3.8 cm (1.5") (solid)
Casters	20.3 cm x 5 cm (8" x 2") (solid)
Suspension	N/A
Drive train	Two-motor, rear-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	50 amp PG VR2 controller
Standard seat	Sling
Specialty controls	No
Battery charger	Off-board, 5 amp
Per-charge range¹	Up to 17.28 km (10.74 miles)
Battery requirements	(2) 12 volt deep-cycle
Battery size (largest)	35 AH
Battery weight³	10.6 kg (23.5 lbs.) each
Base weight	42 kg (92.5 lbs.) (without batteries)
Standard seat weight	N/A (seat is integrated with base)
Std. seating dimensions	45.7 cm x 40.6 cm (18" x 16")
Warranty	2-year limited

COLOUR



Silver

Pride[®]
Mobility Products
Australia Pty. Ltd.

20-24 Apollo Drive Hallam, Victoria 3803
 ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • www.pridemobility.com.au



The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

FEATURES

- High-tensile carbon steel frame
- Folding back
- Height-adjustable armrests
- Adjustable back height
- Swing in/out leg rests
- Manual wheel locks
- 250 watt DC motor

Flash
 Shown with sling seat,
 swing in/out leg rests
 and PG VR2 Controller



FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tyre pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %.
- 3) Battery weight may vary based on manufacturer.



AU_ (R4) FLASH-SS-120117