

Footnotes

1. All shroud panels included.

2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition.
This specification can be subject to a variance of (+ or -) 10%.

3. Battery weight not included.

4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

5. 12 AH battery pack weighs 8.8 kg (19.5 lbs.).

6. Durable solid plastic.

7. Deep cycle (AGM or gel cell type recommended).

8. 18 AH battery pack weighs 13.4 kg (29.5 lbs.).

9. Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance.
Test conducted at maximum weight capacity.
10. Up to 14 miles per charge on 18 AH batteries.

11. Battery weight may vary based on manufacturer.

12. Width of base (not including seat).

13. Does not account for anti-tip clearance.

14. Includes standard seat weight.

15. Requires upgrade to medium or high-back contoured seat; additional charges apply.

16. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.

17. At midframe.

18. Range tested with 31-36 AH Batteries. Standard larger battery size would increase range.

PRIDE®

MOBILITY SCOOTERS

LIVE YOUR BEST®

END-TO-END



	Go-Go® Ultra X	Go-Go Elite Traveller® Plus	Go-Go® LX with CTS Suspension	Quest Folding Scooter	Jazzy® Zero Turn	Revo® 2.0	Pathrider® ES 10	Pathrider® 10	Victory 10® LX with CTS Suspension	Pathrider® 10 DX	Celebrity® DX	Celebrity® XL DX	Pathrider® 130 XL	Pathrider® 140 XL	Sport Rider
FRONT TYRES	5 cm x 20.3 cm (2" x 8") solid (3W) 5 cm x 17.8 cm (2" x 7") solid (4W)	5 cm x 20.3 cm (2" x 8") solid (3W)	5 cm x 17.8 cm (2" x 7") solid	4 cm x 17.8 cm (1.6" x 7") solid	5 cm x 17.8 cm (2" x 7") solid	7.6 cm x 20.3 cm (3" x 8") solid	10 cm x 25.4 cm (4" x 10")	9 cm x 25 cm (3.5" x 10")	9 cm x 25.4 cm (3.5" x 10")	9 cm x 25.4 cm (3.5" x 10")	8.89 cm x 25.40 cm (3.5" x 10")	10.16 cm x 33 cm (4" x 13")	10 cm x 33 cm (4" x 13")	10 cm x 33 cm (4" x 13")	7 cm x 35.5 cm (2.75" x 14")
REAR TYRES	6.35 cm x 20.3 cm (2.5" x 8") solid	7.6 cm x 22.9 cm (3" x 9") solid	6.35 cm x 20.3 cm (2.5" x 8") solid	5 cm x 20.3 cm (2" x 8") solid	6.35 cm x 22.9 cm (2.5" x 9") solid	7.6 cm x 22.9 cm (3" x 9") solid	10 cm x 25.4 cm (4" x 10")	9 cm x 25 cm (3.5" x 10")	9 cm x 25.4 cm (3.5" x 10")	10 cm x 25.4 cm (4" x 10")	10.16 cm x 25.40 cm (4" x 10")	10.16 cm x 33 cm (4" x 13")	10 cm x 33 cm (4" x 13")	10 cm x 33 cm (4" x 13")	7 cm x 35.5 cm (2.75" x 14")
MAXIMUM SPEED ²	Up to 6.4 km/h (4 mph)	Up to 7.2 km/h (4.5 mph)	Up to 7 km.h (4.4 mph)	Up to 6.4 km.h (4 mph)	Up to 6.4 km.h (4 mph)	Up to 8 km.h (5 mph)	Up to 8.85 km/h (5.5 mph)	Up to 8.4 km/h (5.25 mph)	Up to 9 km/h (5.6 mph)	Up to 8.85 km/h (5.5 mph)	Up to 9.65 km/h (6 mph)	Up to 15 km/h (9.3 mph)	Up to 15 km/h (9.32 mph)	Up to 15 km/h (9.32 mph)	Up to 15 km/h (9.32 mph)
GROUND CLEARANCE ^{4,13, 17}	7.6 cm (3") at midframe 3.2 cm (1.25") at motor	9.1 cm (3.6") at frame 4.3 cm (1.7") at motor	6.4 cm (2.5") at frame 3.8 cm (1.5") at motor	6 cm (2.3") at motor	8.9 cm (3.5") at frame 3.3 cm (1.3") at motor	8.9 cm (3.5") at frame 3.2 cm (1.25") at motor	7.0 cm (2.75") at mid-frame	7.6 cm (3") at mid-frame 5 cm (2") at motor	5.6 cm (2.2") at motor 7.4 cm (2.9") at frame	8.25 cm (3.25") at mdframe 6.4 cm (2.5") at motor	8.26 cm (3.25")	10.79 cm (4.25")	9.69 cm (3.81")	13.3 cm (5.25") at mid-frame	12.7 cm (5") at shroud
3-W TURNING RADIUS ⁴	83.8 cm (33")	N/A	N/A	N/A	N/A	101.6 cm (40")	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	158 cm (62.2")
4-W TURNING RADIUS ⁴	111.8 cm (44")	116.1 cm (45.72")	113.6 cm (44.75")	82 cm (32")	97.2 cm (38.25")	132.1 cm (52")	132.08 cm (52")	137 cm (54")	152 cm (59.9")	133.35 cm (52.5")	149.86 cm (59")	168.27 cm (66.25")	170 cm (67")	210.8 cm (83") Right 208.3 cm (82") Left	N/A
3-W LENGTH/WIDTH ^{4,12}	94 cm (37")/49.5 cm (19.5")	N/A	N/A	N/A	N/A	116.9 cm (46")/52.7 cm (20.75")	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
4-W LENGTH/WIDTH ^{4,12}	100.3 cm (39.5")/49.5 cm (19.5")	102.3 cm (40.25")/54 cm (21.25")	106 cm (41.75")/52.1 cm (20.5")	93 cm (37")/48.5 cm (19")	116.2 cm (45.75")/50.8 cm (20")	119.4 cm (47")/52.7 cm (20.75")	115.57 cm (45.5")/57.15 cm (22.5")	115 cm (45.25")/56.5 cm (22.25")	119.4 cm (47")/56.5 cm (22.2")	119.38 cm (47")/56.21 cm (22.13")	118.75 cm (46.75")/61 cm (24")	130.17 cm (51.25")/64.13 cm (25.25")	136 cm (53.54")/63.65 cm (25")	142.8 cm (56.25")/69.2 cm (27.25")	165 cm (65")/78 cm (30.7")
SEAT-TO-GROUND HEIGHT RANGE	46 - 51.13 cm (18.13"-20.13")	52 - 55.88 cm (20.5"-22")	52.7 - 56.5 cm (20.75"-22.25")	52 cm (20.5")	53.3 - 60 cm (21"-23.5")	53.3 - 60 cm (21"-23.5")	51.42 - 55.24 cm (20.25"-21.75")	52 - 55.8 cm (20.5"-22")	60 - 64.8 cm (23.5"-25.5")	60 - 63.8 cm (23.625"-25.125")	58.5 - 61.5 cm (23"-24.25")	53.97 - 64.13 cm (21.24"-25.24")	56.69 cm (22.31")	66.67 cm-70.48 cm (26.25"-27.75")	55.37 cm (21.8")
SEAT TYPE (STANDARD)	COMPACT, FOLDABLE 43.18 x 43.18 cm (17"x17")	COMPACT, FOLDABLE 43.18 x 43.18 cm (17"x17")	COMPACT, FOLDABLE 43.18 x 43.18 cm (17"x17")	COMPACT, FOLDABLE 41 x 33 cm (16"x13")	FOLDABLE MOLDED PLASTIC 43.18 x 43.18 cm (17"x17")	COMPACT, FOLDABLE 50.8 x 45.7 cm (20"x18")	FOLDABLE MOLDED PLASTIC 46 x 43 cm (18"x17")	FOLDABLE MOLDED PLASTIC 46 x 43 cm (18"x17")	RECLINING H/B W/SLIDERS & HEADREST 46 x 46 cm (18"x18")	CRS JUNIOR SEAT 44.25 x 43.18 cm (17.5"x17")	DELUXE MOLDED PLASTIC W/ SLIDERS 48.25 x 40.5 cm (19"x16")	CRS DELUXE RECLINING HIGH-BACK W/ SLIDERS & HEADREST 45.72 x 45.72 cm (18"x18")	PILLOW TOP DELUXE HIGH-BACK W/ SLIDERS 49.5 cm x 44.5 cm (19.5"x17.5")	PILLOW TOP DELUXE HIGH-BACK W/ SLIDERS 49.5 cm x 44.5 cm (19.5"x17.5")	CAPTAINS HIGH-BACK ADJUSTABLE SEAT 51 cm x 56 cm (20"x22")
(3/4-WHEEL) WEIGHT ^{3,4,14}	34.5 kg (76 lbs.)/36.3 kg (80 lbs.)	39.2 kg (86.5 lbs.)/39.5 kg (87 lbs.)	N/A/40.8 kg (90 lbs.)	N/A/23.5 kg (51.4 lbs.)	N/A/54 kg (120 lbs.)	55 kg (120 lbs.)/60 kg (132 lbs.)	53.29 kg (117.5 lbs.)/ 59.42 kg (131 lbs.)	N/A/62.4 kg (137.5 lbs.)	N/A/81.2 kg (179 lbs.)	N/A/84.82 kg (187 lbs.)	N/A/68.5 kg (151 lbs.)	N/A/90.71 kg (200 lbs.)	N/A/65.28 kg (143 lbs.)	N/A/117.70 kg (259.5 lbs.)	N/A/115 kg (253.5 lbs.)
SEAT-TO-DECK HEIGHT RANGE	40 cm - 43.8 cm (15.75"-17.25")	41.275 cm - 45 cm (16.25"-17.75")	41.3 cm - 43.8 cm (16.25"-17.25")	33 cm (13")	40.6 cm - 45.7 cm (16"-18")	38.1 cm - 45.7 cm (15"-18")	43.8 cm-47.62 cm (17.25"-18.75")	44.5 cm-48.3 cm (17.5"-19")	42 cm-47 cm (16.5"-18.5")	44.5 cm-48.3 cm (17.5"-19")	41 cm-44.5 cm (16.25"-17.5")	43.18 cm-53.34 cm (17"-21")	47 cm (18.5")	53.34 cm-57.15 cm (21"-22.5")	47.5 cm (18.7")
BATTERY WEIGHT ¹¹	4 kg (9 lbs.) ea. ⁵	6.35 kg (14 lbs.) ea. ⁸	4 kg (9 lbs.) ea. ⁵ 6.35 kg (14 lbs.) ea. ⁸	2.3 kg (5 lbs.) each	15 kg (32 lbs.) battery pack	11.5 kg (25 lbs.) each	11 kg (24.5 lbs.) each	11 kg (24.5 lbs.) ea./ 15 kg (32.5 lbs.) ea.	14.74 kg (32.5 lbs.) ea.	11 kg (24.5 lbs.) ea./ 14.74 kg (32.5 lbs.) ea.	11 kg (24.5 lbs.) ea./ 14.74 kg (32.5 lbs.) ea.	14.74 kg (32.5 lbs.) ea./ 15.1 kg (33.2 lbs.) ea.	17.16 kg (37.83 lbs.) ea.	(75 AH) 24 kg (53 lbs.) ea. (100AH) 34.24 kg (75.5 lbs.) ea.	(75 AH) 24 kg (53 lbs.) ea. (100AH) 34.24 kg (75.5 lbs.) ea.
BATTERY REQUIREMENTS ⁷	(2)12AH	(2)18AH	(2)12AH or (2)18AH	(1)11.5AH (Lithium-ion)	(2)22AH	(2) U-1 (30-34 AH)	(2) U-1 (30-34 AH)	(2) U-1 (30-34 AH) or (2) 40 AH	(2) 40 AH	(2) 35 AH or (2) 40 AH	(2) 32 AH or (2) 40 AH	(2) 40AH	(2) 55AH	(2) 75AH or (2) 100AH	(2) 75AH or (2) 100AH
PER CHARGE RANGE (MAX.) ^{2,9}	3W (12AH): Up to 11 km (6.9 miles) 4W (12AH): Up to 11.5 km (7.2 miles)	Up to 20.4 km (12.7 miles)	(12AH): Up to 9.5 km (5.9 miles) (18AH): Up to 15.1 km (9.4 miles)	Up to 20 km (12.5 miles)	Up to 16 km (10 miles)	3W: Up to 19 km (12 miles) 4W: Up to 21 km (13 miles)	Up to 25 km (15.5 miles)	Up to 21 km (13 miles)	Up to 26 km (16.16 miles)	Up to 33.38 km (21 miles) with 40AH batteries	Up to 20.63 km (12.8 miles)	Up to 40.23 km (25 miles) ¹⁸	Up to 24.99 km (15.53 miles)	70AH: Up to 35 km (21.75 miles)	100AH: Up to 50 km (31.25 miles)
BATTERY CHARGER	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2.5A	OFF-BOARD, 3.5A	OFF-BOARD, 3.5A	OFF-BOARD, 3A	OFF-BOARD, 3A(32AH) OFF-BOARD, 5A(40AH)	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 8A	OFF-BOARD, 8A	OFF-BOARD, 8A
SUSPENSION	N/A	N/A	Front & Rear	N/A	Front & Rear	Front & Rear	N/A	N/A	FRONT AND REAR	FRONT AND REAR	FRONT AND REAR	FRONT AND REAR	FRONT AND REAR	FRONT AND REAR	FRONT AND REAR
MAXIMUM WEIGHT CAPACITY	118 kg (260 lbs.)	136 kg (300 lbs.)	136 kg (300 lbs.)	115 kg (250 lbs.)	147 kg (325 lbs.)	170 kg (375 lbs.)	158 kg (350 lbs.) max	181.4 kg (400 lbs.) max	181.4 kg (400 lbs.) max	158.75 kg (350 lbs.) max	158.75 kg (350 lbs.) max	181 kg (400 lbs.) max	158.75 kg (350 lbs.) max	181 kg (400 lbs.) max	181 kg (400 lbs.) max
STANDARD BODY COLOURS	RED	RED, BLUE ¹	RED, BLUE ¹	WHITE	FIRE OPAL, GREEN MOONSTONE, ONYX BLACK, PEWTER, PEARL WHITE, SAPPHIRE BLUE	GREY STREET, TRUE BLUE	RED	VIPER BLUE	CANDY APPLE RED	VIPER BLUE	CANDY APPLE RED	CANDY APPLE RED OR SILVER	CANDY APPLE RED OR MERCURY SILVER	CANDY APPLE RED, YELLOW OR SILVER	BLACK
WEIGHT OF HEAVIEST PIECE ⁴	3W 13.4 kg (29.5 lbs.) (front section) 4W 15.4 kg (34 lbs.) (front section)	3W N/A 15.4 kg (34 lbs.) (front section)	3W N/A 4W 16.6 kg (36.6 lbs.) (front section)	3W N/A 4W 23.5 kg (51.4 lbs.) (unit does not disassemble)	3W N/A 4W 21.8 kg (48 lbs.) (rear section)	3W 21.8 kg (48 lbs.) (rear section) 4W 25.9 kg (57 lbs.) (front section)	3W N/A 25.8 kg (57 lbs.) 4W front section	3W N/A 27.7 kg (61 lbs.) 4W rear section	3W N/A 35.2 kg (77.6 lbs.) 4W rear section	3W N/A 29.48 kg (65 lbs.) 4W front section	3W N/A 32 kg (70.5 lbs.) 4W front section	3W N/A 37.64 kg (83 lbs.) 4W front section	3W N/A 34 kg (74.95 lbs.) 4W rear section	3W N/A 93.21 kg (205.5 lbs.) 4W main frame	92.3 kg (203.5 lbs.) 3W main frame 4W N/A
AVAILABLE IN 3-WHEEL FORMAT	YES	YES	NO	NO	NO	YES	NO	NO	NO	NO	NO	NO	NO	NO	YES
AVAILABLE IN 4-WHEEL FORMAT	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	YES	NO
TRAVEL MOBILITY						FULL-SIZE LUXURY						HEAVY-DUTY			

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or – 3 %.



AU_SC E2E_1-9-19



20-24 Apollo Drive Hallam, Victoria 3803
ACN 088 609 661 ABN 66 088 609 661

Tel: 03 8770 9600 • Fax: 03 9703 2720

pridemobility.com.au

