

Single-Fold Ramps from Pride® Lifts and Ramps offer a reliable solution for transitioning curbs, thresholds and other obstacles. They are available in several different lengths and can accommodate a weight capacity of 272 kg (600 lbs.). Plus, lightweight Single-Fold Ramps are easily portable, making them ideal for home and travel.

Single-Fold Ramps

FEATURES

- 272 kg (600 lbs.) weight capacity
- Lightweight & portable
- Integrated, high-traction surface
- Comprehensive product line
- Competitively priced
- 75 cm (29.5") wide
- Available in: 61 cm, 91 cm, 122 cm, 152 cm and 183 cm (2', 3', 4', 5' and 6')

Ramp Length	Rise				
	2' (61 cm)	3' (91 cm)	4' (122 cm)	5' (152 cm)	6' (183 cm)
17" (30%)	11.2" (20%)	8.4" (15%)	6.7" (12%)	5.6" (10%)	4.8" (8%)
14.5" (26%)	9.6" (17%)	7.2" (13%)	5.7" (10%)	4.8" (8%)	4" (7%)
12" (21%)	8" (14%)	6" (10%)	4.8" (8%)	4" (7%)	3.2" (6%)
9.6" (17%)	6.4" (11%)	4.8" (8%)	3.8" (7%)	3.2" (6%)	2.4" (4%)
7.2" (13%)	4.8" (8%)	3.6" (6%)	2.9" (5%)	2.4" (4%)	-
4.8" (8%)	3.2" (6%)	2.4" (4%)	1.9" (3%)	-	-
-	-	-	-	-	-

