

Pride[®] Lift Chairs

Experience the Luxury

Power Stand Assist Chair



A comfortable chair that's stylish and durable

The Power Stand Assist is a seat lift designed for use throughout your entire home. The chair can easily assist with standing up from the dining room table or transferring into a wheelchair in the living room. With its stylish wood design and comfortable fabric cushions, the Power Stand Assist is the perfect addition to any décor.

SPECIFICATIONS

LC-440 Power Stand Assist

Back Width:	53.5 cm (21")
Seat-To-Floor:	49 cm (19")
Seat Depth:	46 cm (18")
Seat Width:	54.5 cm (21.5")
Top of Back to Seat:	69 cm (27")
Weight Capacity:	125 kg (275 lbs.)

FABRIC SELECTION (MICRO-SUEDE)



Cocoa



Fern



Sandal

Note: Measurements may vary +/- 2.5 cm (1").
Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Experience the Luxury