

**Pride® Lift Chairs**

*Experience the Luxury*

# T-3 Lift Chair

## An infinite-position lift chair that delivers comfort and durability

Engineered for comfort and durability, the T-3 Infinite-Position Lift Chair from Pride® reclines to a full Trendelenburg position that is designed to elevate the legs to a point above the head. The stylish T-3 is available in Black Leather or Oatmeal fabric to enhance your home's décor.

### SPECIFICATIONS

#### T-3

**Back Width:** 69 cm (27")

**Seat-to-Floor:** 48 cm (19")

**Seat Depth:** 51 cm (20")


**Seat Width:** 51 cm (20")

**Top of Back to Seat:** 68.5 cm (27")

**Weight Capacity:** 170 kg (375 lbs.)

**Warranty:** 2-year limited

### FABRIC SELECTION



**Black  
Leather**



**Oatmeal**

Note: Measurements may vary +/- 2.5 cm (1").  
Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

**Pride®**  
Mobility Products  
Australia Pty. Ltd.

**PRIDE®**  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST® —

20-24 Apollo Drive  
Hallam, Victoria 3803  
[www.pridemobility.com.au](http://www.pridemobility.com.au)

Tel: 03 8770 9600  
Fax: 03 9703 2720

ACN 088 609 661 ABN 66 088 609 661

*Experience the Luxury*

