

TRU COMFORT™

The new Synergy® TRU-Comfort™ seating system provides the benefits of adjustable positioning and pressure management while maintaining a sporty and stylish look. TRU-Comfort's innovative lumbar and lateral bolsters easily adjust to deliver personalized comfort.

TRU-Comfort offers positioning through the medial thigh support of 31.8mm in height and the lateral thigh support of 49.8mm in height. This support comes from the cushion base design with these features not appearing as dramatic due to the properties of the upper layer of viscoelastic foam which is pulled down at the peaks by the cover. TRU-Comfort also keeps the occupant from slouching forward due to the lower back to front to base design.

FEATURES

- Ultra-soft, leather-like cover
- High density HR polyurethane base for long lasting comfort
- Sensus viscoelastic foam top layer to manage pressure, promote proper circulation and provide comfort
- Easily added to TRU-Balance® and other Synergy® seats
- Adjustable lumbar support and lateral bolsters to dial in positioning comfort
- Modular seating components allow use of alternate seating cushions
- T-slot channel for easy swing-away lateral attachment and adjustment

TRU-Comfort™
seating system
shown on
Quantum® 600



"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".

Adjustable lumbar
support and lateral
bolsters



T-slot channel for
swing-away laterals



Electronics are
concealed to create
a non-medical
appearance



Removable
back



Pride
Mobility Products
Australia Pty. Ltd.



21 Healey Rd.
Dandenong, Victoria 3175
www.pridemobility.com

ACN 088 609 661 ABN 66 088 609 661

Tel: (03) 9706 4611
Fax: (03) 9706 4622

QUANTUM
REHAB
INNOVATIVE REHAB SOLUTIONS