

*Additional Seating Options Available from Quantum Rehab®:*



**Synergy® Seating System**

- Easily adjustable and versatile design that can grow with the client
- Easy customisation
- Several width and depth adjustment packages ranging from 25 cm (10") to 61 cm (24")
- Accepts a wide variety of after-market seats and backs
- *Heavy-Duty* option provides a weight capacity of 295 kg (650 lbs.) and seating up to 81 cm x 71 cm (32" x 28") (model specific)

**Synergy® Manual Tilt**

- 34 kg (75 lbs.) weight capacity
- Adjustable widths of 25 cm - 30.5 cm (10" - 12") and 30.5 cm - 35.5 cm (12" - 14")
- Adjustable depths of 25 cm - 35.5 cm (10" - 14")
- Mercury drive lockout
- Optional paediatric stroller-type push handle



**Synergy® Manual Recline**

- Solid curved back pan and 6.35 cm (2 1/2") raised pivot point to reduce shear forces
- Attendant activated recline seat
- 25 cm - 61 cm (10" - 24") wide / 25 cm - 61 cm (10" - 24") deep
- 136 kg (300 lbs.) weight capacity
- 180° of recline




**TRU BALANCE®**  
POWER POSITIONING SYSTEMS



**SYNERGY®**  
cushions and backs

- An excellent selection of pressure management and positioning solutions
- Cushions available as small as 25 cm x 25 cm (10" x 10") with sizes up to 81 cm x 71 cm (32" x 28") (295 kg (650 lbs.) weight capacity) to accommodate clients of all sizes
- The industry's shortest lead times for custom sizes ensure your clients get the exact protection they need quickly and easily
- Unique design and exclusive materials provide excellent support and protection



Spectrum® Air Contoured Cushion

*For more information, please refer to our Synergy® Cushions and Backs brochure.*

**QUANTUM®**  
LIFE BEYOND LIMITS™

20-24 Apollo Drive Hallam, Victoria 3803  
ACN 088 609 661 ABN 66 088 609 661  
(Tel) 03 8770 9600 • (Fax) 03 9703 2720  
[www.pridemobility.com.au](http://www.pridemobility.com.au)



The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

*Advanced Choices for Complete Multiple Power Solutions*

### TRU-Balance® Power Tilt

- 55° of tilt
- Lower standard seat-to-floor height
- New tilt interface for fast and easy adjustments
- Maintains center of gravity shift



### TRU-Balance® Lift & Tilt

- Delivers 25.5 cm (10") of elevation with 50° of tilt
- Four-point mount lift for maximum stability, low seat-to-floor heights and the use of maximum capacity batteries
- Maximum versatility and ease of adjustment
- Available with power recline, power front riggings and TRU-Comfort seating
- 136 kg (300 lbs.) weight capacity



### TRU-Balance® Power Recline

- 95-168° of recline
- 136 kg (300 lbs.) weight capacity
- Does not increase seat-to-floor height
- Standard shear reduction and reclining armrests
- Available with TRU-Balance Power Tilt and TRU-Balance Lift & Tilt



### TRU-Balance® Power Elevating Seat

- Offers 18 cm (7") of lift with power elevating seat
- Allows access to batteries when seat is elevated
- 136 kg (300 lbs.) weight capacity
- Allows drive capability at ¼ speed when elevated



### TRU-Balance® Paediatric Tilt

- Delivers maximum adjustability and versatility for pediatric applications
- Offers low seat-to-floor heights along with easy width and depth adjustment to accommodate growth
- Available in sizes as small as 25 cm x 25 cm (10" x 10")



### TRU-Balance® HD Power Tilt

- 40 degrees of pressure-relieving tilt
- Available on the Q6000Z HD and the Q600 XL (HD)
- Weight capacities up to 204 kg (450 lbs.) on the Q6000Z HD and 181 kg (400 lbs.) on the Q600 XL (HD)
- Maximum seat sizes of 71 cm x 71 cm (28" x 28") on the Q6000Z HD and 61 cm x 61 cm (24" x 24") on the Q600 XL (HD)



# TRUBALANCE®

POWER POSITIONING SYSTEMS

### TRU-Balance® Power Articulating Foot Platform

- 20 cm (8") of articulation
- 90-170° of elevation
- 6-way adjustable calf pads
- Bi-lateral, angle-adjustable split foot platform



### TRU-Balance® Power Articulating Leg Rests

- 70° or 80° angle
- Actuator is mounted to the outside of leg rest and will not increase the seat-to-floor height of the seating system
- Connected with a ¼" headphone-style jack connector (pull out in order to release)
- Depth and height-adjustable calf pads
- 10 cm (4") or 15 cm (6") articulation
- Drop-in or swing-away



### Q-Logic Expandable Attendant Control

- Full proportional expandable joystick
- Standard on multiple power configurations equipped with Q-Logic Drive Controls
- Allows caregivers the ability to provide patient assistance



### TRU-Comfort

- Ultrasoft, leather-like cover
- High-density HR polyurethane base for long lasting comfort
- Viscoelastic foam top layer to manage pressure, promote proper circulation and provide comfort
- Adjustable lumbar support and lateral bolsters to dial in positioning comfort



### TRU-Comfort Plus

- All of the benefits of TRU-Comfort Seating
- Four-way stretch, fluid resistant cover
- Additional 1¼ cm (½") of viscoelastic foam in seat cushion for maximum comfort and pressure management



### Recline Back Interface for After-Market Backs

- Optional integration system allows for an unprecedented level of compatibility between TRU-Balance® Power Recline and most after-market back rests
- Retains full functionality, including shear reduction



### Vent Tray Solutions

- Numerous vent tray and oxygen holder configurations are available like our compact vent tray (for LTV style ventilators) and articulating vent tray



Compact Vent Tray

**Please Note:** Some options may not be available on all power bases. Please contact Quantum® Customer Service for availability.

Articulating Vent Tray