

# TRU COMFORT

## SEATING AND POSITIONING

**TRU-Comfort Seating and Positioning** products provide adjustable positioning and pressure management and are compatible with both TRU-Balance® Power Positioning Systems and the Synergy® Seat Frame System. TRU-Comfort cushions offer a deep contoured, high-density moulded foam seat cushion with a viscoelastic total top layer, an ergonomically contoured back cushion and a highly adjustable headrest. **TRU-Comfort 2 Cushions** have a deep contoured, high-density moulded foam seat cushion with a viscoelastic total top layer, which makes for superb comfort and support.




### TRU-COMFORT HEADREST

- Available in 20 cm (8") and 25 cm (10") widths  
*(20 cm (8") shown)*
- Choice of TRU-Comfort and TRU-Comfort Plus covers
- Streamlined mounting bracket offers a greater range of adjustment



### TRU-COMFORT PLUS

- All of the benefits of TRU-Comfort Seating
- Four-way stretch, fluid resistant cover
- Additional 1.25 cm (1/2") of viscoelastic foam in seat cushion for maximum comfort and pressure management



### TRU-COMFORT 2

- Four-way stretch, recovery, fluid-resistant, breathable cover
- Additional 2.5 cm (1") of viscoelastic foam in seat cushion for maximum comfort and pressure management
- Pelvic well
- High contour side for femur positioning
- Reduces pressure (ischial tuberosities)



20-24 Apollo Drive Hallam, Victoria 3803  
ACN 088 609 661 ABN 66 088 609 661  
(Tel) 03 8770 9600 • (Fax) 03 9703 2720  
[www.pridemobility.com.au](http://www.pridemobility.com.au)

The information contained herein is correct at the time of publication.  
We reserve the right to alter specifications without prior notice.



QUANTUM REHAB®