

A SMOOTH RIDE for your Active Lifestyle



VICTORY® 10
LX WITH CTS SUSPENSION

CTS SUSPENSION



Front suspension



Rear suspension


The Victory® 10 LX with CTS Suspension offers a smooth ride with performance to match. Comfort-Trac Suspension (CTS) is the next generation in mobility scooter ride comfort. This industry first technology incorporates advanced suspension components for improved absorption over rough surfaces and a more comfortable ride experience over various terrains.

PRIDE® **MOBILITY** SCOOTERS
LIVE YOUR BEST®

| pridemobility.com.au



Specifications 4W:S710LX

Weight Capacity	181 kg (400 lbs.)
Maximum Speed ¹	Up to 9 km/h (5.6 mph)
Ground Clearance ³	5.6 cm (2.2") at motor
Turning Radius ³	165.1 cm (64.9")
Overall Length ³	119.4 cm (47")
Overall Width ³	56.5 cm (22.2")
Seat-to-Ground Height	42 - 47 cm (16.5"-18.5")
Seat-to-Deck Height	59.7 - 57.2 cm (23.5"-22.5")
Front Tyres	9 x 25.4 cm (3.5" x 10") solid
Rear Tyres	9 x 25.4 cm (3.5" x 10") solid
Suspension	Front and rear
Range Per Charge ^{1,2}	Up to 26 km (16.16 miles)
Total Weight Without Batteries ^{3,7}	70.9 kg (156.2 lbs.)
Heaviest Piece When Disassembled ³	35.2 kg (77.6 lbs.) (front section)
Standard Seating	Type: Deluxe reclining high-back Weight: 19 kg (42 lbs.) Material: Black Vinyl Dimensions: 45.7 cm (18") width, 45.7 cm (18") depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Electronic, regenerative and electromechanical
Battery Requirements ^{3,4,5}	Size: (2) 40 AH (12V deep-cycle) Weight: 14.5 kg (32 lbs.) each
Battery Charger	5-amp, off-board (40 AH)
Colours	 Candy Apple Red

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).
- ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).
- ³ AGM or gel cell type required.
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

Features

- Comfort-Trac front and rear suspension
- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Easy drive tiller with wraparound handles
- Backlit battery gauge
- Exclusive stylish, lightweight, non scuffing, black, low-profile tyres
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination with a courtesy light under tiller
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrest
- Storage compartment below tiller
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port
- Auto-connecting front-to-rear harness

Accessories

- Pride[®] rear bag (includes crutch/walking stick holder)
- Walker holder (must have Pride[®] rear bag to attach)
- Oxygen tank holder (must have Pride[®] rear bag to attach)
- Saddlebag
- Mobile phone holder
- Cup holder
- Flag
- Rearview mirror