

Pride® Walkers

PMW460



FEATURES

- Lightweight aluminium frame
- Detachable padded backrest tube
- 5" front wheels and 5" rear twin wheels
- Bolt head-inner tube knob interface design tightens the handle tubes securely
- Push down brake
- Vinyl pouch standard

PMW460

4-Wheeled, Aluminium, Push Down

TECHNICAL SPECIFICATIONS

Seat Height:	55 cm (21.5")
Width:	63.5 cm (25")
Height:	85-95 cm (33.5"-37.5")
Length:	57 cm (22.5")
Seat Width:	35 cm (14")
Seat Depth:	31 cm (12")
Weight:	6.1 kg (13.5 lb)
Weight Capacity:	114 kg (250 lb)
Colour:	Red

PMW462

4-Wheeled, Aluminium, Loop Brake, 6" Wheel

TECHNICAL SPECIFICATIONS

Seat height:	55 cm (21.5")
Width:	63.5 cm (25")
Height:	85-95 cm (33.5"-37.5")
Length:	57 cm (22.5")
Seat Width:	35 cm (14")
Seat Depth:	31 cm (12")
Weight:	6.5 kg (14.5 lb)
Weight Capacity:	114 kg (250 lb)
Colour:	Red

FEATURES

- Lightweight aluminium frame
- Detachable padded backrest tube
- Four 6" x 1" wheels
- Bolt head-inner tube knob interface design tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard

PMW462



PMW463



FEATURES

- Lightweight aluminium frame
- Detachable padded backrest tube
- Four 200 x 38 PU wheels
- Bolt head-inner tube knob interface design tightens the handle tubes securely
- Loop brake is ideal for slow down and permanent stop
- Vinyl pouch standard

PMW463

4-Wheeled, Aluminium, Loop Brake, 8" Wheel

TECHNICAL SPECIFICATIONS

Seat height:	55 cm (21.5")
Width:	63.5 cm (25")
Height:	88-98 cm (33.5"-38.5")
Length:	60 cm (23.5")
Seat width:	35 cm (14")
Seat depth:	31 cm (12")
Weight:	7 kg (15.5 lb)
Weight capacity:	114 kg (250 lb)
Colour:	Red

21 Healey Rd.
Dandenong, Victoria 3175
www.pridemobility.com

Tel: 03 9706 4611
Fax: 03 9706 4622
www.pridemobility.com

Pride
Mobility Products
Australia Pty. Ltd.