

ITEM#

A **L** **B** **R** **1** **6** **S** **A** **Q**

SHELL TYPE	BACK REST	CHAIR WIDTH	HEIGHT	CONTOUR	MOUNTING
AL-AL Series (Scapular Cut Out)	NOTE: SHELL ASSEMBLY WEIGHT AL & AFT SERIES - 1.5lb/680gm - above weights are of a 16" STD - all weights approx.	16"/40.5cm	S Standard (16"/40.5cm)	A Active (2.5"/63.5mm)	Q Q-Back (Re-configured with TB3 compat- ible hole pattern and mounting hardware)
AFT-AFT Series (Flat Top)		18"/45.5cm	T Tall (18"/45.7cm)	D Deep (4"/100mm)	
		20"/51cm	X x-tall (20"/51cm)	6 x-deep (6"/150mm)	
			XX xx-tall (22"/56cm)		



Contact us by e-mail:
info@stealthproducts.com
 or by phone:
 +1 (512) 715-9995
 +1 (800) 965-9229

