



## LUXURY FEATURES AT AN EXCEPTIONAL VALUE

The **Celebrity® Series** of scooters are feature-rich luxury products at an attractive price. Available in a 3 or 4-wheel format, the Celebrity LR offers elegant simplicity, unmatched reliability, impressive performance and affordable pricing. The Celebrity DX Sport 4-wheel is a perfect blend of style, performance technology, convenient standard features and easy-to-use operation, making any of the Celebrity Series of scooters a great overall choice.

### Celebrity® LR 3 & 4-Wheel Scooter



### Celebrity® DX Sport 4-Wheel Scooter



**Pride**  
Mobility Products Co.  
LIVE YOUR BEST®

PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

Beamsville, Ontario • 888-570-1113  
[pridemobility.com](http://pridemobility.com)

# SPECIFICATIONS

	<b>CELEBRITY® LR (3 &amp; 4-wheel)</b>	<b>CELEBRITY® DX SPORT</b>
Max. weight capacity:	350 lbs.	350 lbs.
Maximum speed: <sup>1</sup>	Up to 6.5 mph	Up to 9.5 mph
Ground clearance: <sup>2</sup>	4"	4.25"
Turning radius: <sup>2</sup>	3W: 44.5"/ 4W: 56.5"	66.25"
Overall length: <sup>2</sup>	3W: 47.5"/ 4W: 48"	47.5"
Overall width: <sup>2</sup>	24"	25.75"
Battery charger:	Off-board, 5 amp	Off-board, 5 amp
Per-charge range: <sup>1,5</sup>	3W: Up to 17.99 miles 4W: Up to 18.97 miles	Up to 10.81 miles
Battery requirements: <sup>3</sup>	(2) 12-volt, deep cycle Size: NF-22	(2) 12-volt, deep cycle Size: 40 AH
Battery pack weight: <sup>4</sup>	NF-22: 36.5 lbs. each	40 AH: 32.5 lbs. each
Weight w/o battery pack: <sup>2</sup>	3W: 160 lbs. / 4W: 173 lbs	189 lbs.
Heaviest piece: <sup>2</sup>	3W: 61 lbs. (front section) 4W: 73 lbs. (front section)	63 lbs. (front section)

## FEATURES

- Deluxe, dual-adjustable delta tiller with wraparound handles
- Ultra-simple, auto-latching lockup for quick disassembly/reassembly
- Deluxe light package with turn signals
- Speed adjustable up to 6.5 mph
- 350 lbs. weight capacity
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front frame-mounted seat post for optimal stability
- Jumbo front basket and rearview mirror standard
- Larger battery capacity for optimal range

- State-of-the-art rear suspension
- Larger rear wheels
- Top speed of 9.5 mph
- Maximum weight capacity of 350 lbs.
- Ultra-simple, auto-latching lockup for quick disassembly/reassembly
- Easily accessible tie down points (for transport of unoccupied scooter)
- Includes standard angle-adjustable dual bulb curb light
- Front basket and rearview mirror standard

### COLOR CHOICES (FOR BOTH)



Candy Apple Red



Onyx Black

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup>Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.

<sup>3</sup>AGM or gel cell type recommended.

<sup>4</sup>Battery weight may vary based on manufacturer.

<sup>5</sup>Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

CA\_Celeb Series SS\_9-1-15\_REV A



**PRIDE® MOBILITY SCOOTERS**  
LIVE YOUR BEST®

Beamsville, Ontario • 888-570-1113  
pridemobility.com