


DYNO™



Facebook.com/
LeggeroUSA



- Folds Compactly ◇
- 3 Wheel Configuration ◇
- Swiveling Front Wheel ◇
- ACTIVATOR™ Dynamic Seating ◇
- 30° Hip Angle Adjustment ◇
- Lightweight ◇
- Contoured Seating System ◇



Distributed by Pride Mobility Products Company, Beamsville, ON, Canada
For customer enquiries, email Productorders@pridemobility.com or call us at 888-570-1113.
Web: Leggero.us • facebook.com/LeggeroUSA • © 2019 LEGGERO LLC CADYNOB-1909

DYNO™



The only WC-19
transit tested
compact folding
3-wheeler!

Available with
Activator™
Dynamic Seating



SPECIFICATIONS & ADJUSTMENT RANGE

Weight Limit	80 lbs / 36 kg
Weight of Chair	26.5 lbs / 12 kg
Chair Width	26.5 in / 67 cm
Seat Depth	6-14 in / 15-35 cm
Seat Width	14 in / 35 cm
Back Height	13-21 in / 33-53 cm
Lower Leg Length	8-14 in / 20-35 cm
Chest Width Range	7- 14 in / 18-35 cm
Seat to Floor	19 in / 48 cm
Hip Angle Adjustment	85°, 90°, 105°
Adjustable Seat Tilt	25-30°
Knee Active Range (optional)	10°
Seat to Back Active Range (optional)	30°
Push Handle Height	39-46 in / 99-117 cm
Front Wheel Diameter	10 in / 25 cm
Rear Wheel Diameter	16 in / 41 cm
Folded with Wheels	42 x 17 x 26.5 in 107 x 43 x 67 cm
Folded without Wheels	35 x 17 x 25 in 89 x 43 x 64 cm

Freedom to go where you want and the
positioning options you need.



For more information visit www.Leggero.us