



USB  
Charger  
Port

LC-570M shown in  
Black Cherry fabric

# ELEGANCE COLLECTION

Pride® Elegance Collection Lift Chairs are designed to be the ultimate in lift chair comfort, style and performance. Available in several beautiful fabrics and styles to fit any décor, Pride Elegance Lift Chairs also feature a 7-year prorated warranty and head & arm covers. The state-of-the-art lift system utilized in Pride Elegance Lift Chairs will recline you back, providing wonderful relaxation.

NORTH AMERICA'S #1  
POWER LIFT RECLINER

POWER  
RECLINERS

Protect your  
investment with  
the Pride® Lift  
Chair Care Kit



Beamsville, Ontario | 888-570-1113  
pridemobility.ca | prideliftchairs.com

**Pride**  
Mobility Products Co.  
LIVE YOUR BEST®

# ELEGANCE COLLECTION

NORTH AMERICA'S #1  
POWER LIFT RECLINER

POWER  
LIFT  
RECLINERS



SPECIFICATIONS	Medium	Small	Medium	Large
	LC-450C	LC-570S	LC-570M	LC-570L
	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'9" - 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	20"	19"	20"	21"
Seat Depth	20"	19"	20"	22"
Seat Width	20"	19"	20"	22"
Top of Back to Seat	27"	24.5"	26.5"	27"
Back Style	Blown Fiber Split Back <sup>4</sup>	Blown Fiber Pillow Back	Blown Fiber Waterfall Back	Blown Fiber Waterfall Back
Fully Padded Chaise	Yes <sup>5</sup>	Yes	Yes	Yes
Heat and Massage	Yes	Yes	Yes	Yes
Seating Type	Total Comfort <sup>6</sup>	Total Comfort <sup>6</sup>	Total Comfort <sup>6</sup>	Total Comfort <sup>6</sup>
Distance From Wall	18"	16.5"	19"	24"
Overall Chair Width	37"	35"	36"	38"
Chair Weight	119 lbs.	110 lbs.	113 lbs.	114 lbs.
Warranty	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*



LC-570M shown in  
Cloud 9 Stone fabric

## ELEGANCE FABRIC SELECTION

### CLOUD 9



BLACK CHERRY



PACIFIC



STONE



WALNUT

## BENEFITS

- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation
- Premium furniture grade laminate/hardwood frames for added strength and durability
- Easy-to-use hand control with USB charger port (*not available when adding optional heat and massage*)
- Standard head and arm covers
- Large, dual pockets for extra storage
- Quiet and smooth lift system
- Integrated battery backup ensures chair is operational when power goes out<sup>1</sup>

## OPTIONS\*

- Deluxe heat and massage
- Residential cleaning kit
- Head pillow
- Footrest extension
- Wireless remote (LC-450C)
- Mattress option (LC-450C only)

Additional lead time may  
apply for some options.

\*Additional charges apply for all options.

1. Batteries not included.
2. Available with deluxe heat and massage only. Standard heat and massage not available.
3. Please check for available inventory.
4. Back upholstery is fully removable and adjustable.
5. Seat pad is fully removable and adjustable.
6. Unique spring design and a one-piece soft foam seat for outstanding relaxation.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.  
Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.