

# GO-GO

ELITE TRAVELLER®

## The Quickest Way to Travel.

The Go-Go Elite Traveller®'s compact design allows it to easily maneuver in tight spaces while providing stable outdoor performance. Take the guesswork out of travel with simple, one hand feather-touch disassembly. Go where you want to go, easily, with the Go-Go Elite Traveller.



Go-Go Elite Traveller® 3-wheel



Go-Go Elite Traveller®  
4-wheel

Available as a 3 & 4-wheel scooter



Feather-touch disassembly.

**Pride**  
Mobility Products Co.  
— LIVE YOUR BEST —



## SPECIFICATIONS

Model numbers	SC40E (3-wheel) SC44E (4-wheel)
Max. weight capacity	275 lbs.
Maximum speed <sup>2</sup>	Up to 4.25 mph
Ground clearance <sup>1</sup>	2.5"
Turning radius <sup>1</sup>	3-wheel: 32.5" 4-wheel: 44.125"
Overall length <sup>1</sup>	3-wheel: 37" 4-wheel: 39.5"
Overall width <sup>1</sup>	19.5"
<b>TIRES:</b>	
Front	3-W: 2" x 8" solid 4-W: 2" x 7" solid
Rear	2.5" x 8" solid
Anti-tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded plastic

Drivetrain	24-volt DC motor, sealed mini-transaxle
------------	---

Dual braking system	Regenerative and electromechanical
---------------------	------------------------------------

Battery charger	Off-board, 2 amp
-----------------	------------------

Per-charge range <sup>2,5</sup>	Up to 14 miles
---------------------------------	----------------

Battery requirements <sup>3</sup>	(2) 12V deep-cycle (standard)
-----------------------------------	-------------------------------

Compatible batteries	17 AH (standard)
----------------------	------------------

## WEIGHT:

Batteries (2 required) <sup>4</sup>	33 lbs. battery pack
Weight w/o battery pack	3-wheel: 75 lbs. 4-wheel: 78.5 lbs.
Heaviest piece	3-wheel: 28 lbs. (front section) 4-wheel: 32.5 lbs. (front section)

## Colors



Red



Blue

2 sets of exciting colors with every scooter!


## FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Exclusive black, non-scuffing tires
- Auto-connecting front to rear cable
- Includes 2 sets of easily changeable colored panels in Red and Blue
- Front frame-mounted seat post offers maximum stability
- All-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Standard front tiller-mounted basket included
- Removable deck basket for additional storage (4-wheel only)
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board


## Pride® Lifts & Ramps

These Pride™ Lifts are compatible with the Go-Go Elite Traveller®.

- Backpacker® Plus
- Backpacker® MV
- Power Tote Mini (3-wheel only)
- SilverBoom 250



SilverBoom 250



Outlander

Beamsville, Ontario • 888-570-1113 • [www.pridemobility.com](http://www.pridemobility.com)

<sup>1</sup>Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.

<sup>2</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>3</sup>AGM or gel cell type recommended.

<sup>4</sup>Battery weight may vary based on manufacturer.

<sup>5</sup>Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

