



Shown in Cloud 9
Black Cherry

Grandeur

NORTH AMERICA'S #1
POWER LIFT RECLINER

Pride
POWER
LIFT
RECLINERS

Experience luxury at its finest with the **Grandeur Collection** from Pride® Power Lift Recliners. True infinite-position seating and pocketed coils deliver lavish comfort and maximum performance. With three different fabric selections to choose from, the **Grandeur Collection** offers superior relaxation with contemporary style and elegance.

Beamsville, Ontario | 888-570-1113
pridemobility.com | prideliftchairs.com

Pride
Mobility Products Co.
— LIVE YOUR BEST —

Grandeur

NORTH AMERICA'S #1
POWER LIFT RECLINER

Pride
POWER
LIFT
RECLINERS



LC-515iL SPECIFICATIONS: LARGE

Position	True Infinite Trendelenburg Zero Gravity 5'4" - 6'
Weight Capacity	375 lbs.
Height	43"
Seat to Floor	20"
Seat Depth	20"
Seat Width	21"
Top of Back to Seat	25"
Back Style	Blown fiber, Pillow back
Fully Padded Chaise	Yes
Footrest Extension	Standard
Heat and Massage	Yes
Seat Type	Pocketed coil springs
Distance from the Wall	21.5"
Overall Chair Width	33.5"
Chair Weight	131 lbs.
Warranty	7-Year Pro-rated

STANDARD:

CLOUD 9




OPTIONAL:

ULTRALEATHER™



ARIA² CRYPTON[®]

super fabric



FEATURES AND BENEFITS

- Easy to remove back
- Industry-first hand control with USB charger port (not available with heat and massage)
- Integrated battery backup in case of power failure¹
- Large, dual pockets
- Pocketed coils deliver individualized support and relaxation
- Premium furniture grade laminate/hardwood frames for added strength and durability
- Standard head and arm covers
- Standard lumbar pillow
- Transformer features self diagnostic electronics with 9' cord

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

QS Look for the Quick Ship Indicator for the shortest lead times available.

¹ Batteries not included.
² Crypton[®] Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in fabric.