



LC-358M shown in
Cloud 9 Stone fabric

HERITAGE COLLECTION

NORTH AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS

Comfort and style complement one another perfectly in the **Heritage Collection** from Pride® Power Lift Recliners. With an array of plush fabrics to fit any décor and a superior design that provides exceptional strength, Heritage Collection delivers style and durability that will last for years to come. Pride's Heritage Collection offers numerous standard features, superb comfort and quiet and smooth performance, all at a great value.

Beamsville, Ontario | 888-570-1113
pridemobility.com | prideliftchairs.com

Pride
Mobility Products Co.
LIVE YOUR BEST®



SPECIFICATIONS	Petite Wide	Small	Medium	Large	Heavy-Duty	Very Heavy-Duty
	LC-358PW ^{2,5}	LC-358S ⁵	LC-358M ^{2,5}	LC-358L ^{2,5}	LC-358XL ^{4,6}	LC-358XXL ^{4,6,7}
	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4" - 6'	3-Position Chaise Lounger 5'8" - 6'2"	3-Position Chaise Lounger 5'6" - 6'1"	2-Position Chaise Lounger 5'6" - 6'1"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	19"	19"	20"	21"	20"	20"
Seat Depth	19"	19"	20"	22"	22"	22"
Seat Width	24"	18"	20"	22"	26"	30"
Top of Back to Seat	23.5"	24"	24.5"	25.5"	29"	29"
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage	Optional	Optional	Optional	Optional	Optional	NA
Seating Type	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Total Comfort Seating	Sinusuous Spring	Sinusuous Spring
Distance From Wall	18"	18"	17"	19"	19"	14"
Overall Chair Width	40"	35"	35"	37"	39.5"	43"
Chair Weight	130 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.
Warranty	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*	7-Year Prorated	7-Year Prorated

HERITAGE FABRIC SELECTION

CLOUD 9



BLACK CHERRY PACIFIC STONE WALNUT

*358XXL - Stone and Walnut only available

LEXIS STA-KLEEN®



BLACK MUSHROOM BURGUNDY CHESTNUT

*358XXL - Chestnut and Mushroom only available

ULTRALEATHER™



BUFF FUDGE GARNET PECAN CHARCOAL

ARIA³



ESPRESSO COOL GREY RED LAZULI SAND

FEATURES AND BENEFITS

- Standard footrest extension on 358PW, M and L for additional leg support while reclined
- Easily adjustable leg levelers provides additional support for uneven surfaces
- Easy-to-use hand control with USB charger port (not available when adding optional heat and massage)
- Standard head and arm covers
- Large dual pockets for extra storage
- Industry-first Total Comfort Seating features a unique spring design and one-piece soft foam seat for outstanding relaxation (except LC-358XL and LC-358XXL)
- Quiet and smooth lift and scissor mechanism system
- Integrated battery backup ensures chair is operational when power goes out
- Energy-efficient power box equipped with self-testing electronics and 9' cord
- Premium furniture grade laminate/hardwood frame for added strength and durability

OPTIONS*

- Heat and massage (except XXL)
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Wireless remote (not available on XL or XXL)
- Non-skid leg levelers

Additional lead time may apply for some options.

*Additional charges apply for all options.

1. Batteries not included.

2. Standard footrest extension.

3. Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.

4. Removable, fully padded chaise.

5. Unique spring design and a one-piece soft foam seat for outstanding relaxation.

6. Features Dual Motor Sync System.

7. Only available in Stone, Walnut, Lexis Sta-Kleen, Chestnut and Mushroom.

8. Standard with footrest extension and lumbar pillow.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

*Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



Look for the Quick Ship Indicator for the shortest lead times available.