

# Color Options



Edge 3

Blue Jean	Cherry Bomb	Lemon Crush	Green Machine	Raspberry Beret	Evolution Orange
Electric Blue	Deep Purple	Sugar Plum	Red Corvette	White Satin	Grey Street
Midnight Blue	Back in Black	Night Train	Military Green	Carbon Fiber	Digital Camo



Q6 Edge<sup>®</sup> 2.0



Q6 Edge<sup>®</sup> Z



Q6 Edge<sup>®</sup> HD

Pretty in Pink	Lime Light	Yellow Blaze	Deep Purple	Orange Crush	Electric Blue
Red Corvette	White Satin	Grey Street	Midnight Blue	Silver Bullet	Back in Black

# Specifications

	EDGE 3	Q6 EDGE <sup>®</sup> 2.0	Q6 EDGE <sup>®</sup> Z	Q6 EDGE <sup>®</sup> HD
<b>DRIVE CONFIGURATION</b>	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive
<b>DRIVE WHEELS</b>	14" Solid <i>Optional:</i> 14" Pneumatic; 14" Knobby	14" Solid <i>Optional:</i> 14" Pneumatic; 14" Knobby	14" Run Flat <i>Optional:</i> 14" Pneumatic; 14" Knobby	14" Run Flat <i>Optional:</i> 14" Pneumatic; 14" Knobby
<b>CASTER WHEELS</b>	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A	N/A	N/A	N/A
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 6.25 mph (10 Km/h) 4.5 mph (7.2 Km/h) w/iLevel	Up to 6.25 (10 km/h) mph 3.5 mph (5.6 Km/h) w/iLevel	Up to 6 mph (9.6 km/h); 3.5 mph (5.6 Km/h) w/iLevel; <i>Optional:</i> 8 mph (12.8 Km/h) <sup>2</sup>	Up to 4.7 (7.6 Km/h) mph 3.5 mph (5.6 Km/h) w/iLevel
<b>GROUND CLEARANCE</b>	2.7" (battery tray)	2.7" (battery tray)	3" (frame base)	3" (frame center)
<b>TURNING RADIUS</b>	20.5"	6" casters: 20.5"	19.75"	19.75"
<b>OVERALL LENGTH</b>	32.23"	6" casters: 32.23"	35.6"	35.6"
<b>BASE WIDTH</b>	24.4"	24.4"	26.25"	26.25"
<b>Seating Sizes:</b>				
<b>SYNERGY<sup>®</sup> HD/HD TILT</b>	N/A	N/A	N/A	N/A
<b>TRU-BALANCE<sup>®</sup> 3 STATIC</b>	W: 10-22" D: 10-22"	W: 10-22" D: 10-22"	W: 12-24" D: 12-24"	N/A
<b>TRU-BALANCE<sup>®</sup> 3 POWER TILT</b>	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-24" D: 14-24"	N/A
<b>TRU-BALANCE<sup>®</sup> 3 LIFT &amp; TILT</b>	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-24" D: 12-24"	N/A
<b>TRU-BALANCE<sup>®</sup> 3 HD STATIC</b>	N/A	N/A	N/A	W: 20-32" D: 20-28"
<b>TRU-BALANCE<sup>®</sup> 3 HD POWER TILT</b>	N/A	N/A	N/A	W: 20-32" D: 20-25"
<b>TRU-BALANCE<sup>®</sup> 3 HD LIFT &amp; TILT</b>	N/A	N/A	N/A	W: 20-24" D: 20-25" <sup>6</sup>
<b>TRU-BALANCE<sup>®</sup> FLEX SEAT</b>	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	N/A	N/A
<b>KOZMO<sup>™</sup> PEDIATRIC STATIC SEAT</b>	N/A	N/A	N/A	N/A
<b>Seat-to-Floor Heights:</b>				
<b>SYNERGY<sup>®</sup> HD</b>	N/A	N/A	N/A	N/A
<b>TRU-BALANCE<sup>®</sup> 3 STATIC <sup>3,4</sup></b>	16.5" - 18.5"	16.5 - 18.5"	16.5"-18.5"	17.5 - 19.5"
<b>POWER TILT <sup>3,4</sup></b>	17.5" - 19.5"	17.5-19.5"	18.125"-20.125"	18.5-20.5"
<b>LIFT &amp; TILT <sup>3,4</sup></b>	17.5" - 19.5"	17.5-19.5"	19.375"-21.375"	17.75-19.75"
<b>TRU-BALANCE<sup>®</sup> FLEX SEAT <sup>3,4</sup></b>	16.5" - 19.5"	16.5-19.5"	N/A	N/A
<b>MANUAL RECLINE</b>	Yes	Yes	Yes	No
<b>BATTERY SIZE</b>	NF-22 (2)	NF-22	Group 24	Group 24
<b>BATTERY WEIGHT <sup>7</sup></b>	38 lbs. each	38 lbs.	52 lbs. each	52 lbs.
<b>AVAILABLE ELECTRONICS</b>	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+	100A Q-Logic 3 EX 75A Q-Logic NE+	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+
<b>BATTERY CHARGER</b>	8A Off-board	8A Off-board	8A Off-board	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole Mid-Wheel 6 <sup>®</sup>	2-Pole Mid-Wheel 6 <sup>®</sup>	4-Pole Mid-Wheel 6 <sup>®</sup> or 4-Pole 8 mph	4-Pole 4.7 mph
<b>WEIGHT CAPACITY</b>	300 lbs.	300 lbs.	300 lbs.	450 lbs. <sup>5</sup>
<b>BASE WEIGHT</b>	138.6 lbs.	138.6 lbs. (6" casters) 148 lbs. (8" casters)	143.5 lbs.	143.5 lbs.
<b>WARRANTY</b>	Lifetime limited warranty on frame; 2-year warranty on electronics and drive motor	Lifetime limited warranty on frame; 2-year warranty on electronics and drive motor	Lifetime limited warranty on frame; 2-year warranty on electronics and drive motor	Lifetime limited warranty on frame; 2-year warranty on electronics and drive motor

**FOOTNOTES:**

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) High-speed motor option will reduce range.
- 3) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 4) Power Recline and Power Front Riggings do not increase seat-to-floor height.
- 5) Weight capacity is reduced to 400 lbs. with iLevel<sup>®</sup> power adjustable seat height.
- 6) Required with high-speed motors.
- 7) Battery weight may vary based on manufacturer.
- 8) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4, standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.



Revolutionizing the Power Chair Experience





# 25,000 LIVES EMPOWERED



## 12" LIFT AT 4.5 MPH (7.2 Kmh)

Maximum lift of 10" on Q6 Edge® HD



12" of power adjustable seat height while at walking speed (10" on Q6 Edge® HD)

iLevel's patented Extra Stability Technology® enhances safety for transfers, reach, and numerous other activities of daily living

Up to 19 exciting color options

Retrofittable on the Edge 3 and Q6 Edge® 2.0

300 lbs. weight capacity on Edge 3, Q6 Edge® 2.0, Q6 Edge 2.0 X, and Q6 Edge® Z

400 lbs. weight capacity on Q6 Edge® HD

Smooth and stable ride

Front and rear LED fender lights come standard with iLevel®

iLevel® Power Chairs feature a revolutionary new power chair design for a better experience. iLevel improves daily functionality and quality of life by allowing operation of the power chair with the seat elevated while at walking speed up to 4.5 mph (7.2 Kmh).

## Benefits of iLevel®

**iLevel® Power Wheelchairs** feature a revolutionary new power chair design for a better experience. iLevel from Quantum® improves daily functionality and quality of life by allowing operation of the power chair with the seat fully elevated while at speeds at 4.5 mph (7.2 Kmh).

The world is not designed around seated height, and unfortunately wheelchair users for too long have been excluded from many everyday functional tasks and interactions beyond their physical reach. iLevel is a feature that can dramatically increase the quality of life for power chair users. From in the home to out in the community, a power elevating seat can help make life safer for users, make social interactions more fulfilling, and give people who use motorized wheelchairs access to common areas of their home.

iLevel's patented **Extra Stability Technology®** enhances safety for transfers, reach, and numerous other activities of daily living. iLevel elevating seat technology is unique in that it's user-inspired, where not only does it elevate, but it stabilizes the power chair for safe operation at at 4.5 mph (7.2 Kmh). For users, this means being able to grocery shop while elevated; this means navigating crowds at eye level; this means accessing high-top tables in restaurants; and, this means strolling the walking path arm-in-arm with your partner. Of course, simply being able to adjust the seat height, with added stability, makes transfers safer.

Specifically, iLevel can benefit you if:

- You have difficulty, or need assistance, to transfer from surfaces (beds, toilets, shower chairs, car/van seats, etc.) that are higher than the seat on motorized wheelchairs
- You have difficulty, or need assistance, to reach high enough to access one or more of the following areas at home, school, work, or the community:
  - ✓ Dresser drawers; clothes rod, washer or dryer
  - ✓ Medicine cabinet; bathroom sink/faucet; mirror; shower head/faucet
  - ✓ Freezer/refrigerator; oven/stove (including light/exhaust fan switches); microwave, cupboards/shelves;
  - ✓ counter; sink/faucet; disposal switch
  - ✓ Light switches; thermostat; fire alarm; wall phone; speaker; door eye hole/viewer; doorbell; elevator buttons
- You have difficulty seeing where you are going, moving through crowds, or safely crossing the street from a low (standard) seat height in motorized wheelchairs
- You experience neck pain, have difficulty seeing or hearing those standing nearby when seated in power wheelchairs
- You experience any difficulty sitting in motorized wheelchairs as a result of looking up or reaching overhead



## iLevel® ready units



### Edge 3

The design of the Edge 3 is for the most advanced power chair experience ever! Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. The Edge 3 is available in a total of 18 color options, including seven brand-new, bold colors.

With iLevel® option...

- Up to 12" of seat elevation at walking speed up to 4.5 mph (7.2 Kmh)



### Q6 Edge® 2.0

With features like standard 6 mph motors and advanced responsiveness, the Q6 Edge® 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.

With iLevel® option...

- Up to 12" of seat elevation at walking speed up to 3.5 mph (5.6 Kmh)



### Q6 Edge® HD

With features like standard 4-pole motors, Mid-Wheel 6® Drive Design and ATX Suspension, the Q6 Edge® HD is engineered to meet the performance needs of the most active user. The Q6 Edge HD accepts a complete range of seating and electronics options and offers an increased weight capacity of 450 lbs. to deliver exceptional rehab capability at an excellent value. It is also available with our iLevel® technology (400 lbs. weight capacity).

With iLevel® option...

- Up to 10" of seat elevation at walking speed up to 3.5 mph (5.6 Kmh)



### Q6 Edge® Z

The Q6 Edge® Z utilizes powerful 4-pole motors for maximum torque and Mid-Wheel 6® design for aggressive outdoor performance and precision indoor maneuverability. The Q6 Edge Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability, and is available with our iLevel® technology.

With iLevel® option...

- Up to 12" of seat elevation (in 16 seconds) at walking speed up to 3.5 mph (5.6 Kmh)



## Options

- Q-Logic 3 Advanced Drive Controls
- Trilogy vent tray
- LED full lighting package
- 4" wide flat-free knobby extended traction tires
- Swing-away joystick mount
- USB mobile device charger
- Transfer bar
- Cup holder
- Personal item hook
- Oxygen holder
- Push bars
- Backpack holder
- Phone holder

