



LC-525i shown in Aria® Cool Grey fabric

USB
Charger
Port



INFINITY COLLECTION

The Infinity Collection from Pride® Power Lift Recliners features true infinite-position seating and pocketed coils designed to deliver comfort, style and performance. Available in four different fabric selections, the Infinity Collection models deliver maximum relaxation and will match any décor.

NORTH AMERICA'S #1
POWER LIFT RECLINER

POWER
RECLINERS

Protect your investment with the Pride® Lift Chair Care Kit



Beamsville, Ontario | 888-570-1113
pridemobility.com | prideliftchairs.com

Pride
Mobility Products Co.
LIVE YOUR BEST®

INFINITY COLLECTION

NORTH AMERICA'S #1
POWER LIFT RECLINER

POWER LIFT
RECLINERS



LC-525i shown in
Ultraleather™ Charcoal fabric



	Petite Wide	Small	Medium	Large
SPECIFICATIONS	LC-525iPW²	LC-525iS²	LC-525iM²	LC-525iL²
	Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'2" & Below	Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'3" & Below	Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'4" - 6'0"	Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'9" - 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	19"	19"	20"	21"
Seat Depth	19"	19"	20"	22"
Seat Width	24"	19"	20"	22"
Top of Back to Seat	28"	26"	28"	30"
Back Style	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
Fully Padded Chaise	Yes	Yes	Yes	Yes
Heat and Massage	Optional	Optional	Optional	Optional
Seating Type	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs
Distance From Wall	19"	20"	20"	21"
Overall Chair Width	38"	34"	34.5"	37"
Chair Weight	130 lbs.	118 lbs.	130 lbs.	130 lbs.
Warranty	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*	7-Year Prorated*

INFINITY FABRIC SELECTION

QS

DURASOFT



EMBER



OAT



DEEP SKY



TIMBER

LEXIS STA-KLEEN® sta+kleen. PERFORMANCE FABRIC



BLACK



MUSHROOM



BURGUNDY



CHESTNUT

ULTRALEATHER™ Uf



BUFF



FUDGE



GARNET



PECAN



CHARCOAL

ARIA³ CRYPTON[®] super fabric



ESPRESSO



COOL GREY



RED



LAZULI



SAND

BENEFITS

- Pocketed coils deliver individualized support and relaxation
- Premium furniture grade laminate/hardwood frames for added strength and durability
- Easy-to-use hand control with USB charger port *(not available when adding optional heat and massage)*
- Standard lumbar pillow
- Standard head and arm covers
- Large, dual pockets for extra storage
- Quiet and smooth lift system
- Integrated battery backup ensures chair is operational when power goes out
- Energy-efficient power box equipped with self-tested electronics and 9' cord
- Footrest extension standard (LC-108 only)

OPTIONS*

- Deluxe heat and massage
- Residential cleaning kit
- Head pillow

Additional lead time may apply for some options.

*Additional charges apply for all options.

1. Batteries not included.
2. Available with deluxe heat and massage only. Standard heat and massage not available.
3. Crypton[®] fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.



Look for the Quick Ship Indicator for the shortest lead times available.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.
*Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.