

Jazzy®

ELITE ES PORTABLE

The Jazzy® Elite ES-Portable offers a combination of great style, comfort, performance and portability and disassembles into 5 convenient pieces for your active lifestyle. The in-line, front-wheel drive design allows tight turns around corners, with large front wheels for extra absorption over obstacles. With superb climbing capabilities and front anti-tip wheels, the Jazzy Elite ES Portable transitions over grade changes smoothly and handles various terrains with ease.



*Jazzy® Elite ES Portable
shown disassembled
in Jazzy Red*

SPECIFICATIONS

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.5"
Turning Radius	25.5"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 11 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	24.5 lbs. each
Base Weight	47.6 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

FEATURES

- Easy disassembly for convenient transport
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Height-adjustable footrest bracket
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Swing-away joystick
- Military patches
- XLR USB charger

RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander
- Zeus
- Commander

COLOR CHOICES



Candy
Apple Red



Viper Blue

OPTIONAL SEATING

- Solid seat pan



Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com

Pride[®]
Mobility Products Co.

Beamsville, Ontario • 888-570-1113 • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.