

# Jazzy<sup>®</sup> 2.0

## SELECT 6

The Jazzy Select 6<sup>®</sup> 2.0 offers a truly unique ride experience with a sporty design and superior performance. Pride's patented Active-Trac<sup>®</sup> Suspension combines with Mid-Wheel 6<sup>®</sup> drive design for a smooth, comfortable and stable ride, both indoors and out. The depth-adjustable high-back seat folds flat for greater convenience, while offering maximum comfort and support.

**THIS FOOT PLATFORM OFFERS SUPERIOR ADJUSTABILITY FOR AMPLE FOOT ROOM.**



**NORTH AMERICA'S  
FAVORITE  
POWER CHAIR**

Learn more online.



[gojazzy.com](http://gojazzy.com)

**THIS CONVENIENT FOLDING SEAT IS DEPTH ADJUSTABLE FOR MAXIMUM COMFORT AND OFFERS REPLACEABLE BACK AND SEAT COVERS.**

# Specifications

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed<sup>1</sup></b>	Up to 4.3 mph
<b>Ground Clearance<sup>2</sup></b>	2.5" at footplate, 1.77" at motor bracket
<b>Turning Radius<sup>2</sup></b>	23.76"
<b>Overall Length<sup>2</sup></b>	43.27"
<b>Overall Width</b>	23.27"
<b>Drive Wheels</b>	10" flat free
<b>Front Wheels</b>	6" flat free
<b>Rear Wheels</b>	6" flat free
<b>Suspension</b>	Active-Trac® (ATX) with Mid-Wheel 6® Technology
<b>Drivetrain</b>	Two-motor, Mid-Wheel 6® drive
<b>Braking System</b>	"Intelligent Braking" electronic regenerative, disc park brake
<b>Standard Electronics</b>	40 amp, PG GC 3 controller
<b>Battery Charger</b>	Off-board, 3.5 amp
<b>Per-Charge Range<sup>2,4,6</sup></b>	Up to 19 miles
<b>Battery Requirements<sup>3,5</sup></b>	(2) 12 volt, deep cycle, U-1
<b>Battery Weight</b>	24.5 lbs. (U-1)
<b>Base Weight</b>	170 lbs. without batteries
<b>Standard Seat Weight</b>	50 lbs. (comfort, high-back)
<b>Weight of Heaviest Piece</b>	114 lbs.
<b>Max. Seat Dimensions</b>	20" x 20" (high-back)

Colors:



Sapphire  
Blue



Ruby  
Red



Pearl  
White



Onyx  
Black

\*Additional charges may apply.

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>5</sup> Battery weight may vary based on manufacturer.

<sup>6</sup> Tested with Interceptor U1 batteries.

# Features

- Active-Trac® Suspension combines with Mid-Wheel 6® Technology for a smooth, stable ride over uneven terrain
- Convenient, limited recline, high-back black vinyl folding seat is depth adjustable for maximum comfort and includes replaceable back and seat covers available in sizes: 16" x 16"-18", 18" x 18"-20" and 20" x 18"-20"
- Height-adjustable, flip-up foot platform offers ample foot room for maximum comfort
- Black, flat-free, non-marking tires saves your floors
- In-line joystick mount is equipped with performance-matched electronics for easy operation
- Equipped with 60" lap belt for security and peace of mind

## Optional Seating\*

- Limited recline, high-back black vinyl solid seat pan with depth adjustment available in sizes: 16" x 16"-18" or 18" x 18"-20"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 20" x 18" or 20" x 20"
- Power elevating seat (VR2 6-key required)

## Options

- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away footrests
- Swing-away joystick
- Synergy® Seating

## Accessories

- Cane/crutch holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger