



TRAK™ with Medical Mounting Accessories

TRAK™




Oxygen Tank Holder

Accommodates "D" size Oxygen tank and is easily removed

LTV Bracket

Specifically designed to accommodate the Pulmonetics LTV 950



IV Pole

Equipment Tray

Rigid tray with modified frame for increased equipment clearances.

Folds with accessories in place!*

* Cargo may need to be removed from equipment tray.



Combined weight of accessories, cargo & occupant must not exceed chair capacity. Vent, Oxygen tank, and other cargo must be removed and secured off chair during transit in a motor vehicle.

Equipment Tray Details

	TRAK™ 12	TRAK™ 14
Tray Length	14.5 in. / 37 cm	14.5 in. / 37 cm
Tray Width	14.5 in. / 37 cm	16.5 in. / 37 cm
Clearance Height (0° of Seat Tilt)	12 in. / 30 cm	12 in. / 30 cm
Clearance Height (45° of Seat Tilt)	6 in. / 15 cm	6 in. / 15 cm
Tray Ground Clearance	3 in. / 7.5 cm	3 in. / 7.5 cm
O2 Ground Clearance	2 in. / 5 cm	2 in. / 5 cm

Part Numbers & Descriptions

3002-7016	IV Pole	All items must be factory installed.
3002-7017	Oxygen Tank Holder	
3002-7018	Equipment Tray	
3002-7019	LTV Bracket	
3002-7021	LTV Bracket & Equipment Tray	