

NM-435

FEATURES AND BENEFITS:

- Premium furniture-grade laminate/hardwood frames for added strength and durability
- 3-position full recline ideal for napping and relaxing
- Standard large dual pockets
- Standard footrest extension
- Easy-to-use hand control with USB charger port
- Removable back
- Energy-efficient transformer, features self-diagnostic electronics
- Quiet and smooth lift system
- Heavy-duty lift actuator and scissor mechanism
- Integrated battery backup in case of power failure

OPTIONS:

- Wireless remote*
- Lift chair overlay
- Residential cleaning kit
- Protective pad
- Non-skid leg levelers
- Head pillow
- Lumbar pillow



	Petite Small	Medium	Large/Tall
SPECIFICATIONS	NM-435PS	NM-435M	NM-435LT
	3-Position Chaise Lounger 5'3" and below	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 6'1" - 6'8"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	19"	20"	20"
Seat Depth	19"	20"	21"
Seat Width	20"	20"	22"
Top of Back to Seat	24.5"	26.5"	29"
Back Style	Blown-Fiber Waterfall Back	Blown-Fiber Waterfall Back	Blown-Fiber Waterfall Back
Fully Padded Chaise	Yes	Yes	Yes
Heat and Massage Option	Yes	Yes	Yes
Seating Type	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs
Distance From Wall	17"	19"	21"
Overall Chair Width	34"	34"	35"
Chair Weight	127 lbs.	118 lbs.	125 lbs.
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated



AVAILABLE FABRICS

BISHOP



OPTIONAL FABRICS

CRYPTON ARIA®



LEXIS STA-KLEEN®



WIRELESS REMOTE

*Performs lifted, seated, TV and full recline functions. Optimum range up to 4 feet with continuous operation up to 15 ft., designed to be used by the occupant of the lift chair



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

Look for the Quick Ship Indicator for the shortest lead times available.

Beamsville, Ontario • 888-570-1113
 pridemobility.com |

Pride
 Mobility Products Co.
 -- LIVE YOUR BEST --