



USB
Charger
Port



LC-580i shown in
Ultraleather® Charcoal fabric

Oasis

COLLECTION

The **Oasis** brings a new level of comfort to Pride® Power Lift Recliners. Overstuffed padding provides elegant comfort and design, while contoured memory foam seating and pocketed coils deliver individualized support and relaxation. Features include a unique quad pocket design for extra storage and an easy-to-use hand control with USB charger port. The Oasis is available in true-infinite position and 3-position.

NORTH AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS

Beamsville, Ontario | 888-570-1113
pridemobility.com | prideliftchairs.com

Pride
Mobility Products Co.
LIVE YOUR BEST®

Oasis

COLLECTION

NORTH AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS



LC-580i shown in
Saratoga Charcoal fabric



	Medium	Large
SPECIFICATIONS	LC-580iM	LC-580iL
	Infinite Position Trendelenburg Chaise Lounger 5'4" - 6'0"	Infinite Position Trendelenburg Chaise Lounger 5'9" - 6'2"
Weight Capacity	375 lbs.	375 lbs.
Seat-to-Floor	20"	21"
Seat Depth	21"	21"
Seat Width	20"	21"
Top of Back to Seat	26"	29"
Back Style	Blown Fiber Biscuit Back	Blown Fiber Biscuit Back
Fully Padded Chaise	Yes	Yes
Heat & Massage	Optional	Optional
Seating Type	Pocketed Coil Springs w/ Memory Foam Topper	Pocketed Coil Springs w/ Memory Foam Topper
Distance From Wall	20"	21"
Overall Chair Width	39"	39"
Chair Weight	130 lbs.	132 lbs.
Warranty	7-Year Pro-rated ¹	7-Year Pro-rated ¹

OASIS FABRIC SELECTION

QS

SARATOGA



CASHMERE



CHARCOAL



EGGPLANT



GODIVA



NAVY

LEXIS STA-KLEEN®

sta+kleen.
PERFORMANCE FABRIC



BLACK



MUSHROOM



BURGUNDY



CHESTNUT

ULTRALEATHER™

Uf
ULTRALEATHER™



BUFF



FUDGE



GARNET



PECAN



CHARCOAL

ARIA² CRYPTON®

CRYPTON®
super fabric



ESPRESSO



COOL GREY



RED



LAZULI



SAND

BENEFITS

- Premium furniture grade laminate/hardwood frames for added strength and durability
- Contoured memory foam seating and pocketed coils offer individualized comfort and support
- Sinuous no-sag back springs combine durability with just the right amount of flexibility for personalized comfort
- Overstuffed biscuit back design provides body-conforming comfort
- Specially designed arms afford a luxurious feel and firm support
- Industry-exclusive quad pocket design offers extra storage – easy access top pockets provide convenient storage for your hand control and other devices
- Easy-to-use hand control with USB charger port
(not available when adding optional heat and massage)
- Removable filling allows for custom lumbar support
- Headrest and armrest covers included

OPTIONS*

- Deluxe Heat and Massage (LC-580iM/L)
- Residential cleaning kit
- Head pillow
- Lumbar pillow

Additional lead time may apply for some options.

*Additional charges apply for all options.

1. Batteries not included.
2. Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.



Look for the Quick Ship Indicator for the shortest lead times available.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.
*Lifetime warranty on steel lift frame.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

US: 1-800-800-8586 • CAN: 1-888-570-1113
pridemobility.com |

Pride
Mobility Products Corp.
— LIVE YOUR BEST —