

The Outlander™ Exterior Lift from Pride® Lifts & Ramps was designed to provide a durable, easy-to-use system to take your travel mobility scooter with you. Simple, hitch-based installation requires little to no modification of your vehicle and allows for easy removal.

FEATURES

- 150 lbs. lifting capacity
- Integrated scooter lockdown arm
- Lifts Go-Go® Travel Mobility Scooters and a variety of other manufacturers' travel mobility scooters
- Available for Class II, III, and IV hitches
- On/off switch and toggle switch for up/down
- Electrical housing offset to the driver's side for ease of access
- Manual override
- Puddle light for ease of use at night
- Optional swing-away adaptor

Outlander™

PLATFORM TYPE	Travel Mobility Scooter
BASE WEIGHT ³	76 lbs.
MAXIMUM SCOOTER OVERALL WEIGHT	150 lbs.
MAXIMUM PRODUCT WIDTH	24"
OVERALL PLATFORM LENGTH	37"
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	21"
MOTOR	Fully sealed 12V DC
WARRANTY	3-year
HITCH CLASS ⁴	II, III, IV
VEHICLE PLATFORM	Full-size van, Minivan, Truck, Sedan, SUV and Crossover
LIFTS CURRENT PRIDE® PRODUCTS ^{1,2}	All Pride® Travel Mobility Scooters



1. Lifting weight does not include accessories or certain upgraded seating. Please contact Customer Service for greater detail.
2. Vehicle dependent.
3. Does not include any accessories or hardware.
4. Hitch weight carrying capacity varies by vehicle. Call for details.

Pride® Lifts & Ramps

www.prideliftsandramps.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

GO GO®
TRAVEL MOBILITY