

# EDGE 3

4.5 MPH iLEVEL® READY



Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. The Edge 3 is available in a total of 19 color options, including seven brand-new, bold colors.



RASPBERRY BERET

**Q**  
**QUANTUM®**  
#1 FOR REHAB POWER

 | [quantumrehab.ca](http://quantumrehab.ca)

# Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Solid Optional: 14" Pneumatic; 14" Knobby
<b>CASTER WHEELS</b>	Front: 6" Solid Rear: 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 6.25 mph (10 Kmh) Up to 4.5 mph (7.2 Kmh) w/iLevel activated
<b>GROUND CLEARANCE</b>	2.7" (battery tray)
<b>TURNING RADIUS <sup>2</sup></b>	20.5"
<b>OVERALL LENGTH <sup>2</sup></b>	32.23"
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b>	
<b>TRU-Balance® 3 Static Seat</b>	Width: 10-22" Depth: 10-22"
<b>TRU-Balance® 3 Power Tilt</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® 3 Tilt and iLevel®</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® 3 Flex Seat</b>	Width: 15-22" Depth: 15-20"
<b>SEAT-TO-FLOOR HEIGHTS <sup>10</sup></b>	
<b>TRU-Balance® 3 Static Seat <sup>5</sup></b>	16.5-18.5"
<b>TRU-Balance® 3 Power Tilt <sup>5</sup></b>	17.5-19.5"
<b>TRU-Balance® 3 Tilt &amp; iLevel® <sup>5</sup></b>	17.5-19.5"
<b>TRU-Balance® 3 Flex Seat</b>	16.5-19.5"
<b>MANUAL RECLINE <sup>6</sup></b>	Yes
<b>BATTERY SIZE <sup>3</sup></b>	NF-22 (2)
<b>BATTERY WEIGHT <sup>7</sup></b>	38 lbs. each
<b>AVAILABLE ELECTRONICS <sup>8</sup></b>	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole, Mid-Wheel <sup>6</sup>
<b>WEIGHT CAPACITY <sup>4</sup></b>	300 lbs.
<b>BASE WEIGHT</b>	138.6 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849), Single Power (K0856/K0857), Multiple Power (K0861)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability
- 10) 16.5" seat-to-floor height is not available with power tilt or iLevel®. 19.5" seat-to-floor height is not available with static seat or power recline.

# Features

- Upgraded SRS (Smooth Ride Suspension)
- 4.5 mph (7.2 Kmh) with iLevel®
- Fender lights and USB charger included
- Bold new colors
- Built in Canada

# Options

- TB3 Static and Power Positioning Seating System
- Flex Seating system
- High Back Seat
- Upgradable and specialty drive controls
- Accu-Trac Advanced Tracking Technology
- Manual and Power Elevating Leg rests
- Manual and Power Elevating Center-Mount Foot-platform
- Full light package
- See order form for additional accessories

# Colors

