



Q6 EDGE® Z
shown in *Back in Black*
with iLevel® technology

Q6 edge® Z

The Q6 Edge® Z utilizes powerful 4-pole motors for maximum torque and Mid-Wheel 6® design for aggressive outdoor performance and precision indoor maneuverability. The Q6 Edge® Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability, and is available with our iLevel® technology.

QUANTUM®

Specifications

DRIVE WHEELS	Standard: 14" Run flat Optional: 14" Pneumatic 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Standard: Up to 6 mph Optional: Up to 8 mph ¹⁰
GROUND CLEARANCE	3" at frame base
TURNING RADIUS ²	19.75"
OVERALL LENGTH ²	35.6"
BASE WIDTH	26.25"

SEATING SIZES		
SYNERGY ®	Width: 10-24"	Depth: 10-24"
TRU-BALANCE ® 3 STATIC	Width: 12-24"	Depth: 12-24"
TRU-BALANCE ® 3 POWER TILT	Width: 12-24"	Depth: 14-24"
TRU-BALANCE ® 3 LIFT	Width: 12-24"	Depth: 12-24"

SEAT-TO-FLOOR HEIGHTS	
SYNERGY ®/TRU-BALANCE® 3 STATIC ⁵	16.5-18.5"
TRU-BALANCE ® 3 POWER TILT ⁵	18.5"-20.5"
TRU-BALANCE ® 3 iLEVEL® & TILT ⁵	18.5"-20.5"

MANUAL RECLINE	Yes
BATTERY SIZE ³	Group 24
BATTERY WEIGHT ⁷	52 lbs. each
AVAILABLE ELECTRONICS ⁸	75A Q-Logic 2 NE+ 100A Q-Logic 2 EX
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-Pole Mid-Wheel 6®
WEIGHT CAPACITY	300 lbs.
BASE WEIGHT	143.5 lbs.

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Part of Community Use Package with 6 mph motors. Increases width by 1.5".
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Seat backs for units without TRU-Balance® 3 seating are only available in Back in Black.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability
- 10) iLevel® option not available if 8 mph motors are selected.

Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® 3 Power Positioning Systems
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- Compatible with iLevel®
- Standard switch-operated LED fender lights

Options

- iLevel® option, featuring patented Extra Stability Technology™, enables up to 12" of seat elevation while at walking speed (up to 3.5 mph)
- TRU-Balance® 3 Seating
- Q-Logic 2 Drive Controls
- 8 mph motor package
- Accu-Trac Advanced Tracking Technology⁹
- Trilogly vent tray⁹
- Compact vent tray⁹
- 4" wide flat-free or pneumatic knobby extended traction tires^{9,4}
- Stealth® Seating
- Contoured Seating
- Stealth® Positioning Products
- i-Drive® Electronics

Colors

Available on Q-Logic 2 controllers and units with TRU-Balance® 3 (base and back)⁶

