



ALSO AVAILABLE:
KOZMO™ PORTABLE



AVAILABLE WITH 12" iLEVEL®



AVAILABLE WITH 12" iLEVEL® AT 4.5 MPH



AVAILABLE WITH 12" iLEVEL®



AVAILABLE WITH 10" ELEVATION



AVAILABLE WITH 10" ELEVATION



AVAILABLE WITH 10" HD iLEVEL®



	KOZMO®	J6	Q4	Q6 EDGE® 2.0	EDGE 3	Q6 EDGE® Z	4FRONT™	QUANTUM® RIVAL®	Q6 EDGE® HD	QUANTUM® 1450
DRIVE CONFIGURATION	Rear-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Front-wheel Drive	Rear-wheel Drive	Mid-wheel Drive	Front-wheel Drive
DRIVE WHEELS	9" Solid	10" Solid Optional: 10" Pneumatic	14" Solid Optional: 14" Pneumatic	14" Solid Optional: 14" Pneumatic; 14" Knobby	14" Solid Optional: 14" Pneumatic; 14" Knobby	14" Run Flat Optional: 14" Pneumatic; 14" Knobby	14" Flat Free Optional: 14" Pneumatic; 14" Knobby	14" Solid Optional: 14" Pneumatic; 14" Knobby	14" Run Flat Optional: 14" Pneumatic; 14" Knobby	14" Solid
CASTER WHEELS	Front: 6" Solid	Front: 5" Solid Rear: 6" Solid	Front: 6" Solid Rear: 6" Solid	Front: 6" Solid Rear: 6" Solid	Front: 6" Solid Rear: 6" Solid	Front: 6" Solid Rear: 6" Solid	Rear: 9" Flat Free	Front: 8" Solid; 9" Solid; 9" Pneumatic; 10" Pneumatic	Front: 6" Solid Rear: 6" Solid	9" Solid
ANTI-TIP WHEELS	3" Solid	N/A	N/A	N/A	N/A	N/A	Front: 4"	3" Solid	N/A	4" Solid
MAXIMUM SPEED ³	Up to 3.6 mph	Up to 4 mph	Up to 6 mph	Up to 6.25 mph (3.5 mph w/iLevel)	Up to 6.25 mph (4.5 mph w/iLevel)	Up to 6 mph (3.5 mph w/iLevel) (Optional: 8 mph) ¹⁰	Up to 6 mph (3.2 mph w/10" elevated)	Up to 6 mph (1.5 mph while elevated) (Optional: 8 mph) ¹⁰	Up to 4.7 mph (3.5 mph w/iLevel)	Up to 4.5 mph
GROUND CLEARANCE	1.6" at anti-tip bracket	2.5" (frame)	2.5" (at caster pivot)	2.7" (battery tray)	2.7" (battery tray)	3" (frame base)	3" (motor pivot)	3.25" (frame base)	3" (frame center)	1.5" (anti-tip wheels)
TURNING RADIUS	21.875" without foot rigging	22"	20"	20.5"	20.5"	19.75"	24.5"	35.75"	19.75"	26.5"
OVERALL LENGTH	33.09" with foot platform	34.74"	35.5"	32.23"	32.23"	35.6"	40"	36"	35.6"	38.25"
BASE WIDTH	22.6"	23.19"	24"	24.4"	24.4"	26.25"	24"	25"	26.25"	29.25"
Seating Sizes:										
SYNERGY® HD/HD TILT	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	W: 20-32" D: 20-28"
TRU-BALANCE® 3 STATIC	N/A	W: 10-20" D: 10-20"	W: 16-20" D: 16-20"	W: 10-22" D: 10-22"	W: 10-22" D: 10-22"	W: 12-24" D: 12-24"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	N/A	N/A
TRU-BALANCE® 3 POWER TILT	N/A	W: 14-20" D: 14-20"	W: 16-20" D: 16-20"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-24" D: 14-24"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	N/A	N/A
TRU-BALANCE® 3 LIFT & TILT	N/A	N/A	N/A	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-24" D: 12-24"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	N/A	N/A
TRU-BALANCE® 3 HD STATIC	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	W: 20-28" D: 20-28"	N/A
TRU-BALANCE® 3 HD POWER TILT	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	W: 20-28" D: 20-25"	N/A
TRU-BALANCE® 3 HD LIFT & TILT	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	W: 20-28" D: 20-25" ⁸	N/A
TRU-BALANCE® FLEX SEAT	N/A	N/A	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	W: 15-22" D: 15-20"	N/A	N/A	W: 15-22" D: 15-20"	N/A	N/A
KOZMO™ PEDIATRIC STATIC SEAT	W: 12-16" D: 12-16"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Seat-to-Floor Heights^{6,13,14}:										
SYNERGY® HD	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	16.5-18.5"
TRU-BALANCE® 3 STATIC	N/A	16.75"	16.5-18.5"	16.5-18.5"	16.5" - 18.5"	16.5"-18.5"	15.25-17.25"	16.5-18.1"	17.5-19.5"	N/A
POWER TILT	N/A	16.75"	17.25"-19.25"	17.5-19.5"	17.5" - 19.5"	18.5"-20.5"	17.5-19.5"	18-19.6"	18.5-20.5"	18.125-20.125"
LIFT & TILT	N/A	N/A	N/A	17.5-19.5"	17.5" - 19.5"	18.5"-20.5"	17.75-19.75"	18.5-20.1"	17.75-19.75"	N/A
TRU-BALANCE® FLEX SEAT	N/A	N/A	16.5"-18.25"	16.5"-19.5"	16.5"-19.5"	N/A	N/A	16.5"-18.1"	N/A	N/A
KOZMO™ PEDIATRIC STATIC SEAT	15.75"-16.75"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
MANUAL RECLINE	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	No
BATTERY SIZE ⁴	(2) 18A, 12 V	U-1	NF-22	NF-22	NF-22 (2)	Group 24	Group 24: 51.5 lbs. Group 34: 38.5 lbs.	Group 24	Group 24	Group 24
BATTERY WEIGHT ¹¹	31 lbs.	23.4 lbs.	39 lbs.	38 lbs.	38 lbs. each	52 lbs. each	51.5 lbs.	52 lbs.	52 lbs.	52 lbs.
AVAILABLE ELECTRONICS ⁹	40A, PG GC3 (non-programmable)	70A Q-Logic NE 70A Q-Logic NE+ 70A Q-Logic 2 EX	120A Q-Logic 2 EX 75A Q-Logic NE+	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+	75A Q-Logic 3 NE+ 100A Q-Logic 3 EX	120A Q-Logic 3 EX	75A Q-Logic NE 75A Q-Logic NE+ 100A Q-Logic 2 EX	120A Q-Logic 3 EX 75A Q-Logic NE 75A Q-Logic NE+	90A VR2 6-key 100A Q-Logic 2 EX
BATTERY CHARGER	3.5A Off-board	5A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board
MOTOR PACKAGES	2 motor, rear-wheel drive	2-Pole Mid-Wheel 6®	4-Pole Mid-Wheel 6®	2-Pole Mid-Wheel 6®	2-Pole Mid-Wheel 6®	4-Pole Mid-Wheel 6® or 4-Pole 8 mph ¹	2-pole, front-wheel drive	4-Pole 6 mph ¹ or 4-Pole 8 mph ¹	4-Pole 4.7 mph	4-Pole 4.5 mph
WEIGHT CAPACITY ⁵	200 lbs.	300 lbs. ⁷	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	450 lbs. ¹²	600 lbs.
BASE WEIGHT	67 lbs.	100.2 lbs.	137 lbs.	138.6 lbs. (6" casters)	138.6 lbs.	143.5 lbs.	140 lbs.	129 lbs.	143.5 lbs.	195 lbs.

1) Optional
2) Requires HD Suspension
3) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
4) Battery capacity may change when adding power positioning systems.
5) Weight capacity is dictated by the selection of seating.

6) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
7) 250 lbs. weight capacity with TRU-Balance® Systems
8) A minimum seating size of 20" x 20" is required with iLevel® on the Q6 Edge® HD; maximum seating size is 24" x 25"
9) All Quantum units can be ordered with specialty controls

10) iLevel® not available with 8 mph motor package
11) Battery weight may vary +/- 3 lbs. based on manufacturer.
12) Weight capacity is reduced to 400 lbs. with iLevel® power adjustable seat height
13) Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.

14) 16.5" seat-to-floor height is not available with power tilt or iLevel®. 19.5" seat-to-floor height is not available with static seat or power recline.

