

| | 4FRONT™ <small>AVAILABLE WITH 10" LIFT</small> | Q6 EDGE® 2.0 <small>AVAILABLE WITH iLEVEL®</small> | Q6 EDGE® Z <small>AVAILABLE WITH iLEVEL®</small> | QUANTUM® RIVAL® | J6 | |
|-------------------------------------|---|---|---|--|---|---|
| DRIVE CONFIGURATION | Front-wheel Drive | Mid-wheel Drive | Mid-wheel Drive | Rear-wheel Drive | Mid-wheel Drive | |
| DRIVE WHEELS | 14" Flat Free | 14" Solid Optional: 14" Pneumatic or Knobby | 14" Run Flat Optional: 14" Pneumatic | 14" Solid Optional: 14" Pneumatic | 10" Solid Optional: 10" Pneumatic | |
| CASTER WHEELS | Rear: 9" Flat Free | Front: 6" Solid (Optional: 8") Rear: 6" Solid (Optional: 8") | Front: 6" Solid Rear: 6" Solid | Front: 9" Pneumatic | Front: 5" Solid Rear: 6" Solid | |
| ANTI-TIP WHEELS | Front: 4" | N/A | N/A | N/A | N/A | |
| MAXIMUM SPEED ⁴ | Up to 6 mph | Up to 6 mph | Up to 6 mph (Optional: 8 mph) ⁸ | Up to 6 mph (Optional: 8 mph) | Up to 4 mph | |
| GROUND CLEARANCE | 3" (motor pivot) | 2.7" (battery tray) | 3" (frame base) | 3.25" (frame base) | 2.5" (frame) | |
| TURNING RADIUS ⁵ | 24.5" | 6" casters: 20.5" 8" casters: 23.25" | 19.75" | 35.75" | 22" | |
| OVERALL LENGTH ¹¹ | 40" | 6" casters: 35.5" 8" casters: 40.48" | 35.6" | 36" | 34.74" | |
| BASE WIDTH | 24" | 24.4" | 26.25" | 25" | 23.19" | |
| Seating Sizes: | TRU-Balance® 3 | TRU-Balance® 3 | TRU-Balance® 3 | TRU-Balance® 2 | TRU-Balance® 3 | TRU-Balance® 2 |
| SYNERGY® | N/A | W: 10-22" D: 10-22" | W: 10-24" D: 10-24" | W: 10-22" D: 10-20" | W: 12-22" D: 12-22" | W: 10-20" D: 10-20" |
| STATIC | W: 12-22" D: 12-22" | W: 10-22" D: 10-22" | W: 12-24" D: 12-24" | W: 10-22" D: 10-20" | W: 12-22" D: 12-22" | W: 10-20" D: 10-20" |
| POWER TILT | W: 12-22" D: 12-22" | W: 12-22" D: 12-22" | W: 12-24" D: 14-24" | W: 14-22" D: 14-20" | W: 12-22" ³ D: 12-22" ³ | W: 14-20" D: 14-20" |
| LIFT & TILT | W: 12-22" D: 12-22" | W: 12-22" D: 12-22" | W: 12-24" D: 12-24" | N/A | W: 12-22" ³ D: 12-22" ³ | N/A |
| | | | | | | |
| MANUAL RECLINE | Yes | Yes | Yes | Yes | | Yes |
| BATTERY SIZE ⁶ | Group 24: 51.5 lbs. Group 34: 38.5 lbs. | NF-22 | Group 24 | Group 24 | | U-1 |
| BATTERY WEIGHT ¹⁰ | 51.5 lbs. | 38 lbs. | 52 lbs. each | 52 lbs. | | 23.4 lbs. |
| AVAILABLE ELECTRONICS ¹² | 120A Q-Logic 3 EX | 120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+ | 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX | 75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX | | 70A Q-Logic 2 NE 70A Q-Logic 2 NE+ 70A Q-Logic 2 EX |
| BATTERY CHARGER | 8A Off-board | 8A Off-board | 8A Off-board | 8A Off-board | | 5A Off-board |
| MOTOR PACKAGES | Two-pole, front-wheel drive | 2-Pole Mid-Wheel 6® | 4-Pole Mid-Wheel 6® | 4-Pole 6 mph ¹ or 4-Pole 8 mph ² | | 2-Pole 4 mph |
| WEIGHT CAPACITY ⁷ | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | | 300 lbs. ⁹ |
| BASE WEIGHT | 140 lbs. | 138.6 lbs. (6" casters) 148 lbs. (8" casters) | 143.5 lbs. | 129 lbs. | | 100.2 lbs. |

Footnotes:

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems
- 7) Weight capacity is dictated by the selection of seating.
- 8) iLevel® not available with 8 mph motor package.
- 9) 250 lbs. weight capacity with TRU-Balance® Systems
- 10) Battery weight may vary +/- 3 lbs. based on manufacturer
- 11) Without foot riggings
- 12) All Quantum units can be ordered with specialty controls

Available Colors:

| | Pretty in Pink | Red Corvette | Candy Apple Red | Orange Crush | Yellow Blaze | Lime Light | White | Silver Bullet | Grey Street | Black in Black | Night Train | Military Green | Electric Blue | Viper Blue | Midnight Blue | Deep Purple |
|----------|----------------|--------------|-----------------|--------------|--------------|------------|-------|---------------|-------------|----------------|-------------|----------------|---------------|------------|---------------|-------------|
| Edge 2.0 | ● | ● | | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | | ● | ● |
| 4Front | ● | ● | | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | | ● | ● |
| Edge Z | ● | ● | | ● | ● | ● | ● | ● | ● | ● | | ● | ● | | ● | ● |
| Rival™ | ● | ● | | ● | ● | ● | ● | ● | ● | ● | | | ● | | ● | ● |
| J6 | | | ● | | | | | | | | | | | ● | | |



(US) 866-800-2002 ▪ (Canada) 888-570-1113 ▪ quantumrehab.com

Due to manufacturing tolerances and continual product improvement, specifications can be subject to variance of + or – 3%. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



QUANTUM® FOR KIDS





QUANTUM[®] 4FRONT™

The 4Front™ features standard speeds of up to 6 mph in the lowered position, and speeds up to 3.2 mph while elevated at 10" with our optional safe seat elevation. Standard LED fender lights allow kids to see and be seen.

Features

- CASE (Caster Angle Sensor Encoder) technology
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Available with optional 10" of power adjustable safe seat elevation at speeds of up to 3.2 mph
- Q-Logic 3 EX Advanced Drive Controls
- LED front and rear fender lights (standard)
- Quiet ride technology
- Bold aesthetics with 14 color options



Q6 edge[®] 2.0

With features like standard 6 mph motors and advanced responsiveness, the Q6 Edge® 2.0 is engineered to meet the needs of the most active child. The Q6 Edge® 2.0 is available with our iLevel® technology and features bold aesthetics with 14 color options to fit any child's personality. It accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

Features

- iLevel® now available for seat sizes as small as 12" x 12" (SOF only)
- Standard 6 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics with 14 color options
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- iLevel® is retrofittable on the Q6 Edge® 2.0
- Available with Q-Logic 3 Advanced Drive Control System
- Available with WC-19 crash-test approved occupied transit kit option on tilt & recline and iLevel
- Front and rear LED fender lights (standard with iLevel®)
- Available with Accu-Trac Advanced Tracking
- Optional 8" casters



RIVAL[®]

The Rival® rear-wheel drive power chair features innovative design that provides exceptional outdoor performance and handling while delivering outstanding tight space maneuverability. The Rival accepts a complete range of rehab seating and electronics options.

Features

- 12 exciting color options
- Compact size for outstanding maneuvering in tight spaces
- Maximum speed up to 8 mph with high-speed motor package
- High-performance, high efficiency motor package
- Accepts a wide range of seating and positioning systems
- Two suspension packages
- Compatible with TRU-Balance® 3 Power Positioning Systems
- Sliding battery tray for easy access



Q6 edge[®] Z

The Q6 Edge® Z utilizes powerful 4-pole motors for maximum torque and Mid-Wheel 6® design for aggressive outdoor performance and precision indoor maneuverability. The Q6 Edge® Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability, and is available with our iLevel® technology.

Features

- User weight capacity of 300 lbs.
- Mid-Wheel 6® allows six wheels on the ground for stability
- 8 mph motor package option available
- iLevel® option available with up to 12" of power adjustable seat height at walking speed (3.5 mph)
- 13 exciting color options
- Front and rear LED fender lights (standard)
- Available with Accu-Trac Advanced Tracking Technology



J6

The compact J6 is the ultimate choice for clients who demand the superb stability of Mid-Wheel 6® Design and tight-quarter maneuverability. Highly adaptable, the J6 accepts a wide range of seating and electronics options while its powerful drivetrain delivers high-performance torque.

Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® 2 Tilt
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac® with extra stability) incorporates enhanced performance over more varied terrain
- Easy front access to batteries



Power Adjustable Seat Lift with iLevel[®]

Features

- iLevel® technology enables the unit to drive up to 3.5 mph while fully elevated
- Offers 12" of power adjustable seat lift
- Available in up to 14 back shroud colors
- Available on the Q6 Edge® 2.0 and Q6 Edge® Z
- Fender lights standard with this option



Power Tilt

Features

- 50° of tilt
- Lower standard seat-to-floor height
- Unmatched adjustment and durability
- Maintains center of gravity shift
- Available in up to 14 back shroud colors

TRU BALANCE[®] 3

POWER POSITIONING SYSTEMS



Power Recline

Features

- 85° - 175° of recline
- 300 lbs. weight capacity
- Does not increase seat-to-floor height
- Standard shear reduction and cantilever armrests
- Available in up to 14 back shroud colors
- Available with TRU-Balance® 3 Power Tilt and TRU-Balance® 3 Power Adjustable Seat Lift with iLevel®