



Facebook.com/
LeggeroUSA



3 SIZES AVAILABLE

16



12



14



REACH™



- ◆ Stylish Folding Frame
- ◆ Lightweight
- ◆ 3 Sizes Available
- ◆ Adjustable Fixed Tilt
- ◆ Knee Angle Adjustment
- ◆ Quick Hip Angle Adjustment
- ◆ Free Sunshade & Storage Bag



Distributed by Pride Mobility Products Company, Beamsville, ON, Canada
For customer enquiries, email Productorders@pridemobility.com or call us at 888-570-1113.
Web: Leggero.us • facebook.com/LeggeroUSA • © 2019 LEGGERO LLC CAREACHB-1909

REACH™



Optional REACH RCB™
Back with Swing
In/Out Laterals

Size Options	Max Seat Width
REACH™ 12	12 in / 30 cm
REACH™ 14	14 in / 35 cm
REACH™ 16	16 in / 40 cm

Adjustable Fixed Tilt

10°
20°
30°



SPECIFICATIONS & ADJUSTMENT RANGE

	REACH™ 12	REACH™ 14	REACH™ 16
Weight Limit	80 lbs / 36 kg	110 lbs / 50 kg	150 lbs / 68 kg
Seat Width	8-12 in / 20-30 cm	10-14 in / 25-35 cm	12-16 in / 30-40 cm
Seat Depth	8-13 in / 20-33 cm	11-16 in / 28-40 cm	13-18 in / 33-46 cm
Back Height	11-18 in / 28-46 cm	13-21 in / 33-53 cm	16-24 in / 40-61 cm
Seat to Footplate	8-13 in / 20-33 cm	11-16 in / 28-40 cm	13-18 in / 33-46 cm
Adjustable Fixed Tilt Positions	10°, 20°, or 30°	10°, 20°, or 30°	10°, 20°, or 30°
Seat to Back Angle Range with 10° Fixed Tilt	85°, 90°, 95°, 100°, 105°, 110°, 115°	85°, 90°, 95°, 100°, 105°, 110°, 115°	85°, 90°, 95°, 100°, 105°, 110°, 115°
Seat to Back Angle Range with 20° Fixed Tilt	85°, 90°, 95°, 100°, 105°	85°, 90°, 95°, 100°, 105°	85°, 90°, 95°, 100°, 105°
Seat to Back Angle Range with 30° Fixed Tilt	85°, 90°, 95°	85°, 90°, 95°	85°, 90°, 95°
Weight of Chair	27 lbs / 12 kg	29 lbs / 13 kg	31 lbs / 14 kg
Maximum Safe Slope	5°	5°	5°
Rear Wheel Diameter	12 in / 30 cm	12 in / 30 cm	12 in / 30 cm
Front Caster Diameter	7 in / 18 cm	7 in / 18 cm	7 in / 18 cm
Chair Width	21.5 in / 55 cm	23.5 in / 60 cm	25.5 in / 65 cm
Folded Dimensions	31 x 21.5 x 18 in / 79 x 55 x 46 cm	31 x 23.5 x 18 in / 79 x 60 x 46 cm	33 x 25.5 x 18 in / 84 x 65 x 46 cm

For more information visit www.Leggero.us