

RUGGED WITH FEATHER-TOUCH DISASSEMBLY



4-wheel

REVO™ 2.0



Experience the **Revo™ 2.0**, the most impressive mid-size scooter that combines an array of features into one great package. Tough and built to last, the Revo 2.0 offers rugged dependability you expect from luxury mid-size scooters. Comfort-Trac Suspension (CTS) ensures a smooth and comfortable ride over varied terrain. Experience convenient, feather-touch disassembly for portability that is comparable to our Travel Mobility line. Enjoy standard under seat storage and speeds up to 5 mph. Ride on with the Revo 2.0!



PRIDE®
MOBILITY
SCOOTERS

| pridemobility.ca |





Specifications 3W: S66 4W: S67

Weight Capacity	375 lbs.
Maximum Speed	5 mph
Ground Clearance	3W: 3.25" at mid-frame 4W: 3.5" at mid-frame
Turning Radius	3W: 40" 4W: 52"
Overall Length	3W: 46" 4W: 47"
Overall Width	20.75"
Seat-to-Ground Height	23.5" highest position; 21" lowest position
Seat-to-Deck Height	15"-18"
Front Tires	3W: 3" x 8" 4W: 3" x 9"
Rear Tires	3" x 9"
Suspension	Front and rear
Drive System	24 volt DC sealed transaxle
Dual Braking System	Regenerative and electro-mechanical
Range Per Charge ^{1,2}	3W: Up to 12 miles 4W: Up to 13 miles
Total Weight Without Batteries ^{3,5}	3W: 120 lbs. 4W: 132 lbs.
Heaviest Piece When Disassembled ³	3W: 48 lbs. (<i>rear section</i>) 4W: 57 lbs. (<i>front section</i>)
Standard Seating	<i>Type: Foldable, molded plastic</i> <i>Dimensions: Width: 20"</i> <i>Depth: 18"</i>
Battery Requirements ⁴	(2) U1 batteries; 25 lbs. each
Battery Charger	3.5 amp
Color	 

Grey Street True Blue

Features

- Comfort-Trac front and rear suspension
- Delta tiller with wraparound handles
- Easily accessible tie down points (*for transport of unoccupied scooter*)
- Easily charge the batteries on or off the unit
- ECO mode/sport mode selector switch offers an optional ECO mode at 4 mph with improved range per charge or the sport mode at 5 mph
- Under-tiller puddle light improves visibility to USB charging port and tiller-mounted storage
- Tiller-mounted headlight and rear running lights
- Pride's exclusive black, flat-free, non-scuffing tires
- Seat designed for extreme comfort
- Standard front basket
- Standard under seat storage
- Ultra heavy-duty drivetrain with 70 amp controller
- User-friendly console
- Easy access tiller mount XLR USB charging port for mobile devices
- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system

Accessories

- Cane/crutch holder (*single, double*)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 50", 60", 70" Lap belt
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

⁵ Includes standard seat weight.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

