

Footnotes

- All shroud panels included.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Battery weight not included.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- Battery pack is 18 lbs. with 12AH batteries.
- Durable solid plastic.
- Deep cycle (AGM or gel cell type recommended).
- 18 AH battery pack weighs 28 lbs.
- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Includes pathway illumination light.
- Battery weight may vary based on manufacturer.
- Width of base (not including seat).
- Does not account for anti-tip clearance.
- Includes standard seat weight.
- Option available through most providers. Not available through Pride Mobility®.
- Optional.
- At midframe.
- Light package.

PRIDE® MOBILITY SCOOTERS

LIVE YOUR BEST®

END-TO-END



	Go Chair®	Go-Go® Ultra X	Go-Go Elite Traveller®	Go-Go® Sport ¹⁰	Go-Go® LX with CTS Suspension ¹⁰	Jazzy® Zero Turn	Go-Go® Folding Scooter ¹⁰	Revo® 2.0 ¹⁰	Victory® 10 ¹⁰	Victory® Twin ¹⁸	Victory® LX with CTS Suspension ¹⁸	Celebrity® LR	Celebrity® DX Sport	Pursuit® ¹⁸	Maxima	Pursuit® XL ¹⁸	Raptor® ¹⁸
FRONT TIRES	3" x 6" SOLID	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 7" SOLID	1.6" x 7" SOLID	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	3.5" x 10" SOLID	DUAL 3" x 9" SOLID	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	4" x 10" PNEUMATIC	4" x 10" PNEUMATIC	4" x 13" PNEUMATIC	3" x 10" SOLID	13" PNEUMATIC	3" x 16" PNEUMATIC
REAR TIRES	3" x 9" SOLID	2.5" x 8" SOLID	2.5" x 8" SOLID	3" x 9" SOLID	2.5" x 8" SOLID	2.5" x 9" SOLID	2" x 8" SOLID	3" x 9" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	4" x 10" PNEUMATIC	4" x 13" PNEUMATIC	4" x 13" PNEUMATIC	4" x 10" SOLID	13" PNEUMATIC	3" x 16" PNEUMATIC
MAXIMUM SPEED²	UP TO 3.7 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4.7 MPH	UP TO 4.4 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 5 mph	UP TO 5.25 MPH	UP TO 5.25 MPH	UP TO 5.6 MPH	UP TO 6.5 MPH	UP TO 9.5 MPH	UP TO 8 MPH HD: UP TO 6.25 MPH	UP TO 5.25 MPH	UP TO 8.5 MPH	UP TO 14 MPH
GROUND CLEARANCE^{4,13,17}	1.6" at anti-tip bracket 1.9" at frame	2.25" at mid-frame 1.25" at motor	2"	3" at midframe 1.5" at motor	1.5" at motor	1.3" at motor	2.3"	1.25" at transaxle	3" at mid-frame	3" at mid-frame	2.2" at motor	2"	4.25"	4.5"	4.25" 3W 3.25" 4W	5.25"	3.5" at center of frame
3-W TURNING RADIUS⁴	N/A	33"	33"	36"	33"	N/A	N/A	40"	42"	N/A	N/A	44.5"	N/A	N/A	50.38"	N/A	60"
4-W TURNING RADIUS⁴	25.4" WITH FOOT PLATFORM	44"	44"	50.5"	44.75"	38.25"	32"	52"	54"	42"	59.9"	56.5"	66.25"	67.75" LEFT/66.75" RIGHT	54"	82.5"	N/A
3-W LENGTH/WIDTH^{4,12}	N/A	37"/19.5"	37"/19.5"	42"/21.25"	39.75"/20.5"	N/A	N/A	46"/21"	45.25"/22.25"	N/A	N/A	47.5"/24"	N/A	N/A	46"/25.5"	N/A	59"/24.75"
4-W LENGTH/WIDTH^{4,12}	28.5" WITHOUT FOOT PLATFORM/22.6"	39.5"/19.5"	39.5"/19.5"	43"/21.25"	41"/20.5"	45.75"/20"	37"/19"	47"/21"	47"/22.25"	46"/26"	47"/22.2"	48"/24"	47.5"/25.25"	53.25"/25.25"	46.9"/25.5"	56"/27.25"	N/A
SEAT-TO-GROUND HEIGHT RANGE	22"-23"	18.13"-20.13"	18.13"-20.13"	20.5"-22"	20.75"-22.25"	N/A	20.5"	21"-23.5"	20.5"-22"	23"-24.5"	24"-26"	22.5"-24"	22.75"-24.75"	24.75"-26.25"	22.5"-24.5"	26.25"-27.75"	28"
SEAT TYPE (STANDARD)	COMPACT, FOLDABLE (18"X17") 5-POINT SWIVEL WITH LATCH RELEASE	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	FOLDABLE, MOLDED PLASTIC (17"X17")	COMPACT, FOLDABLE (16"X13")	COMPACT, FOLDABLE (20"X18")	FOLDABLE, MOLDED PLASTIC (18"X17")	FOLDABLE, MOLDED PLASTIC (18"X17")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	DELUXE CONTOURED MEDIUM-BACK (22"X18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH BACK (20" X 18") (100" LIMITED RECLINE)
(3/4-WHEEL) WEIGHT^{3,4,14}	67 LBS.	76 LBS./80 LBS.	76 LBS./80 LBS.	92.8 LBS./99.2 LBS.	85.4 LBS./87.4 LBS.	120 LBS.	51.4 LBS.	120 LBS./131 LBS.	133.5 LBS./137.5 LBS.	N/A/150 LBS.	N/A/179 LBS.	160 LBS./173 LBS.	N/A/200 LBS.	N/A/225 LBS.	171 LBS./176 LBS.	259.9 LBS.	208 LBS.
SEAT-TO-DECK HEIGHT RANGE	N/A	15.75"-17.25"	15.75"-17.25"	17"-18"	16.25"-17.25"	17.5"	13"	18"	17.5"-19"	17.5"-19"	17"-19"	16.5"-18.5"	16.75"-18.75"	17"-18.75"	3W: 22.3"-22.6" 4W: 24.3"-24.6"	21"-22.5"	19"
BATTERY WEIGHT¹¹	31 LBS. BATTERY PACK	9 LBS. EA. ⁵	14 LBS. EA. ⁵	14 LBS. EA. ⁵	14 LBS. EA. ⁵	13.5 LBS. EA.	9 LBS. EA. (SEALED LEAD ACID) /5 LBS. EA. (LITHIUM-ION)	24 LBS. EA.	24.5 LBS. EA./ 32.5 LBS. EA.	24.5 LBS. EA./ 32.5 LBS. EA.	24.5 LBS. EA./ 32.5 LBS. EA.	36.5 LBS. EA.	24.5 LBS. EA./ 32.5 LBS. EA.	38 LBS. EA./ 53 LBS. EA.	38 LBS. EA. (NF-22 (55AH)) 24.5 LBS. EA. (U-1)	70 LBS. EA. (G31) 60 LBS. EA. (G27) 53 LBS. EA. (G24)	53 LBS. EA.
BATTERY REQUIREMENTS⁷	(2)18AH	(2)12AH	(2)18AH (INCLUDED)	(2)18AH (INCLUDED)	(2)18AH	(2) 22AH	(2) 12AH OR (1) 11.5 AH (LITHIUM-ION)	(2) U1	(2) U1 (INCLUDED) OR (2) 40AH ¹⁶	(2) U1 (INCLUDED) OR (2) 40AH ¹⁶	(2) U1 OR (2) 40AH	(2) NF-22	(2) U1 OR (2) 40AH	(2) NF-22 OR (2) GROUP 24 ¹⁶	(2) NF-22 (INCLUDED) (2) U-1 (POWER ELEVATING SEAT)	(2) GROUP 31 ¹⁵ , (110AH) (2) GROUP 27 (100AH) OR (2) GROUP 24 (75 AH)	(2) 75AH
PER CHARGE RANGE (MAX.)^{3,5}	UP TO 8.74 MILES	3W: UP TO 6.9 MILES 4W: UP TO 7.2 MILES	3W: UP TO 10.4 MILES 4W: UP TO 9.7 MILES	3W: UP TO 12.7 MILES 4W: UP TO 12 MILES	3W: UP TO 9.44 MILES 4W: UP TO 10.5 MILES	UP TO 10 MILES	4W: UP TO 12.5 MILES (LITHIUM-ION) 4W: UP TO 9.3 MILES (SEALED LEAD ACID)	3W: UP TO 12 MILES 4W: UP TO 13 MILES	UP TO 15.5 MILES	UP TO 11.1 MILES (WITH U-1 BATTERIES)	UP TO 16.16 MILES	UP TO 17.99 MILES (3-W) UP TO 18.97 MILES (4-W)	UP TO 11 MILES	UP TO 25 MILES	3W: UP TO 18 MILES, U-1 (POWER ELEVATING SEAT) UP TO 9.6 MILES 4W: UP TO 16 MILES	G31: UP TO 31.83 MILES G27: UP TO 29.51 MILES G24: UP TO 22.13 MILES	19 MILES (UP TO 400 LBS. AND FULL SPEED) 31 MILES (UP TO 220 LBS. AND AT LOW SPEED)
BATTERY CHARGER	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 3.5A	OFF-BOARD, 2A (SEALED LEAD ACID) OFF-BOARD, 2.5A (LITHIUM-ION)	OFF-BOARD, 3.5A	OFF-BOARD, 3A (U-1) OFF-BOARD, 5A (40 AH)	OFF-BOARD, 3A (U-1) OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 8A	OFF-BOARD, 5A	OFF-BOARD, 8A	OFF-BOARD, 8A
SUSPENSION	LIMITED	N/A	N/A	N/A	FRONT & REAR	FRONT & REAR	N/A	FRONT & REAR	N/A	N/A	FRONT & REAR	3-W REAR/ 4-W FRONT & REAR	REAR	FRONT & REAR	3-W REAR/ 4-W FRONT & REAR	FRONT & REAR	FRONT & REAR
MAXIMUM WEIGHT CAPACITY	300 LBS.	260 LBS.	300 LBS.	325 LBS.	300 LBS.	300 LBS.	250 LBS.	375 LBS.	400 LBS.	350 LBS.	400 LBS.	350 LBS.	350 LBS.	400 LBS.	500 LBS.	400 LBS.	400 LBS.
STANDARD BODY COLORS	RUBY RED, SAPPHIRE BLUE, ROSE QUARTZ, CITRINE YELLOW, AMBER ORANGE, PEARL WHITE	RED ⁶	RED, BLUE ¹	RED, BLUE ¹	RED, BLUE ¹	FIRE OPAL, GREEN MOONSTONE, ONYX BLACK, PEWTER, PEARL WHITE, SAPPHIRE BLUE	WHITE	GREY STREET, TRUE BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, ONYX BLACK	CANDY APPLE RED, ONYX BLACK	CANDY APPLE RED, ONYX BLACK	CANDY APPLE RED	CANDY APPLE RED, SILVER	SILVER
WEIGHT OF HEAVIEST PIECE⁴	36 LBS. (REAR)	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	35 LBS. (3-WHEEL) (FRONT SECTION) 44 LBS. (4-WHEEL) (FRONT SECTION)	34 LBS. (3-WHEEL) (FRONT SECTION) 36.6 LBS. (4-WHEEL) (FRONT SECTION)	48 LBS. (REAR SECTION)	51.4 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK)	48 LBS. (3-WHEEL) (REAR SECTION) 57 LBS. (4-WHEEL) (FRONT SECTION)	54.5 LBS. (3-WHEEL) (FRONT SECTION) 61 LBS. (4-WHEEL) (FRONT SECTION)	55 LBS. (FRONT SECTION)	77.6 LBS. (FRONT SECTION)	61 LBS. (3-WHEEL) (FRONT SECTION) 73 LBS. (4-WHEEL) (FRONT SECTION)	83 LBS. (FRONT SECTION)	95 LBS. (FRONT SECTION)	76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)	205.5 LBS. (MAIN FRAME)	N/A (UNIT DOES NOT DISASSEMBLE)
	TRAVEL MOBILITY							VALUE SERIES			MID-SIZE LUXURY			HEAVY-DUTY		FULL-SIZE LUXURY	