

**Footnotes**

- All shroud panels included.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Battery weight not included.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- Battery pack is 18 lbs. with 12AH batteries.
- Durable solid plastic.
- Deep cycle (AGM or gel cell type recommended).
- 18 AH battery pack weighs 28 lbs.
- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

- Includes pathway illumination light.
- Battery weight may vary based on manufacturer.
- Width of base (not including seat).
- Does not account for anti-tip clearance.
- Includes standard seat weight.
- Option available through most providers. Not available through Pride Mobility®.
- Optional.
- At midframe.
- Light package.

# PRIDE® MOBILITY SCOOTERS

## LIVE YOUR BEST®

# END-TO-END



	Go Chair®	Go-Go® Ultra X	Go-Go Elite Traveller®	Go-Go® Sport <sup>10</sup>	Go-Go® LX with CTS Suspension <sup>10</sup>	Go-Go® Folding Scooter <sup>10</sup>	Victory® 10 <sup>10</sup>	Victory® Twin <sup>18</sup>	Victory® LX with CTS Suspension <sup>18</sup>	Celebrity® LR	Celebrity® DX Sport	Pursuit® <sup>18</sup>	Maxima	Pursuit® XL <sup>18</sup>	
<b>FRONT TIRES</b>	3" x 6" SOLID	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	1.6" x 7" SOLID	3.5" x 10" SOLID	DUAL 3" x 9" SOLID	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	4" x 10" PNEUMATIC	4" x 10" PNEUMATIC	4" x 13" PNEUMATIC	3" x 10" SOLID	13" PNEUMATIC	
<b>REAR TIRES</b>	3" x 9" SOLID	2.5" x 8" SOLID	2.5" x 8" SOLID	3" x 9" SOLID	2.5" x 8" SOLID	2" x 8" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	4" x 10" PNEUMATIC	4" x 13" PNEUMATIC	4" x 13" PNEUMATIC	4" x 10" SOLID	13" PNEUMATIC	
<b>MAXIMUM SPEED<sup>2</sup></b>	UP TO 3.7 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4.7 MPH	UP TO 4.4 MPH	UP TO 4 MPH	UP TO 5.25 MPH	UP TO 5.25 MPH	UP TO 5.6 MPH	UP TO 6.5 MPH	UP TO 9.5 MPH	UP TO 8 MPH HD: UP TO 6.25 MPH	UP TO 5.25 MPH	UP TO 8.5 MPH	
<b>GROUND CLEARANCE<sup>4,13,17</sup></b>	1.6" at anti-tip bracket 1.9" at frame	2.25" at mid-frame 1.25" at motor	2"	3" at midframe 1.5" at motor	1.5" at motor	2.3"	3" at mid-frame	3" at mid-frame	2.2" at motor	2"	4.25"	4.5"	4.25" 3W 3.25" 4W	5.25"	
<b>3-W TURNING RADIUS<sup>4</sup></b>	N/A	33"	33"	36"	33"	N/A	42"	N/A	N/A	44.5"	N/A	N/A	50.38"	N/A	
<b>4-W TURNING RADIUS<sup>4</sup></b>	25.4" WITH FOOT PLATFORM	44"	44"	50.5"	44.75"	32"	54"	42"	59.9"	56.5"	66.25"	67.75" LEFT/66.75" RIGHT	54"	82.5"	
<b>3-W LENGTH/WIDTH<sup>4,12</sup></b>	N/A	37"/19.5"	37"/19.5"	42"/21.25"	39.75"/20.5"	N/A	45.25"/22.25"	N/A	N/A	47.5"/24"	N/A	N/A	46"/25.5"	N/A	
<b>4-W LENGTH/WIDTH<sup>4,12</sup></b>	28.5" WITHOUT FOOT PLATFORM/22.6"	39.5"/19.5"	39.5"/19.5"	43"/21.25"	41"/20.5"	37"/19"	47"/22.25"	46"/26"	47"/22.2"	48"/24"	47.5"/25.25"	53.25"/25.25"	46.9"/25.5"	56"/27.25"	
<b>SEAT-TO-GROUND HEIGHT RANGE</b>	22"-23"	18.13"-20.13"	18.13"-20.13"	20.5"-22"	20.75"-22.25"	20.5"	20.5"-22"	23"-24.5"	24"-26"	22.5"-24"	22.75"-24.75"	24.75"-26.25"	22.5"-24.5"	26.25"-27.75"	
<b>SEAT TYPE (STANDARD)</b>	COMPACT, FOLDABLE (18"X17") 5-POINT SWIVEL WITH LATCH RELEASE	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (16"X13")	FOLDABLE, MOLDED PLASTIC (18"X17")	FOLDABLE, MOLDED PLASTIC (18"X17")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	DELUXE CONTOURED MEDIUM-BACK (22"X18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	
<b>(3/4-WHEEL) WEIGHT<sup>3,4,14</sup></b>	67 LBS.	76 LBS./80 LBS.	76 LBS./80 LBS.	92.8 LBS./99.2 LBS.	85.4 LBS./87.4 LBS.	51.4 LBS.	133.5 LBS./137.5 LBS.	N/A/150 LBS.	N/A/179 LBS.	160 LBS./173 LBS.	N/A/200 LBS.	N/A/225 LBS.	171 LBS./176 LBS.	259.9 LBS.	
<b>SEAT-TO-DECK HEIGHT RANGE</b>	N/A	15.75"-17.25"	15.75"-17.25"	17"-18"	16.25"-17.25"	13"	17.5"-19"	17.5"-19"	17"-19"	16.5"-18.5"	16.75"-18.75"	17"-18.75"	3W: 22.3"-22.6" 4W: 24.3"-24.6"	21"-22.5"	
<b>BATTERY WEIGHT<sup>11</sup></b>	31 LBS. BATTERY PACK	9 LBS. EA. <sup>5</sup>	14 LBS. EA. <sup>5</sup>	14 LBS. EA. <sup>5</sup>	14 LBS. EA. <sup>5</sup>	9 LBS. EA. (SEALED LEAD ACID)/5 LBS. EA. (LITHIUM-ION)	24.5 LBS. EA./32.5 LBS. EA.	24.5 LBS. EA./32.5 LBS. EA.	24.5 LBS. EA./32.5 LBS. EA.	36.5 LBS. EA.	24.5 LBS. EA./32.5 LBS. EA.	38 LBS. EA./53 LBS. EA.	38 LBS. EA. (NF-22 (55AH))/24.5 LBS. EA. (U-1)	70 LBS. EA. (G31) 60 LBS. EA. (G27) 53 LBS. EA. (G24)	
<b>BATTERY REQUIREMENTS<sup>7</sup></b>	(2)18AH	(2)12AH	(2)18AH (INCLUDED)	(2)18AH (INCLUDED)	(2)18AH	(2) 12AH OR (1) 11.5 AH (LITHIUM-ION)	(2) U1 (INCLUDED) OR (2) 40AH <sup>16</sup>	(2) U1 (INCLUDED) OR (2) 40AH <sup>16</sup>	(2) U1 OR (2) 40AH	(2) NF-22	(2) U1 OR (2) 40AH	(2) NF-22 OR (2) GROUP 24 <sup>16</sup>	(2) NF-22 (INCLUDED) (2) U-1 (POWER ELEVATING SEAT)	(2) GROUP 31 <sup>16</sup> , (110AH) (2) GROUP 27 (100AH) OR (2) GROUP 24 (75 AH)	
<b>PER CHARGE RANGE (MAX.)<sup>2,8</sup></b>	UP TO 8.74 MILES	3W: UP TO 6.9 MILES 4W: UP TO 7.2 MILES	3W: UP TO 10.4 MILES 4W: UP TO 9.7 MILES	3W: UP TO 12.7 MILES 4W: UP TO 12 MILES	3W: UP TO 9.44 MILES 4W: UP TO 10.5 MILES	4W: UP TO 12.5 MILES (LITHIUM-ION) 4W: UP TO 9.3 MILES (SEALED LEAD ACID)	UP TO 15.5 MILES	UP TO 11.1 MILES (WITH U-1 BATTERIES)	UP TO 16.16 MILES	UP TO 17.99 MILES (3-W) UP TO 18.97 MILES (4-W)	UP TO 11 MILES	UP TO 25 MILES	3W: UP TO 18 MILES, U-1 (POWER ELEVATING SEAT) UP TO 9.6 MILES 4W: UP TO 16 MILES	G31: UP TO 31.83 MILES G27: UP TO 29.51 MILES G24: UP TO 22.13 MILES	
<b>BATTERY CHARGER</b>	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A (SEALED LEAD ACID) OFF-BOARD, 2.5A (LITHIUM-ION)	OFF-BOARD, 3A (U-1) OFF-BOARD, 5A (40 AH)	OFF-BOARD, 3A (U-1) OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 8A	OFF-BOARD, 5A	OFF-BOARD, 8A	
<b>SUSPENSION</b>	LIMITED	N/A	N/A	N/A	FRONT & REAR	N/A	N/A	N/A	FRONT & REAR	3-W REAR/ 4-W FRONT & REAR	REAR	FRONT & REAR	3-W REAR/ 4-W FRONT & REAR	FRONT & REAR	
<b>MAXIMUM WEIGHT CAPACITY</b>	300 LBS.	260 LBS.	300 LBS.	325 LBS.	300 LBS.	250 LBS.	400 LBS.	350 LBS.	400 LBS.	350 LBS.	350 LBS.	400 LBS.	500 LBS.	400 LBS.	
<b>STANDARD BODY COLORS</b>	RUBY RED, SAPPHIRE BLUE, ROSE QUARTZ, CITRINE YELLOW, AMBER ORANGE, PEARL WHITE	RED <sup>9</sup>	RED, BLUE <sup>1</sup>	RED, BLUE <sup>1</sup>	RED, BLUE <sup>1</sup>	WHITE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, ONYX BLACK	CANDY APPLE RED, ONYX BLACK	CANDY APPLE RED, ONYX BLACK	CANDY APPLE RED	CANDY APPLE RED, SILVER	
<b>WEIGHT OF HEAVIEST PIECE<sup>4</sup></b>	36 LBS. (REAR)	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	35 LBS. (3-WHEEL) (FRONT SECTION) 44 LBS. (4-WHEEL) (FRONT SECTION)	34 LBS. (3-WHEEL) (FRONT SECTION) 38.6 LBS. (4-WHEEL) (FRONT SECTION)	51.4 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK)	54.5 LBS. (3-WHEEL) (FRONT SECTION) 61 LBS. (4-WHEEL) (FRONT SECTION)	55 LBS. (FRONT SECTION)	77.6 LBS. (FRONT SECTION)	61 LBS. (3-WHEEL) (FRONT SECTION) 73 LBS. (4-WHEEL) (FRONT SECTION)	83 LBS. (FRONT SECTION)	95 LBS. (FRONT SECTION)	76 LBS. (3-WHEEL) (FRONT SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)	205.5 LBS. (MAIN FRAME)	
	<b>TRAVEL MOBILITY</b>						<b>VALUE SERIES</b>			<b>MID-SIZE LUXURY</b>			<b>HEAVY-DUTY</b>		<b>FULL-SIZE LUXURY</b>

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.