

Footnotes

1. All shroud panels included.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
3. Battery weight not included.
4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
5. Battery pack is 18 lbs. with 12AH batteries.
6. Durable solid plastic.
7. Deep cycle (AGM or gel cell type recommended).
8. 18 AH battery pack weighs 28 lbs.

9. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
10. Includes pathway illumination light.
11. Battery weight may vary based on manufacturer.
12. Width of base (not including seat).
13. Does not account for anti-tip clearance.
14. Includes standard seat weight.

15. Option available through most providers. Not available through Pride Mobility®.
16. Optional.
17. At midframe.
18. Light package.

PRIDE®

MOBILITY
SCOOTERS

END-TO-END



	Go Chair®	Go-Go® Ultra X	Go-Go Elite Traveller®	Go-Go® Sport ¹⁰	Go-Go® LX with CTS Suspension ¹⁰	Jazzy® Zero Turn	Go-Go® Folding Scooter ¹⁰	Revo® 2.0 ¹⁰	Jazzy® Passport	Victory® 10.2 ¹⁰	Victory® Twin ¹⁸	Victory® LX with CTS Suspension ¹⁸	Victory® LX Sport ¹⁸	Pursuit® ES ¹⁸	Maxima	Wrangler™	Pursuit® XL ¹⁸	Raptor® ¹⁸
FRONT TIRES	3" x 6" SOLID	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 7" SOLID	1.6" x 7" SOLID	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	8" solid (casters)	3.5" x 10" SOLID	DUAL 3" x 9" SOLID	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	4" x 13" PNEUMATIC	3" x 10" SOLID	14.5" TUBELESS PNEUMATIC	13" PNEUMATIC	3" x 16" PNEUMATIC
REAR TIRES	3" x 9" SOLID	2.5" x 8" SOLID	2.5" x 8" SOLID	3" x 9" SOLID	2.5" x 8" SOLID	2.5" x 9" SOLID	2" x 8" SOLID	3" x 9" SOLID	12" pneumatic drive	3.5" x 10" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	4" x 13" PNEUMATIC	4" x 10" SOLID	14.5" TUBELESS PNEUMATIC	13" PNEUMATIC	3" x 16" PNEUMATIC
MAXIMUM SPEED²	UP TO 3.7 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4.7 MPH	UP TO 4.4 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 5 MPH	UP TO 3.6 MPH	UP TO 5.2 MPH	UP TO 5.25 MPH	UP TO 5.6 MPH	UP TO 8 MPH	UP TO 6 MPH HS: UP TO 9.5 MPH	UP TO 5.3 MPH	UP TO 11.4 MPH @ 200 LBS. UP TO 11 MPH @ 350 LBS.	UP TO 7.5 MPH	UP TO 14 MPH
GROUND CLEARANCE^{4,13,17}	1.8" at motor 1.7" at anti-tips 2.7" at footplate bracket	1.3" at motor	1.38" at transaxle	3" at midframe 1.5" at motor	1.5" at motor	1.4" at motor	2.3"	3W: 3.25" 4W: 3.5" at mid-frame	3" at motor	2.4" at transaxle 2.1" at anti-tip 3.0" at side deck	3" at mid-frame	2.4" at transaxle	2.2" at motor 3.15" at mid-frame	3.5"	4.5" 3W 3.25" 4W	5.5" at motor	5.25" at mid-frame	3.5" at mid-frame
3-W TURNING RADIUS⁴	N/A	33"	33"	36"	33"	N/A	N/A	40"	N/A	42"	N/A	N/A	N/A	N/A	50.5"	N/A	N/A	60"
4-W TURNING RADIUS⁴	25.4" WITH FOOT PLATFORM	44"	44"	50.5"	44.75"	38.25"	32"	52"	24.25"	54.25"	42"	60"	61.25"	68.5"	54"	94.75"	82.5"	N/A
3-W LENGTH/WIDTH^{4,12}	N/A	37"/19.5"	37"/19.5"	42"/21.25"	39.75"/20.5"	N/A	N/A	46"/20.75"	N/A	45.25"/22.25"	N/A	N/A	N/A	N/A	46"/25.5"	N/A	N/A	59"/24.75"
4-W LENGTH/WIDTH^{4,12}	33"/25.5"	39.5"/19.5"	39.5"/19.5"	43"/21.25"	41"/20.5"	45.75"/20"	37"/19"	47"/20.75"	40"/23.5"	47"/22.25"	46"/26"	47"/22.25"	47"/25.62"	53.5"/25"	46.9"/25.5"	61.25"/32.5"	56"/27.25"	N/A
SEAT-TO-GROUND HEIGHT RANGE	22"-23"	18.25"-20.25"	18.13"-20.13"	20.5"-22"	20.75"-22.25"	21.5"-23.5"	20.5"	21"-23.5"	19"	22"-25"	23"-24.5"	23.5"-25.5"	23"-25"	24.75"-26.25"	22.5"-24.5"	30"-31.5"	26.25"-27.75"	28"
SEAT TYPE (STANDARD)	COMPACT, FOLDABLE (18"x17") 5-POINT SWIVEL WITH LATCH RELEASE	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (17"x17")	COMPACT FOLDABLE, MOLDED PLASTIC (17"x17")	COMPACT, FOLDABLE (16"x13")	COMPACT, FOLDABLE (20"x18")	GRAY FOAM (17.5" X 16.5")	FOLDABLE, MOLDED PLASTIC (18"x17")	FOLDABLE, MOLDED PLASTIC (18"x17")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	DELUXE CONTOURED MEDIUM-BACK (22"x18")	LIMITED RECLINING H/B PILLOW TOP (18"x16.5")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	HIGH BACK (20" X 18") (100° LIMITED RECLINE)
(3/4-WHEEL) WEIGHT^{3,4,14}	67 LBS.	76 LBS./80 LBS.	76 LBS./80 LBS.	93 LBS./99 LBS.	85 LBS./87 LBS.	120 LBS.	51 LBS.	120 LBS./132 LBS.	N/A/54 LBS. (w/ seat)	130.25 LBS./139 LBS.	N/A/150 LBS.	N/A/179 LBS.	N/A/160 LBS.	N/A/161 LBS.	171 LBS./176 LBS.	274 LBS.	259.9 LBS.	208 LBS.
SEAT-TO-DECK HEIGHT RANGE	16.25-17.25"	15.75"-17.25"	15.75"-17.25"	17"-18"	16.25"-17.25"	16"-18"	13"	15"-18"	14.5" (no adjustment)	15.25" - 18.25"	17.5"-19"	16.5"-18.5"	16.5"-18.5"	17.75"-18.75"	16"-18"	18.75"-20.25"	21"-22.5"	19"
BATTERY WEIGHT¹¹	14 LBS. EA.	9.3 LBS. EA. ⁵	13.3 LBS. EA. ⁸	13.3 LBS. EA. ⁸	13.3 LBS. EA. ⁸	13.5 LBS. EA.	9 LBS. EA. (SEALED LEAD ACID)/5 LBS. EA. (LITHIUM-ION)	25 LBS. EA.	7 LBS.	24 LBS. EA./32 LBS. EA.	24 LBS. EA./32 LBS. EA.	24 LBS. EA./32 LBS. EA.	30.2 LBS. EA.	53 LBS. EA.	38 LBS. EA. (NF-22 (55 AH))/24.5 LBS. EA. (U-1)	48 LBS. EA. (75 AH)/60 LBS. EA. (100 AH)	70 LBS. EA. (G31)/60 LBS. EA. (G27)/53 LBS. EA. (G24)	53 LBS. EA.
BATTERY REQUIREMENTS⁷	(2)18AH (INCLUDED)	(2) 12AH INTERCEPTOR (INCLUDED)	(2) 22AH INTERCEPTOR (INCLUDED)	(2) 22AH INTERCEPTOR (INCLUDED)	(2) 22AH INTERCEPTOR (INCLUDED)	(2) 22AH (INCLUDED)	(2) 12AH SEALED LEAD ACID (STANDARD) OR (1) 11.5 AH LITHIUM ION ¹⁶	(2) U-1 (INCLUDED)	(1) 18 AH LITHIUM-ION (INCLUDED) OR 10.4 AH LITHIUM ION TRAVEL BATTERY ¹⁴	(2) U-1 (INCLUDED) OR (2) 40AH ¹⁵	(2) U-1 (INCLUDED) OR (2) 40AH ¹⁵	(2) U-1 (INCLUDED) OR (2) 40AH ¹⁶	(2) 40AH (INCLUDED)	(2) GROUP 24 (INCLUDED)	(2) NF-22 INTERCEPTOR (INCLUDED) OR (2) U-1 (POWER ELEVATING SEAT) ¹⁶	(2) 75AH (STANDARD) (2) 100 AH (OPTIONAL)	(2) GROUP 31 ¹⁵ , (1) 10AH) (2) GROUP 27 (100AH) OR (2) GROUP 24 (75 AH)	(2) 75AH (INCLUDED)
PER CHARGE RANGE (MAX.)^{2,9}	UP TO 8.7 MILES (10.9 MILES WITH EXTENDED RANGE BATTERIES)	3W: UP TO 6.9 MILES 4W: UP TO 7.2 MILES	3W: UP TO 13.9 MILES 4W: UP TO 13.9 MILES	3W: UP TO 15.9 MILES 4W: UP TO 16.2 MILES	3W: UP TO 11.7 MILES 4W: UP TO 14.6 MILES	UP TO 10.6 MILES	4W: UP TO 13 MILES (LITHIUM-ION) 4W: UP TO 9.3 MILES (SEALED LEAD ACID)	3W: UP TO 12 MILES 4W: UP TO 13 MILES	UP TO 9.6 MILES	3W/4W: UP TO 11 MILES (U-1) 3W: UP TO 13 MILES (40AH) 4W: UP TO 13.5 MILES (40 AH)	UP TO 11.1 MILES (WITH U-1 BATTERIES)	UP TO 16.16 MILES	UP TO 19.7 MILES (200 LBS.) UP TO 13 MILES (400 LBS.)	UP TO 33.6 MILES HS: UP TO 22 MILES	3W: UP TO 18 MILES, U-1 (POWER ELEVATING SEAT) UP TO 9.6 MILES 4W: UP TO 16 MILES	UP TO 18.6 MILES @ 350 LBS. (75AH) UP TO 24.2 MILES AT 350 LBS. (100AH) UP TO 22.6 MILES @ 200 LBS. (75AH) UP TO 29.4 MILES AT 200 LBS. (100AH)	G31: UP TO 31.83 MILES G27: UP TO 29.51 MILES G24: UP TO 22.13 MILES	19 MILES (UP TO 400 LBS. AND FULL SPEED) 31 MILES (UP TO 220 LBS. AND AT LOW SPEED)
BATTERY CHARGER	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 3.5A	OFF-BOARD, 3.5A	OFF-BOARD, 3.5A	OFF-BOARD, 3.5A	OFF-BOARD, 2A (SEALED LEAD ACID) OFF-BOARD, 2.5A (LITHIUM-ION)	OFF-BOARD, 3.5A	OFF-BOARD, 2A	OFF-BOARD, 3.5A (U-1) OFF-BOARD, 5A (40 AH)	OFF-BOARD, 3A (U-1) OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 8A	OFF-BOARD, 5A	OFF-BOARD, 8A	OFF-BOARD, 8A	OFF-BOARD, 8A
SUSPENSION	LIMITED	N/A	N/A	N/A	FRONT & REAR	FRONT & REAR	N/A	FRONT & REAR	LIMITED	N/A	N/A	FRONT & REAR	FRONT & REAR	FRONT & REAR	3-W REAR/ 4-W FRONT & REAR	FRONT & REAR	FRONT & REAR	FRONT & REAR
MAXIMUM WEIGHT CAPACITY	300 LBS.	260 LBS.	300 LBS.	325 LBS.	300 LBS.	325 LBS.	250 LBS.	375 LBS.	250 LBS.	400 LBS.	350 LBS.	400 LBS.	400 LBS.	350 LBS.	500 LBS.	350 LBS.	400 LBS.	400 LBS.
STANDARD BODY COLORS	RUBY RED, SAPPHIRE BLUE, ROSE QUARTZ, CITRINE YELLOW, AMBER ORANGE, PEARL WHITE	RED ⁶	RED, BLUE ¹	RED, BLUE ¹	RED, BLUE ¹	FIRE OPAL, GREEN MOONSTONE, ONYX BLACK, PEWTER, PEARL WHITE, SAPPHIRE BLUE	WHITE	GREY STREET, TRUE BLUE	SILVER	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED	DESERT SAND, STONE GREY	CANDY APPLE RED, SILVER	SILVER
WEIGHT OF HEAVIEST PIECE⁴	36 LBS. (REAR)	29 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	29 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	35 LBS. (3-WHEEL) (FRONT SECTION) 44 LBS. (4-WHEEL) (FRONT SECTION)	34 LBS. (3-WHEEL) (FRONT SECTION) 37 LBS. (4-WHEEL) (FRONT SECTION)	48 LBS. (REAR SECTION)	51 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK)	48 LBS. (3-WHEEL) (REAR SECTION) 57 LBS. (4-WHEEL) (FRONT SECTION)	N/A	56.5 LBS. (3-WHEEL) (FRONT SECTION) 64.5 LBS. (4-WHEEL) (FRONT SECTION)	55 LBS. (FRONT SECTION)	77.6 LBS. (FRONT SECTION)	77 LBS. (FRONT SECTION)	118 LBS. (MAIN FRAME)	76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)	219 LBS. (UNIT DOES NOT DISASSEMBLE)	205.5 LBS. (MAIN FRAME)	N/A (UNIT DOES NOT DISASSEMBLE)

TRAVEL MOBILITY

VALUE SERIES

MID-SIZE LUXURY

HEAVY-DUTY

FULL-SIZE LUXURY