

NEW From

mo-vis
move with vision



IDM-ARLITE-9

- Mo-Vis All-Round Light is a compact version of the standard joystick module but with less required force (120gr).
- Fully proportional with reduced amplitude and power requirement •Built in a compact and discreet housing with the possibility of 1 or 2 satellite switches.
- Switches for on/off function or mode button.
- All-round use: can be used as a chin joystick or standard joystick.
- Free downloadable software to configure joystick via PC (equipped with a mini USB port)

MSRP - \$2365.00 (E2373)



IDM-HD

- Fully proportional with increased amplitude and high power requirement (650gr)
- Mo-Vis All-Round Heavy Duty Joystick with adapter for Gatlin Mount System
- All-Round Heavy Duty use: can be used as joystick for hand or foot control
- Switches for on/off function or mode button.
- Fully Adjustable to individual possibilities and needs.
- Easily configurable with the free Mo-Vis Configuration software.

MSRP - \$3280.00 (E2373)



IDM-HW

- The Mo-Vis Hand warmer improves the flexibility of cold body parts.
- Always warm hand with optimal finger mobility and functionality
- Heated air is recycled
- Reduced power consumption and more power efficiency.
- A twister switch can be plugged in for on/off or power level selection.
- Auto off, power levels, temperature and fan speed are fully programmable.

MSRP - \$1,180.00



IDM-DOME

- To better maintain the air temperature around the hands.
- MSRP - \$525.00

