



THE MOST COMFORTABLE RECLINER IN THE WORLD!





COLLECTION	ESCAPE		ATLAS PLUS	ELEGANCE		URBANA	LEGACY		PERFECTA	TRANQUIL	TRANQUIL	TRANQUIL	TRANQUIL	METRO
MODEL NUMBER	PLR990iM	PLR990iL	PLR2985M	PLR975M	PLR975L	PLR965M	PLR958M	PLR958L	PLR945M	PLR935PW	PLR935S	PLR935M	PLR935LT	PLR925M
POSITION	Infinite Position Trendelenburg Chaise Lounger Relaxed Reclined Upward Lift	Infinite Position Trendelenburg Chaise Lounger Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift	Infinite Position Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
SEAT-TO-FLOOR	20"	20"	19.5"	20"	20.5"	20"	19.5"	20.5"	19"	18"	18.5"	19.5"	21.5"	19"
SEAT DEPTH	20"	22"	21"	21"	22"	20"	21.5"	22"	20"	18.5"	18.5"	20.5"	22.5"	21"
SEAT WIDTH	20"	23"	22"	21"	22"	21"	22"	23"	22"	24"	20"	22"	23.5"	21"
TOP OF BACK TO SEAT	26"	29"	26"	27"	28"	25"	25"	26"	25"	24.5"	25"	29"	30"	26.5"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard
SEATING TYPE	Pocketed Coil Springs	Pocketed Coil Springs	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	44.5"	47"	42"	42"	43"	42"	42"	43.5"	43"	38"	40"	42"	47"	41"
DISTANCE FROM WALL	18"	21"	16"	14"	16"	14"	16"	16"	15"	12.5"	13"	16"	16"	14"
OVERALL WIDTH	33"	36.5"	39"	35.5"	37.5"	33"	35"	36.5"	31"	36"	33.75"	37.5"	38"	33"
CHAIR WEIGHT	148 lbs.	150 lbs.	164 lbs.	152 lbs.	157 lbs.	145 lbs.	153 lbs.	160 lbs.	142 lbs.	143 lbs.	137 lbs.	145 lbs.	155 lbs.	142 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated


STANDARD FABRIC SELECTION GUIDE	<u>DURASOFT</u>		<u>BADLANDS</u>	<u>BADLANDS</u>	<u>STONEWASH</u>	<u>SAVILLE</u>	<u>MERINO</u>	<u>ASTRO</u>	<u>SAVILLE</u>
									
									
									
									
									
	BUFF	FUDGE							



POWER HEADREST


POWER LUMBAR


USB REMOTE



LITHIUM BATTERY BACKUP



FOOTREST EXTENSION



POWER BACKREST


 = Indicates shortest available lead time.


POSITION GUIDE*:

 RELAXED

 RECLINED

 UPWARD LIFT

 INFINITE POSITION

 TRENDELENBURG

*Not all units achieve all positions. Please refer to specific product.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



POWER HEADREST

1.



POWER LUMBAR



FOOTREST EXTENSION

4.



INFINITE POSITION

6.



2.

UPWARD LIFT







BACKREST

5.



THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

INDEPENDENT FUNCTIONS




-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.

STANDARD FEATURES

-  7. The USB remote adjusts all functions independently and allows you to charge your smart device
-  8. Lithium battery backup ensures your recliner is operational when the power goes out.



USB REMOTE



PATENT PENDING HEADREST/LUMBAR MECHANISM WITH MOTORS

Independently operates the headrest and lumbar functions for enhanced head, neck and back support. Anti-entrapment on headrest for safety.

STRONG WOOD FRAME

Premium furniture-grade hardwood is robust, durable and long-lasting.



STEEL SEAT BOX

Springs are attached to heavy-duty steel frame.

TOTAL COMFORT SEATING

Sinuous springs with plastic-lined clips are attached to morley hinge.

LIFT AND RECLINE MOTORS

Independently adjusts the recliner to achieve relaxed, reclined, lay flat or upward lift positioning.

STEEL FRAME

Heavy-duty metal framework forged to last.

PRIDE
MOBILITY
LIVE YOUR BEST®

The information contained herein is correct at time of publication.
We reserve the right to change specifications without prior notice.

VivaLIFT!
POWER RECLINER COLLECTION



Beamsville, Ontario | 888-570-1113

pridemobility.ca | prideliftchairs.com |    