

END-TO-END POWER CHAIR SOLUTIONS





























| | | | a Go Ta | Or Just 4 | | OF THE CO | | | | | | | | W Table | |
|---|--|--|--|--|--|--|--|--|--|---|--|--|--|--|--|
| | JAZZY AIR® MED¹ | JAZZY AIR® 2 ¹ | JAZZY® ULTRA LIGHT¹ | JAZZY® CARBON¹ | GO CHAIR®1 | GO CHAIR® MED ¹ | JAZZY® EVO 613 SERIES¹ | JAZZY® EVO 614¹ | JAZZY® EVO 614HD¹ | JAZZY® ELITE ES¹ | JAZZY® SELECT ¹ | JAZZY® ELITE 14 ¹ | JAZZY® 600 ES ¹ | JAZZY® ELITE HD¹ | JAZZY® 1450¹ |
| Model Number | JAZZY AIR MED 2SEC/2SE-SS | JAZZYAIR2 | JZULIGHTMS | JAZZY CARBON | GO-CHAIR | GO CHAIR MED GOCHAIRMED | JAZZY EVO 613 2S-C/2S-SS JAZZY EVO 613 Li 2S-C/2S-SS | JAZZY EVO 614 2S-C/2S-SS | JAZZY EVO 614HD 2SHD-C/2SHD-SS | JAZZY ELITE ES JELITEES 2S-C/2S-SS | JAZZY SELECT 2S-C/2S-SS | JAZZY ELITE 14 JELITE14 2S-C/2S-SS | JAZZY 600 ES 2S-C/2S-SS | JAZZY ELITE HD JELITEHD 2HD-C/2HD-SS | JAZZY 1450 J1450 2VHD-0 2VHD-SS/2SPVHD-SS |
| Weight Capacity | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 450 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 300 lbs. | 450 lbs.7 | 600 lbs. |
| Speed (Maximum) ³ | Up to 3.7 mph (elevated and non-elevated) | Up to 4 mph (elevated and non-elevated) | Up to 3.7 mph | Up to 3.7 mph | Up to 3.8 mph | Up to 3.7 mph | Up to 5.1 mph | Up to 4.4 mph | Up to 4.88 mph | Up to 4 mph | Up to 4.3 mph | Up to 4 mph | Up to 4 mph | Up to 4 mph | Up to 4.9 mph |
| Ground Clearance ¹¹ | 1" at motor bracket | 1" at motor bracket | 2.25" at motor | 2" | 1.75" at motor | 1.75" at motor | 3" at motor bracket | 3" at motor bracket | 3" at motor bracket | 1.85" at motors | 0.75" at motor bracket | 3" at battery tray | 4" at motor | 1.5" at anti-tip wheels | 1.5" at anti-tip wheels |
| Turning Radius ¹¹ | 24.5" | 24.25" | 33.25" | 24" | 25.5" | 25.5" | 22" | 22" | 22" | 24.75" | 22.5" | 24" | 20.5" | 24" | 26.5" |
| Overall Length ^{11,2} | 41.5" (with footplate) | 33.25"/42" (with footplate) | 41" | 37" (with footplate) | 28.5"/ 33" (with footplate) | 28.5"/ 33" (with footplate) | 35" (with footplate) | 40.5" (with footplate) | 40.5" (with footplate) | 34.75"/41" (with footplate) | 32.99"/44" (with footplate) | 35.5"/42.5" (with footplate) | 35.5"/40.25" (with footplate) | 35.5" | 38.25" |
| Overall Width ¹¹ | 23.25" | 23.25" | 24.5" | 24" | 22.5" | 23" | 22" | 25.4" | 25.4" | 22.75" | 23.75" | 24.25" | 24.5" | 24.25" | 29.25" |
| Front Wheels | 6" solid casters | 6" solid casters | 8" solid | 7" solid casters | 6" solid casters | 6" solid casters | 5" solid | 6" solid | 6" solid | 3" solid anti-tips | 5" solid anti-tips | 3" solid anti-tips | 6" solid OMNI-casters | 3" solid anti-tips | 4" solid anti-tips |
| Drive Wheels | 10" solid | 10" solid | 10.75" solid | 8.5" solid | 9" solid | 9" solid | 13" solid | 13.75" perforated | 13.75" perforated | 9" solid | 10" solid | 14" foam filled | 14" foam filled | 14" foam filled | 14" solid |
| Rear Wheels | 6" solid casters | 6" solid casters Active-Trac® (ATX) with | N/A | N/A Front | 3" solid anti-tips Limited | 3" solid anti-tips Limited | 5" solid Active-Trac® (ATX) | 6" solid | 6" solid Full independent front | 6" solid casters | 6" solid casters Active-Trac® (ATX) | 6" solid casters Limited | 6" solid OMNI-casters | 6" solid casters | 9" solid casters |
| Suspension Type | Active-Trac® (ATX) with Mid-Wheel 6® Technology | Mid-Wheel 6® Technology | Suspension Caster Forks | Up to 9.3 miles at 300 lbs. | | | EVO 613 Li: Up to 21 miles at 200 lbs. | Full independent front & rear Active-Trac® (ATX) Up to 21.5 miles at 200 lbs. | & rear Active-Trac® (ATX) Up to 21.5 miles at 200 lbs. | Up to 15 miles at 300 lbs. | , | | Active-Trac® (ATX) with Mid-wheel 6® technology | | Swing Arm Up to 16.5 miles at a 600 lbs |
| Range per Charge ^{3,6} | Up to 14.1 miles at 220 lbs. Up to 13 miles at 300 lbs. | Up to 16.6 miles at 220 lbs. Up to 15 miles at 300 lbs. | Up to 7 miles at 200 lbs. (10Ah) Up to 16.4 miles at 200 lbs. (20Ah) Up to 14.4 miles at 300 lbs. (20Ah) | | Up to 10.9 miles at 200 lbs. Up to 8.8 miles at 300 lbs. | Up to 12.1 miles at 200 lbs. Up to 9.8 miles at 300 lbs. | Up to 17.6 miles at 300 lbs. EVO 613 : Up to 17.1 miles at 200 lbs. Up to 14 miles at 300 lbs. | Up to 19 miles at 300 lbs. | Up to 13 miles at 450 lbs. | · | Up to 19.4 miles at 300 lbs. | Up to 18 miles at 300 lbs. | Up to 29.5 miles at 200 lbs. Up to 24.8 miles at 300 lbs. | Up to 14.4 miles at 450 lbs. | |
| Weight of Heaviest Piece ⁵ | 197.2 lbs. (with seat) ⁹ | 202 lbs. (with seat) 9 | 32.6 lbs. (unit doesn't disassemble) ¹⁰ | 39.6 lbs. (unit doesn't disassemble)10 | 36 lbs. | 36 lbs. | 94 lbs. (base) | 125.7 lbs. (base) | 125.7 lbs. (base) | 69 lbs. | 94 lbs. | 99 lbs. | 124 lbs. | 99 lbs. | 195 lbs. |
| Drivetrain | Two-motor, Mid-Wheel 6 [®] drive | Two-motor, Mid-Wheel 6 [®] drive | Two motor, rear-wheel drive | 24v/250W brushless motors | Two-motor, rear-wheel drive | Two-motor, rear-wheel drive | Two-motor, in-line, mid-wheel drive | Two in-line HD motors | Two in-line HD motors | Two-motor, in-line, front-wheel drive | Two-motor, mid-wheel drive | Two-motor, front-wheel drive | Two-motor, mid-wheel drive | Two-motor, front-wheel drive | Two-motor, front-wheel drive |
| Brakes | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical | Regenerative and electro-mechanical |
| Batteries | (2) 12 volt, U-1 (std.) | (2) 12 volt, U-1 (std.), 40Ah (optional) | 25.2V, 10Ah Lithium-ion 25.2V, 20Ah Lithium-ion | (1) 24V 12Ah (288WH) lithium ion | (2) 20Ah, 12 volt, nano-silica | (2) 20Ah, 12 volt, nano-silica | EVO 613 Li: (1) 25.6V, 32.4Ah lithium iron phosphate EVO 613: (2) 12 volt, U-1 | (2) NF-22 | (2) NF-22 | (2) 12 volt, U-1 | (2) 12 volt, U-1 | (2) 12 volt, NF-22 | (2) 12 volt, NF-22 | (2) 12 volt, NF-22 | (2) 12 volt, Group 24 |
| Battery Weight ⁸ | 22.2 lbs. each | 23.5 lbs. (U-1) 32.5 lbs. (40Ah) | 5 lbs. (10Ah) 7.5 lbs. (20Ah) | 4 lbs. | 16 lbs. each (20Ah) | 16 lbs. each (20Ah) | EVO 613 Li: 18.5 lbs. EVO 613: 23.5 lbs. each | 36 lbs. each | 36 lbs. each | 24.5 lbs. each | 25.2 lbs. each | 37 lbs. each | 37 lbs. each | 37 lbs. each | 52 lbs. each |
| Battery Charger | 5A, off-board | 5A, off-board | 3.5A, off-board | 2A, off-board | 3.5A, off-board | 3.5A, off-board | 3.5A, off-board | 5A, off-board | 5A, off-board | 3.5A, off-board | 3.5A, off-board | 5A, off-board | 5A, off-board | 5A, off-board | 8A, off-board |
| Standard Seating | 45 lbs. memory high back | 41 lbs. sport, high-back | N/A | N/A | 30 lbs. (18") 32 lbs. (20") compact, foldable | 18" wide: 29 lbs.4 20" wide: 32 lbs.4 compact, foldable | 46 lbs. (16") 48 lbs. (18") 50 lbs. (20") memory high-back | 57.8 lbs. (20") memory high-back | 57.8 lbs. (20") memory high-back | 50.2 lbs comfort, high-back | 50.2 lbs. comfort, high-back | 18" wide: 51.8 lbs. 20" wide: 54 lbs. 22" wide: 57.6 lbs. high-back | 18" wide: 51.8 lbs. 20" wide: 54 lbs. high-back | 18" wide: 51.8 lbs. 20" wide: 54 lbs. 22" wide: 57.6 lbs. high-back | 52 lbs. high-back |
| Standard Seat Size | 16" x 16-18", 18" x 18-20", 20" x 18-20" | 16" x 16-18", 18" x 18-20", 20" x 20-22" | 18" x 18" | 17" x 17" | 18" x 17" 20" x 18" | 18" x 17", 20" x 18" | 16" x 16-18", 18" x 18-20", 20" x 18-20" | 16" x 16-18", 18" x 18-20", 20" x 18-20" | 18" x 18-20", 20" x 18-20", 22" x 20-22", 24" x 22-24" (option) | 16" x 16-18", 18" x 18-20", 20" x 18-20" | 16" x 16-18", 18" x 18-20", 20" x 18-20" | 18" x 18-20", 20" x 18-20", 22" x 20-22" | 18" x 18-20", 20" x 18-20" | 18" x 18-20", 20" x 18-20", 22" x 20-22", 24" x 24" | 24" x 24" |
| High-Back Seat (max. size) | 20"W x 18-20"D | 20"W x 20-22"D | N/A | N/A | N/A | N/A | 20"W x 18-20"D | 20"W x 18-20"D | 24"W x 22-24"D | 20"W x 20"D | 20"W x 20"D | 22"W x 22"D | 20"W x 20"D | 22"W x 22"D | 32"W x 28"D |
| Seat to Ground Height | 21.06"-32.28" | 21.25"-33" | 21.26" | 18.75/20.25"12 | N/A | N/A | 19.75-21.25" | 22-24" | 22-24" | 21.5-23.5" | 21-22" | 20-22" | 20.5-22.5" | 20-22" | 20.5-22.5" |
| Synergy [®] Seating | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | 20"W x 20"D | 20"W x 20"D | 20"W x 20"D | 24"W x 24"D | 32"W x 28"D |
| Synergy [®] Seat to Ground Height | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | 17.25-19.25" | 15.75-17.75" | 15.875-17.875" | 15.75-17.75" | 17.25-19.25" |
| Power Elevating Seat | Standard | Standard | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A |
| Standard Electronics | 40A, PG VR2 | 40A, PG VR2 | 40A, PG nVR2 | 35A, Brushless Micon | 40A, PG GC3 | 40A, PG GC3 | 40A, PG GC3 | 60A Dynamic LiNX | 60A Dynamic LiNX | 40A, PG GC3 | 40A, PG GC3 | 50A, Dynamic LiNX | 50A, Dynamic LiNX | 50A, Dynamic LiNX | 90A, PG VR2 |
| Specialty Controls | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | N/A | Yes |
| HCPCS Codes | K0830/K0831 Group 2: Standard, Seat Elevator | K0899 Power Mobility Device not coded by DME PDAC. | K0899 Power Mobility Device not coded by DME PDAC. | K0899 Power Mobility Device not coded by DME PDAC. | K0899 Power Mobility Device not coded by DME PDAC. | K0821 Group 2: Standard, Portable | K0822/K0823 Group 2: Standard | K0822/K0823 Group 2: Standard | K0824/K0825 Group 2: Heavy Duty | K0822/K0823 Group 2: Standard | K0822/K0823 Group 2: Standard | K0822/K0823 Group 2: Standard | K0822/K0823 Group 2: Standard | K0824/K0825 Group 2: Heavy Duty | K0826/K0827 Group 2: Very Heavy Duty K0839 Single Power Very Heavy Dut |
| Warranty | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | 5-year limited on frame; ; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | Lifetime limited on frame; 1-year limited on electrical control/lever system, motor/drive system, bearing sleeve and battery | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month limited on SLA batteries; 3-year limited on lithium batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month limited on batteries | 5-year limited on frame; 13-month limited on electronics 13-month on drive motors; 13-month on batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | 5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries | K0826/K0827 Models: Lifetime limited on frame; 13-month limited on drive mote 13-month limited on drive mote 13-month on batteries K0839 Model: Lifetime limited on frame; 2-year limited on electronics; 18-month limited |

^{1.} Pride FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

BATTERY HCPCS CODES

U-1 E2365 NF-22 E2361 Group 24 E2363 Lithium E2397 18 AH K0733 HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.





Without front riggings.

^{3.} Varies with user weight, terrain, battery type, battery condition and tire pressure.

Seat weight includes joystick.
 Does not include seat or battery weights.

Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.

^{7.} Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18" deep seat configurations.

Battery weight may vary based on manufacturer.
 Seat weight included (because tools are required to remove seat). Battery weight not included.

 ^{10.} Seat weight included. Unit does not disassemble.
 11. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3 %. Dimensions listed are for power base only. Overall

measurements will vary based on seating and accessory selections.

12. With and without standard cushion.