



	JAZZY AIR® 2¹	GO CHAIR®¹	JAZZY® ULTRA LIGHT¹	JAZZY® CARBON¹	JAZZY® CARBON 27¹	JAZZY® CARBON 27X¹	JAZZY® CARBON HD¹
Model Number	JAZZYAIR2	GO-CHAIR	JULTRALIGHT	JAZZY CARBON	JAZZY CARBON 27	JAZZY CARBON 27X	JAZZY CARBON HD
Weight Capacity	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	275 lbs.	400 lbs.
Speed (Maximum)²	Up to 4 mph (elevated and non-elevated)	Up to 3.8 mph	Up to 3.7 mph	Up to 3.7 mph	Up to 3.7 mph	Up to 3.7 mph	Up to 3.7 mph
Ground Clearance⁸	1" at motor bracket	1.75" at motor	2.25" at motor	2"	1" at motor	3"	1" at motor
Turning Radius⁸	24.25"	25.5"	33.25"	24"	33.4"	35.4"	36.4"
Overall Length⁸	33.25"/42" (with footplate)	28.5"/33" (with footplate)	31"/41" (with footplate)	28.3"/37" (with footplate)	27.5"/35.04" (with footplate)	34"/41" (with footplate)	31.3"/36.81" (with footplate)
Overall Width⁸	23.25"	22.5"	24.5"	24"	23.43"	23"	26.4"
Front Wheels	6" solid casters	6" solid casters	8" solid casters	7" solid casters	6.9" solid casters	7" solid casters	6.5" solid casters
Drive Wheels	10" solid	9" solid	10.75" solid	8.5" solid	8.3" solid	10" solid	10.6" solid
Rear Wheels	6" solid casters	3" solid anti-tips	2" solid anti-tips	2" solid anti-tips	2" solid anti-tips	2" solid anti-tips	2" solid anti-tips
Suspension Type	Active-Trac® (ATX) w/Mid-Wheel 6® Technology	Limited	Suspension Caster Forks	Full Front and Rear Suspension	Full Front Suspension	N/A	Full Front and Rear Suspension
Range per Charge²,⁴	Up to 16.6 miles at 220 lbs. Up to 15 miles at 300 lbs.	Up to 10.9 miles at 200 lbs. Up to 8.8 miles at 300 lbs.	Up to 8 miles at 200 lbs. (10Ah) Up to 7 miles at 300 lbs. (10Ah) Up to 16.4 miles at 200 lbs. (20Ah) Up to 14.4 miles at 300 lbs. (20Ah)	Up to 9.3 miles at 300 lbs.	Up to 9.3 miles at 220 lbs.	Up to 8.1 miles at 220 lbs.	Up to 12.4 miles at 220 lbs.
Weight of Heaviest Piece³	202 lbs. (with seat)⁵	36 lbs.	32.6 lbs. (unit doesn't disassemble)⁷	39.6 lbs. (unit doesn't disassemble)⁷	27.78 lbs. (unit doesn't disassemble)⁷	27.8 lbs. (unit doesn't disassemble)⁷	32.2 lbs. (unit doesn't disassemble)⁷
Drivetrain	Two-motor, Mid-Wheel 6® drive	Two-motor, rear-wheel drive	Two motor, rear-wheel drive	24v/250W brushless motors	24V/250W	24V/150W	24V/250W Brushless motors
Brakes	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical
Batteries	(2) 12V, U-1 (std.), 40Ah (optional)	(2) 20Ah, 12V, nano-silica	25.2V, 10Ah Lithium-ion 25.2V, 20Ah Lithium-ion	(1) 24V, 12Ah (288WH) Lithium-ion	(1) 24V, 12Ah, Lithium-ion	(1) 24V, 10Ah, Lithium-ion	(1) 24V, 20Ah, Lithium-ion
Battery Weight⁵	23.5 lbs. (U-1) 32.5 lbs. (40Ah)	16 lbs. each (20Ah)	5 lbs. (10Ah) 7.5 lbs. (20Ah)	4 lbs.	4.6 lbs.	3.5 lbs.	7.5 lbs.
Battery Charger	5A, off-board	3.5A, off-board	3.5A, off-board	2A, off-board	2A, off-board	3.5A, off-board	2A, off-board
Standard Seating Weight	41 lbs. sport, high-back	30 lbs. (18") 32 lbs. (20") compact, foldable	N/A	N/A	N/A	N/A	N/A
Standard Seat Size	16" x 16-18", 18" x 18-20", 20" x 20-22"	18" x 17" 20" x 18"	18" x 18"	17" x 17"	17" x 17"	18" x 18"	20" x 20"
High-Back Seat (max. size)	20"W x 20-22"D	N/A	N/A	N/A	N/A	N/A	N/A
Seat to Ground Height	21.25"-33"	N/A	21.26"	18.75/20.25"⁹	19.1"	19.6"	19.7"
Power Elevating Seat	Standard	N/A	N/A	N/A	N/A	N/A	N/A
Standard Electronics	40A, PG VR2	40A, PG GC3	40A, PG nVR2	35A, Brushless	35A, Brushless	35A, Brushless	35A, Brushless
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	Lifetime limited on frame; 1-year limited on electrical control/lever system, motor/drive system, bearing sleeve and battery	5-year limited on frame; 1-year limited on electronics; 1-year limited on drive motors; 1-year limited on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month on batteries	5-year limited on frame; 1-year limited on electronics; 1-year limited on drive motors; 1-year limited on batteries

1. Pride FDA Class II Medical Devices are designed to aid individuals with mobility impairments.  
2. Varies with user weight, terrain, battery type, battery condition and tire pressure.  
3. Does not include seat or battery weights.  
4. Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities.This specification can be subject to a variance of +10%, -5%.  
5. Battery weight may vary based on manufacturer.

6. Seat weight included (because tools are required to remove seat). Battery weight not included.  
7. Seat weight included. Unit does not disassemble.  
8. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3 %.  
Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.  
9. With and without standard cushion.