

Footnotes

- All shroud panels included.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of + or - 10%.
- Battery weight not included.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- Battery pack is 18 lbs. with 12AH batteries.
- Deep cycle (AGM or gel cell type recommended).
- 18 AH battery pack weighs 28 lbs.
- Tested in accordance with ANSIRENSA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Battery weight may vary based on manufacturer.
- Width of base (not including seat).
- Does not account for anti-tip clearance.
- Includes standard seat weight.
- Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.
- At midframe.
- Light package.

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST® **END-TO-END**



	Jazzy® Zero Turn	Revo™ 2.0	Jazzy® Passport	Go Chair®	Go-Go® ES 2	Go-Go® Ultra X	Go-Go Elite Traveller®	Go-Go Elite Traveller® Plus	Go-Go® Sport	Go-Go® LX w/ CTS Suspension	Go-Go® Folding Scooter	Victory® 9	Victory® 10	Victory® 9.2	Victory® 10.2	Victory 10® LX w/ CTS Suspension ¹⁵	Victory® LX Sport ¹⁵	Maxima	Pursuit® ¹⁵	Pursuit® XL ¹⁵	Raptor ¹⁵	
PART NUMBER	JAZZYZTS	S66 (3W) S67 (4W)	FJPT1700	Go Chair	S81	SC40X (3W) SC44X (4W)	SC40E (3W) SC44E (4W)	SC53 (3W) SC54 (4W)	S73 (3W) S74 (4W)	SC50LX (3W) SC54LX (4W)	S19WH1001 (SLA) S19WH1002 (Lithium Ion)	SC609 (3W) SC709 (4W)	SC610 (3W) SC710 (4W)	S6092 (3W) S7092 (4W)	S6102 (3W) S7102 (4W)	S710LX	S710LXW	SC900 (3W) SC940 (4W)	SC713	SC714	R3-1700	
FRONT TIRES	2" x 7" SOLID	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	8" SOLID (CASTER)	6" SOLID (CASTER)	2" x 7.5" SOLID	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	1.5" x 7" SOLID	3W: 3" x 8" SOLID 4W: 3" x 9" SOLID	3.5" x 10" SOLID	3W: 3" x 8" SOLID 4W: 3" x 9" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3" x 10" SOLID	3" x 13" SOLID	13" PNEUMATIC	3" x 16" PNEUMATIC	
REAR TIRES	2.5" x 9" SOLID	3" x 9" SOLID	12" (PNEUMATIC DRIVE)	9" SOLID (DRIVE)	3" x 7.5" SOLID	2.5" x 8" SOLID	2.5" x 8" SOLID	3" x 9" SOLID	3" x 9" SOLID	2.5" x 8" SOLID	2" x 8" SOLID	3" x 9" SOLID	3.5" x 10" SOLID	3" x 9" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	4" x 10" SOLID	13" SOLID	13" PNEUMATIC	3" x 16" PNEUMATIC	
MAXIMUM SPEED²	UP TO 4 MPH	UP TO 5 MPH	UP TO 3.6 MPH	UP TO 3.7 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4.5 MPH	UP TO 4.7 MPH	UP TO 4.4 MPH	UP TO 4 MPH	UP TO 5.3 MPH	UP TO 5.3 MPH	UP TO 5.6 MPH	UP TO 5.2 MPH	UP TO 5.6 MPH	UP TO 8 MPH	UP TO 5.3 MPH	UP TO 8 MPH	UP TO 7.5 MPH	UP TO 14 MPH	
GROUND CLEARANCE^{4,11,14}	1.4" at motor 1.5" at anti-tip 3.5" at side deck	1.25" at transaxle 3.25" at mid-frame (3W) 3.5" at mid-frame (4W)	3" at motor	1.8" at motor 1.7" at anti-tip 2.7" at foot plate bracket	1.5" at motor 0.75" at anti-tip	1.3" at transaxle 0.67" at anti-tip 3.1" at deck (3W) 3" at deck (4W)	1.38" 0.75" at anti-tip 2.25" at mid-frame	1.7" at transaxle (3W) 0.83" at anti-tip (3W) 2.6" at side deck (3W)	1.5" at motor 0.75" at anti-tip 1.25" at rear frame lock (3W) 3" at controller (4W)	1.5" at motor 0.75" at anti-tip	1.5" at motor 0.75" at anti-tip 2.5" at frame	2.3"	1.8" at transaxle 1.4" at anti-tip 2.7" at side deck (3W) 2.1" at motor bracket bolt (4W)	2.2" at motor (3W) 1.26" at anti-tip (3W) 3.4" at side deck (3W) 2" at motor (4W) 1.5" at anti-tip (4W) 3.1 at side deck (4W)	1.8" at transaxle 1.4" at anti-tip 2.7" at side deck	2.4" at transaxle 2.1" at anti-tip 3.0" at side deck	2.4" at transaxle 0.9" at anti-tip 2.9" at front frame	2.2" at motor 3.15" at rear mid-frame	2.5" at motor 2" at anti-tip (3W) 4.5" at side deck (3W) 1.5" at anti-tip (4W) 3.25" (4W)	3.7" at motor 1.8" at anti-tip 4.5" at side deck	2.75" at motor 5.25" at mid-frame	3.5" at center of frame 1.37" at anti-tip 6.18"
3-W TURNING RADIUS⁴	N/A	40"	N/A	N/A	31"	33"	33"	34.5"	36"	33"	32"	40.25"	42"	40.25"	42"	N/A	N/A	50.5"	N/A	N/A	60"	
4-W TURNING RADIUS⁴	38.25"	52"	24.25"	25.5" with foot platform	N/A	44"	44"	45.75"	50.5"	44.75"	N/A	51.5"	54"	52"	52.25"	60"	61.25"	54"	67.75" LEFT 66.75" RIGHT	82.5"	N/A	
3-W LENGTH/WIDTH^{4,10}	N/A	46"/20.75"	N/A	N/A	36.75"/19.5"	37"/19.5"	37"/19.5"	40"/21.25"	42"/21.25"	39.75"/20.5"	37"/19"	43.5"/22.25"	45.25"/22.25"	43.5"/22.25"	45.25"/22.25"	N/A	N/A	46"/25.5"	N/A	N/A	59"/24.75"	
4-W LENGTH/WIDTH^{4,10}	45.75"/20"	47"/20.75"	40"/23.5"	28.5" without footplate 33"/25.5"	N/A	39.5"/19.5"	39.5"/19.5"	40.25"/21.25"	43"/21.25"	41"/20.5"	N/A	45.5"/22.5"	47"/22.25"	46"/22.25"	47"/22.25"	47"/22.25"	47"/25.62"	47"/25.5"	53.25"/25.25"	56"/27.25"	N/A	
SEAT-TO-GROUND HEIGHT RANGE	21.5"-23.5"	21"-23.5"	19"	22"-23"	19"-20"	18.25"-20.25"	18.25"-20.25"	20.5"-22"	20.5"-22"	20.75"-22.25"	20.5"	19.25"-20.75"	20.5"-22"	21.5"-24.5"	22"-25"	23.5"-25.5"	23"-25"	22.5"-24.5"	24.75"-26.25"	26.25"-27.75"	28"	
SEAT TYPE (STANDARD)	COMPACT, FOLDABLE (17"x17")	FOLDABLE, MOLDED PLASTIC (20"x18")	COMPACT, FOLDABLE (17.5"x16.5")	COMPACT, FOLDABLE (18"x17")	COMPACT, FOLDABLE (17"x15")	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (17"x17")	COMPACT, FOLDABLE (16"x13")	FOLDABLE, MOLDED PLASTIC (18"x17")	FOLDABLE, MOLDED PLASTIC (18"x17")	FOLDABLE, MOLDED PLASTIC (18"x17")	FOLDABLE, MOLDED PLASTIC (18"x17")	RECLINING H/B W/SLIDERS & HEADREST (18"x18")	RECLINING H/B W/SLIDERS & HEADREST (18"x18")	DELUXE CONTOURED MEDIUM-BACK (22"x18")	RECLINING H/B W/SLIDERS & HEADREST (18"x18")	HIGH-BACK PILLOW TOP W/SLIDERS & HEADREST (19.5"x16")	HIGH BACK (20" X 18") (100" LIMITED RECLINE)	
(3/4-WHEEL) WEIGHT^{3,4,12}	120 LBS.	120 LBS./132 LBS.	54 LBS.	67 LBS.	58 LBS.	76 LBS./80 LBS.	76 LBS./80 LBS.	86 LBS./87 LBS.	93 LBS./99 LBS.	85 LBS./87 LBS.	51 LBS.	117 LBS./131 LBS.	133 LBS./137 LBS.	125.75LBS./136 LBS.	130.25 LBS./139 LBS.	179 LBS.	160 LBS.	171 LBS./176 LBS.	225 LBS.	254 LBS.	208 LBS.	
SEAT-TO-DECK HEIGHT RANGE	16"-18"	15"-18"	14.5"	16.25"-17.25"	16.25"-17.25"	15.75"-17.25"	15.75"-17.25"	16.25"-17.75"	17"-18"	16.25"-17.25"	13"	17.25"-18.75"	17.5"-19"	15.25"-18.25"	15.25" - 18.25"	16.5"-18.5"	16.5"-18.5"	3W: 16"-18" 4W: 16"-18"	17"-18.75"	21"-22.5"	19"	
BATTERY WEIGHT³	13.5 LBS. EA.	25 LBS. EA.	7 LBS. EA.	14 LBS. EA.	9 LBS. EA.	9 LBS. EA. ⁵	9 LBS. EA. ⁵ / 14 LBS. EA. ⁷	14 LBS. EA. ⁷	14 LBS. EA. ⁷	9 LBS. EA. ⁵ / 14 LBS. EA. ⁷	9 LBS. EA. (SEALED LEAD ACID / 5 LBS. EA. (LITHIUM-ION)	24 LBS. EA.	24 LBS. EA./ 32 LBS. EA.	24.5 LBS. EA.	24 LBS. EA./ 32 LBS. EA.	32 LBS. EA.	30.2 LBS. EA.	38 LBS. EA. (NF-22 OR 55 AH) 24 LBS. EA. (U-1)	38 LBS. EA. (NF-22 OR 55 AH) 53 LBS. EA. (75 AH)	53 LBS. EA. (75 AH) 60 LBS. EA. (100 AH)	53 LBS. EA.	
BATTERY REQUIREMENTS⁶	(2) 22AH	(2) U1	(1) 24V, 18AH (LITHIUM ION)	(2) 18AH	(2) 12AH	(2) 12AH	(2) 12AH OR (2) 18AH	(2) 18AH	(2) 18AH	(2) 12AH OR (2) 18AH	(2) 12AH OR (1) 12AH (LITHIUM-ION)	(2) U1	(2) U1 OR (2) 40AH	(2) U1	(2) U1 OR (2) 40AH	(2) 40AH	(2) 50AH	(2) NF-22 (2) U-1 (POWER ELEVATING SEAT)	(2) NF-22 OR (2) GROUP 24	(2) GROUP 27 OR (2) GROUP 24	(2) 75AH	
PER CHARGE RANGE (MAX.)^{2,8}	UP TO 10 MILES	3W: UP TO 12 MILES 4W: UP TO 13 MILES	UP TO 9.6 MILES	UP TO 8.7 MILES	UP TO 6.8 MILES	3W: UP TO 6.9 MILES 4W: UP TO 7.2 MILES	3W: UP TO 8.5 MILES (12AH) UP TO 12.4 MILES (18AH) 4W: UP TO 8.2 MILES (12AH) UP TO 12 MILES (18AH)	3W: UP TO 14.5 MILES 4W: UP TO 14.5 MILES	3W: UP TO 15.9 MILES 4W: UP TO 14.5 MILES	3W: UP TO 8 MILES (12AH) UP TO 11.7 MILES (18AH) UP TO 9.3 MILES (SEALED LEAD ACID)	UP TO 13 MILES (LITHIUM-ION)	UP TO 13 MILES	UP TO 16 MILES	3W: UP TO 17.5 MILES 4W: UP TO 16.5 MILES	3W/4W: UP TO 11 MILES (U-1) 3W: UP TO 13 MILES (40AH) 4W: UP TO 13.5 MILES (40AH)	UP TO 16 MILES	UP TO 19.7 MILES (200 LBS.) UP TO 13 MILES (400 LBS.)	3W: UP TO 18 MILES, U-1 (POWER ELEVATING SEAT) UP TO 9.6 MILES 4W: UP TO 16 MILES	G24: UP TO 19 MILES NF-22: UP TO 14 MILES	G27: UP TO 30 MILES G24: UP TO 22 MILES	19 MILES (UP TO 400 LBS. AND FULL SPEED) 31 MILES (UP TO 220 LBS. AND AT LOW SPEED)	
BATTERY CHARGER	OFF-BOARD, 3.5A	OFF-BOARD, 3.5A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 3.5A	OFF-BOARD, 3.5A	OFF-BOARD, 2A	OFF-BOARD, 2.5A (LITHIUM-ION)	OFF-BOARD, 3.5A (U-1)	OFF-BOARD, 3.5A (U-1) / OFF-BOARD, 5A (40 AH)	OFF-BOARD, 3.5A (U-1) / OFF-BOARD, 5A (40 AH)	OFF-BOARD, 3.5A (U-1) / OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 5A (NF-22) OFF-BOARD, 8A (GROUP 24)	OFF-BOARD, 8A	OFF-BOARD, 8A	
SUSPENSION	FRONT & REAR	FRONT & REAR	LIMITED	LIMITED	N/A	N/A	N/A	N/A	N/A	FRONT & REAR	N/A	N/A	N/A	N/A	N/A	FRONT & REAR	FRONT & REAR	N/A	FRONT & REAR	FRONT & REAR	FRONT & REAR	
MAXIMUM WEIGHT CAPACITY	325 LBS.	375 LBS.	250 LBS.	300 LBS.	250 LBS.	260 LBS.	300 LBS.	300 LBS.	325 LBS.	300 LBS.	250 LBS.	300 LBS.	400 LBS.	300 LBS.	400 LBS.	400 LBS.	400 LBS.	500 LBS.	400 LBS.	400 LBS.	400 LBS.	
STANDARD BODY COLORS	FIRE OPAL, GREEN MOONSTONE, ONYX, BLACK, PEWTER, PEARL WHITE, SAPPHIRE BLUE	GREY STREET, TRUE BLUE	PEWTER	RUBY RED, SAPPHIRE BLUE, CITRINE YELLOW, AMBER ORANGE, PEARL WHITE, ROSE QUARTZ	BLACK	RED	RED, BLUE ¹	RED, BLUE ¹	RED, BLUE ¹	RED, BLUE ¹	WHITE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, OCEAN BLUE	CANDY APPLE RED, OCEAN BLUE	CANDY APPLE RED, TRUE BLUE	CANDY APPLE RED, TRUE BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, SILVER	CANDY APPLE RED, SILVER	SILVER	
WEIGHT OF HEAVIEST PIECE⁴	48 LBS. (REAR SECTION)	48 LBS. (3-WHEEL) (REAR SECTION) 57 LBS. (4-WHEEL) (FRONT SECTION)	N/A	36 LBS. (REAR SECTION)	41 LBS. (3-WHEEL) (SCOOTER BASE)	29 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	29 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	32 LBS. (3-WHEEL) (REAR SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	35 LBS. (3-WHEEL) (FRONT SECTION) 44 LBS. (4-WHEEL) (FRONT SECTION)	34 LBS. (3-WHEEL) (FRONT SECTION) 37 LBS. (4-WHEEL) (FRONT SECTION)	51 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK)	48 LBS. (3-WHEEL) (FRONT SECTION) 59 LBS. (4-WHEEL) (FRONT SECTION)	54 LBS. (3-WHEEL) (FRONT SECTION) 61 LBS. (4-WHEEL) (FRONT SECTION)	51.75 LBS. (3-WHEEL) (FRONT SECTION) 61.5 LBS. (4-WHEEL) (FRONT SECTION)	56.5 LBS. (3-WHEEL) (FRONT SECTION) 64.5 LBS. (4-WHEEL) (FRONT SECTION)	78 LBS. (4-WHEEL) (FRONT SECTION)	77 LBS. (4-WHEEL) (FRONT SECTION)	76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)	95 LBS. (FRONT SECTION)	205 LBS. (MAIN FRAME)	N/A (UNIT DOES NOT DISASSEMBLE)	
HCPCS	K0899*	K0899*	K0899*	K0899*	K0800 - GROUP 1 STANDARD (3W)	K0899*	K0800 - GROUP 1 STANDARD	K0806 - GROUP 2 STANDARD ¹³	K0801 - GROUP 1 HEAVY DUTY	K0899*	K0899*	K0800 - GROUP 1 STANDARD	K0801 - GROUP 1 HEAVY DUTY	K0800 - GROUP 1 STANDARD	K0801 - GROUP 1 HEAVY DUTY	K0899*	K0899*	K0802 - GROUP 1 VERY HEAVY DUTY (3W) K0806 - GROUP 2 VERY HEAVY DUTY (4W)	K0899*	K0899*	K0899*	
WARRANTY	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	5-YEAR LIMITED ON FRAME; 13-MONTH LIMITED ON ELECTRONICS; 13-MONTH LIMITED ON DRIVE MOTORS	5-YEAR LIMITED ON FRAME; 13-MONTH LIMITED ON ELECTRONICS; 13-MONTH LIMITED ON DRIVE MOTORS	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	2-YEAR LIMITED ON FRAME; 18-MONTH LIMITED ON CONTROLLERS; 1-YEAR LIMITED ON ELEC. COMPONENTS & CHARGE	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN

TRAVEL MOBILITY

MID-SIZE LUXURY

HEAVY-DUTY

FULL-SIZE LUXURY