

PRIDE

MOBILITY

Travel Scooter

End-to-End



	i-Go™ Non-Medical Device*	Go Go® Super Portable Non-Medical Device*	Go Go® Carbon FDA Class II Medical Device¹	Go Go® Endurance Li FDA Class II Medical Device¹	Go Go® Ultra X FDA Class II Medical Device¹	Go Go Elite Traveller® FDA Class II Medical Device¹	Go Go Elite Traveller® 2 FDA Class II Medical Device¹	Go Go Elite Traveller® 2 Platinum FDA Class II Medical Device¹	Go Go® Sport FDA Class II Medical Device¹	Revo® 2.0 FDA Class II Medical Device¹
MODEL NUMBER	SC20	SC15	SC25M	S54LXLIT	S39 (3-Wheel) S49 (4-Wheel)	SC40E	SC442E	SC442EMS	S73 (3-Wheel) S74 (4-Wheel)	S66 (3-Wheel) S67 (4-Wheel)
WEIGHT CAPACITY	275 lbs.	300 lbs.	300 lbs.	325 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	325 lbs.	400 lbs.
SPEED (MAXIMUM)²	Up to 4 mph	Up to 4.5 mph	Up to 3.7 mph	Up to 4.1 mph	3-Wheel: Up to 4.3 mph 4-Wheel: Up to 4.1 mph	Up to 4 mph	Up to 4.2 mph	Up to 4.2 mph	Up to 4.7 mph	Up to 5.2 mph
GROUND CLEARANCE⁴	4.2" at mid-frame	1.85" at motor 3.7" at mid-deck	1.5" at motor	1.5" at motor 3.25" at side deck	1.25" at motor 3.2" at mid-frame	0.75" at anti-tip 2.25" at mid-frame	1.34" at transaxle	1.34" at transaxle	1.5" at motor 3" at mid-frame	1.25" at motor
TURNING RADIUS⁴	45.9"	51.62"	53.15	48.25"	36.5"/52.25"	33"	37"	37"	3-Wheel: 36" 4-Wheel: 41.5"	3-Wheel: 40" 4-Wheel: 52"
LENGTH/WIDTH⁴,¹⁰	44"/20.3"	37.2"/21.54"	40.35"/17.52"	42"/20.5"	3-Wheel: 37"/19.25" 4-Wheel: 40.25"/19.25"	37"/19.5"	39.5"/19.5"	39.5"/19.5"	3-Wheel: 42"/21.25" 4-Wheel: 43"/22.25"	3-Wheel: 46"/20.75" 4-Wheel: 47"/20.75"
FRONT TIRES⁴	8" hub motor	7" solid	7" solid	8" solid	8" solid	8" solid	7" solid	7" solid	3-Wheel: 8" solid 4-Wheel: 9" solid	3-Wheel: 8" solid 4-Wheel: 9" solid
REAR TIRES⁴	9" perforated	8" solid	8" solid	8" solid	8" solid	8" solid	8" solid	8" solid	9" solid	9" solid
SUSPENSION	N/A	N/A	N/A	Front & Rear	N/A	N/A	N/A	N/A	N/A	Front & Rear
RANGE PER CHARGE²,⁸	Up to 9.3 miles (220 lbs.)¹⁵ Up to 7 miles (275 lbs.)¹⁵	Up to 10.1 miles (220 lbs., 11.7Ah) Up to 8.8 miles (300 lbs., 11.7Ah)	Up to 10 miles	Up to 6.6 miles (8Ah, 200 lbs.) Up to 5 miles (8Ah, 325 lbs.)¹⁵ Up to 9.9 miles (16Ah, 200 lbs.) Up to 13.2 miles (16Ah, 325 lbs.)¹⁵	3-Wheel: Up to 8 miles (200 lbs.) Up to 6.2 miles (300 lbs.) 4-Wheel: Up to 8.9 miles (200 lbs.) Up to 6.6 miles (300 lbs.)	Up to 7.3 miles (12Ah, 300 lbs.) Up to 10.5 miles (18Ah, 300 lbs.)	Up to 6.6 miles (12Ah, 300 lbs.) Up to 10.4 miles (18Ah, 300 l bs.)	Up to 6.6 miles (12Ah, 300 lbs.) Up to 10.4 miles (18Ah, 300 l bs.)	3-Wheel: Up to 12.75 miles (325 lbs.) 4-Wheel: Up to 12 miles (325 lbs.)	3-Wheel: Up to 17.4 miles (200 lbs.) Up to 12.8 miles (400 lbs.) 4-Wheel: Up to 17.8 miles (200 lbs.) Up to 14.8 miles (400 lbs.)
TOTAL WEIGHT WITHOUT BATTERIES³,⁴,¹¹	60.6 lbs.	41.6 lbs.	35 lbs.	93 lbs.	3-Wheel: 76.3 lbs. 4-Wheel: 81.5 lbs.	75.5 lbs.	83 lbs.	83 lbs.	3-Wheel: 92.75 lbs. 4-Wheel: 100 lbs.	3-Wheel: 120 lbs. 4-Wheel: 132 lbs.
WEIGHT OF HEAVIEST PIECE⁴	32.8 lbs. (Rear Base with Seat Back)	N/A (unit does not disassemble)	N/A (unit does not disassemble)	40.75 lbs. (Front Section)	3-Wheel: 29.5 lbs. (Front Section) 4-Wheel: 34.5 lbs. (Front Section)	29.5 lbs. (Front Section)	35 lbs. (Front Section)	35 lbs. (Front Section)	3-Wheel: 35.75 lbs. (Front Section) 4-Wheel: 44 lbs. (Front Section)	3-Wheel: 48 lbs. (Rear Section) 4-Wheel: 57 lbs. (Front Section)
STANDARD SEATING	Compact, Foldable with Removable Back (16" x 16")	Foldable, (15" x 15")	Foldable, 14" x 11.5"	Mesh Back (17" x 17")	Compact, Foldable (17" x 17")	Compact, Foldable (17" x 17")	Mesh Back (17" x 17")	Mesh Back (17" x 17")	Compact, Foldable (17" x 17")	Compact, Foldable (20" x 18")
BATTERY REQUIREMENTS	25.2V, 9.6Ah (Lithium-ion)	11.7Ah (Lithium-ion) or 16Ah (Lithium-ion)	12Ah (Lithium-ion)	8Ah (Lithium-ion) or 16Ah (Lithium-ion)	(2) 12Ah (SLA/AGM)	(2) 12Ah (SLA/AGM) or (2) 18Ah (SLA/AGM)	(2) 12Ah (SLA/AGM) or (2) 18Ah (SLA/AGM)	(2) 12Ah (SLA/AGM) or (2) 18Ah (SLA/AGM)	(2) 18Ah (SLA/AGM)	(2) U-1 (SLA/AGM)
BATTERY WEIGHT⁹	4.8 lbs.	8 lbs. (11.7Ah) 9.8 lbs. (16Ah)	4 lbs.	8 lbs. (8Ah)/ 12 lbs. (16Ah)	9 lbs. ea.⁵	9 lbs. ea.⁵/ 14 lbs. ea.⁷	9 lbs. ea.⁵/ 14 lbs. ea.⁷	9 lbs. ea.⁵/ 14 lbs. ea.⁷	14 lbs. ea.⁷	25 lbs. ea.
BATTERY CHARGER	Off-board, 3.5A	Off-board, 3.5A	Off-board, 2A	Off-board, 3.5A	Off-board, 2A	Off-board, 2A	Off-board, 2A	Off-board, 2A	Off-board, 2A	Off-board, 3.5A
SEAT-TO-DECK HEIGHT RANGE	14.6"	16.61"	16.14"	16.25"-17.75"	16.5"-18"	15.75"-17.25"	16.38"-17.83"	16.38"-17.83"	17"-18"	16.5" - 18"
SEAT-TO-GROUND HEIGHT RANGE⁴	22.25"	21.42"	20.87"	21.25"-22.75"	21.5"-23"	18.25"-20.25"	21.38"-22.83"	21.38"-22.83"	20.5"-22"	22"-23.5"
STANDARD BASE COLORS	Sugar Plum (matte) Iceberg Blue (matte) Black (matte) White (matte)	Black Mint	Black Carbon Fiber	Garnet Red (Matte) Euro Grey ¹⁴	Black	Red Blue¹⁴	Cardinal Red, Black, Iceberg Blue, White Robin's Egg Blue, Chell-O, Sugar Plum, Tangerine	Cardinal Red, Black, Iceberg Blue, White Robin's Egg Blue,Chell-O, Sugar Plum, Tangerine	Red Blue¹⁴	Grey Street True Blue
HCPCS CODE	Non-Medical*	Non-Medical*	K0899**	K0899**	K0800 - Group 1 Standard	K0800 - Group 1 Standard	K0800 - Group 1 Standard	K0899**	K0801 - Group 1 Heavy Duty	K0899**
WARRANTY	3-yr limited on Frame 2-yr limited on Electronics 1-yr limited on Drivetrain, Motor, Battery, & Charger	3-yr limited on Frame 2-yr limited on Electronics 1-yr limited on Drivetrain, Motor, Battery, and Charger	Limited Lifetime on Frame 1-yr limited on Electronics, Drivetrain and Battery	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 2-yr limited on Battery	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime Limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

2. Varies with user weight, terrain type, battery amp-hour (Ah), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

3. Battery weight not included.

4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

5. Battery pack is 19.5 lbs. with 12AH batteries.

6. Deep cycle (AGM or gel cell type recommended).

7. 18 AH battery pack weighs 29.5 lbs.

8. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 71764 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

9. Battery weight may vary based on manufacturer.

10. Width of base (not including seat).

11. Includes standard seat weight.

12. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.

13. Light package.

14. Interchangeable shrouds (all shroud panels included).

15. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

•Non-Medical Devices: This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.

•K0899: Power mobility device, not coded by DME PDAC.

DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or – 3 %.