

THE MOST COMFORTABLE RECLINER IN THE WORLD!



COLLECTION	PREMIER ^{1,2,3}	ULTRA ^{1,2,3}	ULTRA ^{1,2,3}	ULTRA ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	TRANQUIL 2¹	TRANQUIL 2¹	TRANQUIL 2¹	METRO 2¹
MODEL	PLR-5900	PLR-4955S	PLR-4955M	PLR-4955L	PLR-3955PW	PLR-3955S	PLR-3955M	PLR-3955LT	PLR-935S	PLR-935M	PLR-935LT	PLR-925M
SIZE	Medium	Small	Medium	Large	Petite Wide	Small	Medium	Large/Tall	Small	Medium	Large/Tall	Medium
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.
POSITION**	True-Infinite with Tilt 5'5" - 6'1"	True-Infinite with Tilt 5'4" and below	True-Infinite with Tilt 5'5"-6'	True-Infinite with Tilt 6' and above	True-Infinite 5'4" and below	True-Infinite 5'4" and below	True-Infinite 5'5"-6'1"	True-Infinite 6'2" and above	True-Infinite 5'3" and below	True-Infinite 5'4"-5'9"	True-Infinite 5'10" and above	True-Infinite 5'3" and below
OVERALL WIDTH	39.5"	37"	39"	39"	41"	37"	39"	40"	33.75"	37.5"	38"	33"
CHAIR WEIGHT	232 lbs.	158 lbs.	158 lbs.	176 lbs.	161 lbs.	156 lbs.	158 lbs.	178 lbs.	137 lbs.	145 lbs.	155 lbs.	142 lbs.
CHAIR HEIGHT	43"	41"	43"	47"	40"	41"	43"	49"	40"	42"	47"	41"
SEATING TYPE	Pocketed Coiled Springs	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort
SEAT WIDTH	23"	20"	22.5"	23"	24"	20"	22.5"	23.5"	20"	22"	23.5"	21"
SEAT DEPTH	21"	18"	21"	22"	18"	18"	21"	22"	18.5"	20.5"	22.5"	21"
SEAT TO FLOOR HEIGHT	20.5"	18"	19.5"	21.5"	19.5"	18"	19.5"	21.5"	18"	19.5"	21.5"	19"
TOP OF BACK TO SEAT	26"	25"	28"	30"	23"	25"	28"	31"	25"	27"	30"	26.5"
DISTANCE FROM WALL	17"	15"	15"	21.5"	14"	15"	15"	21.5"	13"	16"	16"	14"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard
HEAT	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	NA	NA	NA	NA
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

STANDARD FABRIC SELECTION GUIDE

SORRENTO

ECLIPSE

COFFEE

CAPRICCIO

SLATE

CAPPUCCINO

SMOKE

DOVE

SORRENTO

COFFEE

TRENTO⁴

BROWN

CHARCOAL

CANYON

OCEAN

SILT

STEEL

WALNUT

ASTRO

GREY

BROWN

MUSHROOM

SAVILLE

GREY

BROWN

STANDARD FEATURES:

POWER HEADREST

POWER LUMBAR

USB REMOTE

LITHIUM BATTERY BACKUP

FOOTREST EXTENSION

POSITION GUIDE:

RELAXED

RECLINED

UPWARD LIFT

TRUE INFINITE

= Indicates shortest available lead time.

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.
3. Unit has three different heating pads in the shoulders, lumbar and seat.
4. Available for the PLR-4955M and PLR-4955L only.

*Not all units achieve all positions. Please refer to specific product.
**Height ranges should be used as a guide only due to varying factors of each individual user.
Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.
The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.